

Table S2. Gene-specific primers used in genotyping experiments:

Gene	Forward primer (5'-3')	Reverse primer (5'-3')
<i>Otub1</i> WT allele	ACTTGACCTTGGACAAGATGGAG	TGAAGCCCTAGCAACCTCAAG
<i>Otub1</i> KO allele	GAGATGGCGCAACGCAATTAATG	TACCTGTGTGACCACCTGAAG
<i>NIK</i>	ATCAAGCTGGCCCTTAACCT	CAAGGAGTTCTTGTTTCCCAG
<i>Lym1</i>	CTGCTGAAAGGTCCTGAGAC	CATTCAGTGCACCTGAGGCTG
CD19-Cre WT	AGAGGGAGGCAATGTTGTGC	GTCCAGGTCCCTGACGTCTG
CD19-Cre KO	AGAGGGAGGCAATGTTGTGC	GACGATGAAGCATGTTTAGCTGG
Rosa26-CreER WT	AAAGTCGCTCTGAGTTGTTAT	GGAGCGGGAGAAATGGATATG
Rosa26-CreER KO	AAAGTCGCTCTGAGTTGTTAT	CCTGATCCTGGCAATTTTCG