

Supplementary Information for

Paste manuscript title here

Paste the full author list here

Paste corresponding author name here Email: ruvkun@molbio.mgh.harvard.edu

This PDF file includes:

Supplementary Information Figures S1 to S7 Tables S1 to S4 Legends for Datasets S1 SI References

Other supplementary materials for this manuscript include the following:

Datasets S1

Supplementary Information Text

Supplemental Materials and methods

C.elegans strains used in this study

The following strains were used in the study:

Wildtype	N2
SJ4100	zcls13[hsp-6::GFP]
SJ4058	zcls9[hsp-60::GFP]
SJ4005	zcls4 [hsp-4::GFP] V
AU133	agIs17 [myo-2p::mCherry + irg-1p::GFP] IV
YL162	atfs-1(tm4525); zcls9[hsp-60::GFP]
QC115	atfs-1(et15) V
DA2249	spg-7 (ad2249) I
MQ130	clk-1(qm30) III
CW152	gas-1(fc21) X
MQ887	isp-1(qm150) IV
TK22	mev-1(kn1) III
CU5991	fzo-1(tm1133)
CU6372	drp-1(tm1108)
QC134	nduf-7(et19)
VC1789	nduf-2.2(ok2397)
MQ1333	nuo-6(qm200)
NL1832	ucr-2.3(pk732)

Genomewide keio screen for developmental delay or *hsp-6::gfp* induction phenotype

Keio mutants which represent ~3800 genes were grown in LB media containing kanamycin at 37°C overnight. 100 μ l of the overnight culture were dropped onto NGM media containing 24-well plates and grown overnight at 37°C. Synchronized L1-stage transgenic *hsp-6::gfp* animals prepared by bleach-prep were dropped onto these 24-well Keio mutant in NGM media plates. The wells were screened visually at day 3 under a fluorescent microscope for the induction of *hsp-6::gfp*. All the wells were screened visually at day 5 under the DIC microscope for wells that contained animals growing significantly slower compared to control wells grown on wildtype BW25113 *E. coli*. The primary hits were rescreened again at least three times.

Removal of Kanamycin cassette in ∆fes::Kan Keio mutant and construction of double mutants

Because bacterial insertion drug resistant cassette mutations in one gene in an operon are expected to be polar on downstream genes, the drug resistant cassette alleles are best for identifying relevant operons rather than individual genes. But because genes in operons often act in the same pathway, the identification of the operon is informative (Table S1). Mutations that disrupt the *E. coli* cytochrome oxidase *cyoABCDE* operon were previously shown to induce *C. elegans hsp-6:GFP* (1) ; these were among the strongest hits of the 24 mutants that caused developmental delay (Table 1). None of the other *E. coli* genes identified in the screen of the Keio mutant collection act in the *E. coli* electron transport chain, consistent with our previous finding that inactivation of only specific steps in bacterial respiration causes defects in *C. elegans* growth or mitochondrial function (1). Keio mutants have a Kanamycin cassette flanked by FLP recombinase sites, which allows for the excision of the cassette in the

3

place of gene of interest (2). Kanamycin cassette incision in $\Delta fes::Kan$ was achieved by introducing the plasmid pCP20 by transformation, which contains an Flp recombinase. Transformants were selected for ampicillin resistance at 30°C and the raised at 42°C to remove the pCP20 plasmid. Colonies were verified by inability to grow on LB plates with Kanamycin as well as by sequencing for the recombination site junction.

Double mutants, $\Delta fes entA::kan$, $\Delta fes entB::kan$, $\Delta fes cirA::kan$, $\Delta fes fiu::kan$, or $\Delta fes fepA::kan$ were generated as follows: P1 lysates of Keio collection mutants, entA::kan, entB::kan, cirA::kan, fiu::kan, or fepA::kan strains were generated by growing overnight cultures of each strain in LB media containing 2mM MgSO₄,10mM CaCl₂, 0.2% glucose at 37°C until the OD₆₀₀ is around 0.5. 100 µl of Phage P1 obtained form CGSC was added to 1 ml of each culture. After 2 hours incubation at 37°C until complete lysis, 0.5 ml of chloroform was added the culture. The lysates were centrifuged at 13K RPM for 1 min and the lysates were stored in 4°C until use. 100 µl of P1 lysate from each of the Keio mutant was added to overnight cultures of Δfes mutants grown in LB was suspended in 1ml LB media containing 5mM MgSO₄, 50mM CaCl₂ and incubated for 30 min. The tubes were centrifuged at maximal speed for 1 min and plated on kanamycin containing LB Media agar plates. The colonies were verified for the presence of Δfes and Keio mutant insertions using sequencing.

In vitro ATP measurement

ATP measurements were conducted as per the protocol with the following modifications (3). 1000 animals in 100 μ l of M9 were frozen in liquid nitrogen and stored at -80°C until analysis. The samples were placed immediately in boiling water bath for 15 min. After adding twice amount of distilled water, centrifugation was done at 14,000 rpm for 5 min and supernatant is transferred into a second tube. Equal volume of luciferase reagent

4

(Promega[™]) was added and luminescence was measured using a plate reader (SpectraMax Microplate reader). Protein content was measured using BCA method. ATP concentrations were normalized to protein content of samples. All the values were normalized to wildtype grown on *E. coli* BW25113.

In vivo ATP measurement

The transgenic reporter strain PE254 (*fels4[Psur-5::luc+::gfp; rol-6 (su1006)] V*) (4) was used as a surrogate for determining the intracellular ATP content. ~1000 synchronized PE254 animals grown in *E. coli* BW25113 or Δfes mutant until the L4-larval stage were collected from NGM media plates and washed at least 4 times in M9 buffer to remove bacteria. The animals were transferred to white flat-bottom, 96-well plate wells containing 100 μ M D-luciferin in M9 buffer and the luminescence was measured within 5 min in SpectraMax microplate reader. GFP fluorescence was quantified in SpectraMax microplate reader. GFP fluorescence values were normalized to GFP fluorescence of samples. All the values were normalized to wildtype grown on *E. coli* BW25113.

Measuring endogenous hydrogen peroxide levels using HyPer strain

To measure endogenous hydrogen peroxide levels, transgenic *jrls1*[P*rpl-17*::HyPer] animals expressing a HyPer as described (5). Synchronized L1-stage transgenic *jrls1*[P*rpl-17*::HyPer] animals were grown in *E. coli* BW25113 or Δfes mutant and about 1000 L4-stage animals were harvested in 96 microtiter well plates. For H₂O₂ treatment experiment, H₂O₂ was added to the animals in the 96-well itself and measurements were taken within 15 minutes. An excitation wavelength of either 490 nm or 405 nm was used to measure oxidized or reduced HyPer probe fluorescence respectively with an emission filter of 535 nm. The absorbance at 620 nm was used to

5

normalize for animal numbers. The statistical significance of differences between conditions was determined by using unpaired *t*-test. GraphPad Prism 8.0 was used for these calculations.

Staining of mitochondria

Animals were stained with 10μ M of 10-nonyl acridine orange bromide (NAO) in M9 solution for 1 hour and destained by washing in M9 atleast 5 times.

Antioxidant treatments

Stock solutions of trans-resveratrol, L-ascorbic acid, cis-resveratrol, oxyresveratrol, cistrimethoxyresveratrol, trolox, or N-acetyl cysteine were diluted in M9 to a working concentration of 10μ g/ml and dropped onto NGM media plates seeded with overnight cultures of Δfes mutant *E. coli*. Synchronized L1-stage *spg-7(ad2249)* animals were added to the plates and scored after three days for the suppression of arrest phenotype. The negative control NGM plates contained the solvent and M9.

Microscopy and imaging

Nematodes were mounted onto agar pads and images were taken using a Zeiss AXIO Imager Z1 microscope fitted with a Zeiss AxioCam HRm camera and Axiovision 4.6 (Zeiss) software. Student's t test was used determine statistical significance. Low-magnification bright-field images were acquired using a Zeiss AxioZoom V16 Microscopre, equipped with a Hammamatsu Orca flash 4.0 digital camera camera, and using Axiovision software. Some Immunofluorescence images were acquired with an IX-70 microscope (Olympus, Waltham, MA) fitted with a cooled CCD camera (CH350; Roper Scientific) driven by the Delta Vision system (Applied Precision, Pittsburgh, PA). Images were deconvolved using the SoftWoRx 3.3.6 software (Applied Precision).

Quantification and Statistical analysis

Statistical analysis was performed with GraphPad Prism 7 software. Student's t-test was used for analysis.



Fig. S1: Screen methodology

Individual Keio mutant *E. coli* strains were grown in individual wells of 24 well NGM agar plates, which supports bacterial growth on the surface, and the bacteria that grew on the agar surface were inoculated the next day with ~30 synchronized L1 larval stage *hsp-6::gfp* animals. Each individual well was then screened at two different time-points: 48 hours at 20°C for the Keio mutant strains that induced *hsp-6::gfp* and at 96 hours at 20°C for *E. coli* mutant strains that dramatically slowed *C. elegans* development.

Fig.	S2
------	----



Fig. S2: E. coli mutants induce hsp-6::gfp

While *hsp-6::gfp* is not induced in animals grown on wildtype *E. coli* BW25113, animals grown on mutant Keio *E. coli* strains showed *hsp-6::gfp* activation.



Fig. S3: Diagrammatic representation of bacterial enterobactin content

- A) Diagrammatic representation of ∠fes::kan, ∠fepD::kan, ∠fepG::kan, ∠fepC::kan mutant accumulating ferric enterobactin
- B) Diagrammatic representation of *∆fes* mutant accumulating ferric-enterobactin and lack of enterobactin in *∆fes ∆entA::kan* double mutants.



Fig. S4: Determining the arrest stage of *spg-7(ad2249)* animals grown on Δfes -722 mutant

A) The size and appearance of L1-larval stage wildtype or *spg-7(ad2249)* animals from egg-prep looks very similar to that of the arrested *spg-7(ad2249)* animals grown on *fes*-722 mutant *E. coli* for 4 days. One-day old wildtype animals used for comparison was grown on *fes*-722 mutant *E. coli* from L1-stage for 4 days. Scale bar, 100μM.

- B) Quantification of the animal length in L1-larval stage wildtype or *spg-7(ad2249)* animals from egg-prep and arrested *spg-7(ad2249)* animals grown on *fes-722* mutant *E. coli* for 4 days. One-day old wildtype animals used for comparison were grown on *fes-722* mutant *E. coli* from L1-stage for 4 days. Statistical significance was determined using one-way ANOVA. **** P<0.0001. ns denotes not significantly different. The number of animals used is shown above each bar in the graph. The error bars indicate standard deviation. Data represent one out of two independent experiments.</p>
- C) Whole-mount DAPI staining of animals from L1-larval stage wildtype or *spg-*7(*ad2249*) animals from egg-prep and arrested *spg-7*(*ad2249*) animals grown on *fes-722* mutant *E. coli*. One-day old wildtype animals were used for comparison. The germline in the L1-larval stage wildtype or *spg-7*(*ad2249*) animals from egg-prep and arrested *spg-7*(*ad2249*) animals grown on *fes-722* mutant *E. coli* were not apparently visible at this stage because at this stage there are only 2 germline progenitor cells. If the arrested worms were in their L2 or L3 stages they would have apparently visible germ cells. The germline from adult wildtype animal (red box) is shown for comparison.



Fig. S5: Interaction between mitochondrial mutants and ∆fes-722 mutant

E. coli mutants cause development arrest phenotype in *nduf-7(et19)*, *isp-1(qm150)*, and *clk-1(qm30)* mutants. No development arrest phenotype was observed in *gas(fc21)* and *ucr-2.3(pk732)* mutants grown on Δ *fes-722*.

Fig. S6



Figure S6: Measuring hydrogen peroxide levels

Endogenous hydrogen peroxide levels were not significantly different between animals grown on *E. coli* BW25113 or Δfes -722 *E. coli*; however, addition of exogenous peroxide results in significantly higher level of endogenous peroxide levels in animals grown on Δfes -722 *E. coli* compared to animals grown on *E. coli* BW25113.

Fig. S7

Α



в

E. coli 1 KCIYVAIGQKASTISNVVRKL Gammaproteobacteria 1 YCIYVAIGQKRSTVAQIVKRL C. elegans 1 FCIYVAVGQKRSTVAQIVKRL Kordiimonas gwangyangensis 1 Kordiimonas lipolytica 1 Altererythrobacter sp. ZODW24 1

Figure S7: ATP Synthase protein sequence analysis

- A) C. elegans ATP-1 is highly homologous to *E. coli* ATP synthase. The black boxes shows the 21 amino acid peptide sequence that was shown to bind enterobactin
- B) The 21 amino acid peptide sequence shown to bind enterobactin is highly

conserved

Table S1: Table showing the *spg-7(ad2249)*, *nduf-7(et19)*, *isp-1(qm150)*, *clk-1(qm30)* tested on N-Acetyl Cysteine antioxidant for suppression of arrest phenotype

	arrest (N)	
Animal Strain	Fed on <i>E. coli ∆fes-</i> 722; Solvent control	Fed on <i>E. coli ∆f</i> es- 722; NAC
Wildtype	0±0 (142)	0±0 (161)
spg-7(ad2249)	100±0 (147)	0.69±1.2 (151)
nduf-7(et19)	100±0 (142)	1±1.8 (129)
isp-1(qm150)	100±0 (121)	0±0 (108)
clk-1(qm30)	100±0 (155)	0.7±1.3 (144)

Data from three independent trials were collected; the average of percentage of L1 animals and the standard deviation are shown. The total number of animals counted in all the three independent trials is shown in the parentheses.

Table S2: *spg-7(ad2249)* animals grown on *∆fes-722* mutant *E. coli* killed with

gentamicin or tetracycline arrest as L1-larvae.

	spg-7(ad2249) % Larval arrest (N)
Fed on BW25113	0±0 (161)
Fed on <i>∆fes-722::Kan</i>	100±0 (188)
Fed on BW25113, Ampicillin-treated	0.5±0.8 (176)
Fed on <i>∆fes-722::Kan</i> , Ampicillin-treated	0±0 (146)
Fed on BW25113, Tetracycline-treated	0±0 (144)
Fed on <i>∆fes-722::Kan</i> , Tetracycline-treated	98.6±1.1 (155)
Fed on BW25113, Gentamicin-treated	0.7±1.2 (144)
Fed on <i>∆fes-722::Kan</i> , Gentamicin-treated	98.4±1.7 (195)
Fed on BW25113, Kanamycin-treated	0±0 (134)
Fed on <i>∆fes-722::Kan</i> , Kanamycin-treated	97.7±2.2 (175)
Fed on ∆ <i>fes-722::Kan</i> + BW25113 (75:25)	0±0 (119)
Fed on <i>∆fes-722::Kan</i> + BW25113 (75:25), Kanamycin-treated	98±1.9 (151)
Fed on <i>∆fes-722::Kan</i> + BW25113 (50:50)	0.65±1.1 (146)
Fed on <i>∆fes-722::Kan</i> + BW25113 (50:50), Kanamycin-treated	97.8±2.4 (142)
Fed on <i>∆fes-722::Kan</i> + BW25113 (25:75)	0±0 (138)
Fed on <i>∆fes-722::Kan</i> + BW25113 (25:75), Kanamycin-treated	98.7±2.1 (166)
Fed on <i>∆fes-722::Kan</i> + BW25113 (100:1)	0±0 (141)
Fed on <i>∆fes-722::Kan</i> + BW25113 (100:1), Kanamycin-treated	99±1.6 (166)
Fed on <i>∆fes-722::Kan</i> + BW25113 (1000:1)	0±0 (152)
Fed on ∆fes-722::Kan + BW25113 (1000:1), Kanamycin-treated	98.5±2.4 (140)
Fed on <i>∆fes-722::Kan</i> + BW25113 (10,000:1)	0±0 (124)
Fed on <i>∆fes-722::Kan</i> + BW25113 (10,000:1), Kanamycin- treated	99.4±0.9 (163)

Data from three independent trials were collected; the average of percentage of L1 animals and the standard deviation are shown. The total number of animals counted in all the three independent trials is shown in the parentheses.

Table S3: Measurement of CFU/mI of wildtype *E. coli* BW25113 and *E. coli* ∆fes-

722::Kan in NGM plates

	BW25113 (CFU/ml)	fes::Kan	Ratio
		(CFU/ml)	
Day 0	1x10 ⁴	1x10 ⁸	10,000
Day 1	2.6x10 ⁷	5.6x10 ⁸	~200
Day 2	6x10 ⁸	4.6x10 ⁸	~1.3

Table S4: Trans-resveratrol treatment acts through antagonizing animal ROS pathways

spg-7(ad2249) animals fed on Δfes -722	% Larval arrest (N)	
solvent control	100±0 (105)	
Trans-Resveratrol, 10µg/ml	0±0 (123)	
Trans-Resveratrol, 10µg/ml, washed	100±0 (143)	

spg-7(ad2249) animals grown on Δfes -722 mutant grown in trans-resveratrol arrest do not arrest while animals grown on Δfes -722 mutant *E. coli* pre-grown in the presence of trans-resveratrol and washed to remove trans-resveratrol arrest as L1-larvae.

Supplementary Dataset 1

Table showing list of the hits from the keio screen for mutants that induces *hsp-6::gfp*. All these mutant were found to induce *hsp-60::gfp*. None of these mutantsactivated *hsp-4::gfp* or *pgp-5::gfp*.

References

- 1. Govindan JA, Jayamani E, Zhang X, Mylonakis E, Ruvkun G (2015) Dialogue between E. coli free radical pathways and the mitochondria of C. elegans. *Proc Natl Acad Sci U S A*. doi:10.1073/pnas.1517448112.
- 2. Baba T, et al. (2006) Construction of Escherichia coli K-12 in-frame, singlegene knockout mutants: the Keio collection. *Mol Syst Biol* 2:2006–0008.
- 3. Palikaras K, Tavernarakis N (2016) Intracellular Assessment of ATP Levels in Caenorhabditis elegans. *Bio Protoc* 6(23). doi:10.21769/BioProtoc.2048.
- 4. Lagido C, McLaggan D, Glover LA (2015) A Screenable In Vivo Assay for Mitochondrial Modulators Using Transgenic Bioluminescent Caenorhabditis elegans. *J Vis Exp* (105):e53083–e53083.
- 5. Back P, et al. (2012) Exploring real-time in vivo redox biology of developing and aging Caenorhabditis elegans. *Free Radic Biol Med* 52(5):850–859.