

## Supplementary Material to “Fibrin biopolymer as scaffold candidate to treat bone defects in rats”

**Additional file 1** – Surface cell markers used on mesenchymal stem cell characterization by flow cytometry.

<b>MARKER</b>	<b>EXPRESSION</b>	<b>TRADEMARK</b>	<b>REF.CODE</b>	<b>DILUTION</b>	<b>FLUOROCHROME</b>
CD 45	Negative	ABCAM	AB123522	1/100	FITC
CD 11b	Negative	ABCAM	AB33816	1/100	PE
MHC- II	Negative	ABCAM	AB111825	1/100	PE
RT1- Aw2	Negative	ABCAM	AB119770	1/100	FITC
CD 90	Positive	ABCAM	AB226	1/100	FITC
ICAM1	Positive	ABCAM	AB23835	1/100	FITC