



Supplemental Fig. S3. Expression and purification of the recombinant His-Blp1₂₆₅₂₋₃₃₆₂ C-terminal fragment. Fractions of flow-through, column washing, elution and dialysis in PBS buffer were fractionated in 12% SDS-PAGE. The position of His-Blp1₂₆₅₂₋₃₃₆₂ C-terminal fragment is indicated by the arrow. The positions of standard molecular mass markers (PageRuler™ Prestained Protein Ladder (Thermo Fisher Scientific)) are shown on the left.