

Supporting Information Table S2. Antibodies used in WB of this study.

Antibody	Company	Dilution
TAK1 (#4505)	Cell Signalling Technology	1/ 1,000 in Immuno-enhancer
Phospho-TAK1 (Thr184 / 187) (#4531)	Cell Signalling Technology	1/ 1,000 in Immuno-enhancer
YAP1 (G-6, sc-376830)	Santa Cruz Biotechnology, Inc.	1/ 1,000 in Immuno-enhancer
Phospho-YAP1 (Ser397) (D1E7Y, #13619)	Cell Signalling Technology	1/ 1,000 in Immuno-enhancer
TAZ (M2-616, 560235)	BD biosciences	1/ 1,000 in Immuno-enhancer
α Tubulin (B-7, sc- 5286)	Santa Cruz Biotechnology, Inc.	1/ 1,000 in Immuno-enhancer
TBP (#8515)	Cell Signalling Technology	1/ 1,000 in Immuno-enhancer
HA (F-7, sc-7392)	Santa Cruz Biotechnology, Inc.	1/ 5,000 in 0.2% Twenn-TBS containing 10% Block Ace
Ub (P4D1, sc-8017)	Santa Cruz Biotechnology, Inc.	1/ 1,000 in Immuno-enhancer
Caspase-3 (#9662)	Cell Signalling Technology	1/ 1,000 in Immuno-enhancer
Gapdh (FL-335, sc- 25778)	Santa Cruz Biotechnology, Inc.	1/ 5,000 in 0.2% Twenn-TBS containing 10% Block Ace