

Supplementary Table 1. Primer sets used for real-time PCR.

Gene	Sense Primer (5' - 3')	Antisense Primer (5' - 3')
Cyclin D1 (Mouse)	GCGTACCCTGACACCAATCTC	ACTTGAAGTAAGATACGGAGGGC
p15ink4b (Mouse)	ATGTTGGGCGGCAGCAGTGACG	ATCTCCAGTGGCAGCGTGACG
p18ink4c (Mouse)	GATTTGGGAGAACTGCGCTG	GCCTGGA ACTCCAGCAAAGC
α -SMA (Mouse)	AGAACACGGCATCATCAC	TCCAGAGTCCAGCACAAT
Collagen-III (Mouse)	TCCTGAAGATGTCGTTGATGTG	GTCTCATTGCCTTGCGTGTT
p19af ^r (Mouse)	GGGTCGCAGGTTCTTGGTC	GTGCGGCCCTCTTCTCAA
TGF- β 1 (Rat)	ACCGCAACAACGCAATCTA	AATGGGGGTTCTGGCACT
GAPDH (Mouse)	CATGGCCTTCCGTGTTCCCTA	GCCTGCTTCACCACCTTCTT
GAPDH (Rat)	GGCAAGTTCAACGGCACA	TCTCGCTCCTGGAAGATGG