


Supplemental information

Surmiak et. Al. LTB4 and 5-oxo-ETE from extracellular vesicles stimulate neutrophils in granulomatosis with polyangiitis



Supplemental Figure S1. Gating strategy for flow cytometry analysis of EV coated aldehyde beads. Based on forward and side scatter (FS/SS), beads were first gated for singlets (Panel A). Histograms (FS) of aldehyde -beads before gating (Panel B) and after gating (Panel C). Panel D – CD81 expression of EV coated aldehyde-beads.