

S1 Table. Fractionation conditions of the GELFREE 10% Mass Cartridge Kit

Step	1	2	3	4	5	6	7	8	9	10	11	12	13
Voltage (V)	50	50	50	50	50	100	100	100	100	100	100	100	100
Fraction interval (min)	16	36	2	3	4	3	4	5	7	10	15	20	35
Total elapsed time (min)	16	52	54	57	61	64	68	73	80	90	105	125	160
Fraction No.	---	1	2	3	4	5	6	7	8	9	10	11	12