YMTHE, Volume 28

Supplemental Information

HucMSC-Derived Exosomes Mitigate

the Age-Related Retardation of Fertility

in Female Mice

Weijie Yang, Jing Zhang, Boqun Xu, Yuanlin He, Wei Liu, Jiazhao Li, Songying Zhang, Xiaona Lin, Dongming Su, Tinghe Wu, and Jing Li

Supplementary Figure Legend 1


Figure S1. Specific effects of HucMSC-exos on follicular activation and 2 development. (A) Western blot showing the expressions of p-Akt (Ser473), p-mTOR 3 (Ser2448), p-rpS6 (Ser235/236) in ovaries treated with or without HucMSC-exos for 4 24 h. Ovaries co-cultured with HucMSC were used as positive control. The 5 expressions of Akt, mTOR, rpS6, and β-actin were used as internal controls. Ctrl, 6 PBS; Exo, HucMSC-exos 20 µg/mL; MSC, HucMSC. (B) Western blot of the 7 expression of p-Akt (Ser473), p-mTOR (Ser2448), p-rpS6 (Ser235/236) in isolated 8 primordial oocytes treated with HucMSC-exos with/without Akt inhibitor, Akt VIII 9 (AI) and mTOR inhibitor, rapamycin (Ra) for 24 h. (C and D) Blocking effects of 10 mTOR knockdown on HucMSC-exos induced activation of PI3K/mTOR signaling 11 pathway. Newborn ovaries were pretreated with mTOR siRNAs for 60 h, following 12 with 24 h of treatment with or without HucMSC-exos. RT-PCR showing the effective 13 knockdown of mTOR mRNAs (C); Western blot showing the decreased expression of 14 p-Akt (Ser473), mTOR and p-mTOR (Ser2448), p-rpS6 (Ser235/236) in treated 15 16 newborn ovaries (D). The expressions of Akt, rpS6, and β -actin were used as internal controls. siNC, control siRNA; siNC+E, control siRNA+HucMSC-exos; siM, mTOR 17 siRNAs; siM+E, mTOR siRNAs+HucMSC-exos. (E) TUNEL assay of apoptosis in 18 ovaries treated with or without HucMSC-exos. Paired ovaries were separated and 19 treated with or without HucMSC-exos for 24 h and transplanted under kidney 20 capsules of the same recipient mice. Ovaries were collected 48 h after transplantation. 21 Green, TUNEL positive cells; Blue, nuclear staining with Hoechst 33342. Bar=20 22 μm. (F) Ovarian PCNA staining in nuclei of oocytes and granulosa cells. Bar=50 μm. 23 Figure S2. Histology of follicular development between paired ovaries treated 24 with or without HucMSC-exos. Mice at 10 month of age were intra-bursa injected 25 with 20 µL HucMSC-exos in one lateral ovary and same volume of PBS in the other


lateral ovary. Paired ovaries were collected 3 weeks later for morphological analysis. 27 Representative three pairs of ovaries were shown. In each group, right panels are 28 magnifications of left panels. Bar=50µm. 29

30


26

Α






С




D


В






Ε







Supplementary Table SI				
Antibody	Company	Catalog Number	Dilution	
PE-conjugatedCD11b Antibody	BD	561689	1 test/10^6 cells	
FITC-conjugatedCD19 Antibody	BD	560994	1 test/10^6 cells	
FITC-conjugatedCD34 Antibody	BD	560942	1 test/10^6 cells	
FITC-conjugatedCD44 Antibody	BD	560977	1 test/10^6 cells	
FITC-conjugatedCD45 Antibody	BD	560976	1 test/10^6 cells	
PE-conjugatedCD73 Antibody	BD	561014	1 test/10^6 cells	
PE-conjugatedCD90 Antibody	BD	561970	1 test/10^6 cells	
APC-conjugatedCD105 Antibody	BD	562408	1 test/10^6 cells	
FITC-conjugatedHLA-DR/DQ Antibody	BD	555563	1 test/10^6 cells	
β-Actin Antibody	CST	3700	WB 1:2000	
GAPDH Antibody	CST	5174	WB 1:2000	
Alix Antibody	Proteintech	12422-1-AP	WB 1:1000	
Tsg101 Antibody	Abclonal	A2216	WB 1:1000	
CD9 Antibody	Proteintech	20597-1-AP	WB 1:1000	
Gm130 Antibody	BD	610822	WB 1:1000	
rpS6 Antibody	CST	2217	WB 1:1000	
p-rpS6 (S235/236) Antibody	CST	2211	WB 1:1000	
mTOR Antibody	CST	2983	WB 1:1000	
p-mTOR Antibody	CST	5536	WB 1:1000	
AKT Antibody	CST	2920	WB 1:1000	
p-AKT Antibody	CST	4060	WB 1:1000	
Foxo3a Antibody	CST	12829	IHC 1:200	
PCNA Antibody	CST	13110	IHC 1:1000	
Cx37 Antibody	Abcam	Ab181701	IF 1:200 WB 1:1000	

Supplementary Table S1

Gene	Forward primer	Reverse primer	
Actin	5'-CCGTAAAGACCTCTATGCC-3'	5'-CTCAGTAACAGTCCGCCTA-3'	
Gdf9	5'-TCTTAGTAGCCTTAGCTCTCAGG-3'	5'-TGTCAGTCCCATCTACAGGCA-3'	
Zp3	5'-CATCTCAAAGTCGCGCCAG-3'	5'-GCCTGCGGTTTCGAGAAAC-3'	
Bmp15	5'-TCCTTGCTGACGACCCTACAT-3'	5'-TACCTCAGGGGGATAGCCTTGG-3'	
Amhr	5'-GCAGCACAAGTATCCCCAAAC-3'	5'-GTCTCGGCATCCTTGCATCTC-3'	
Kit	5'-CTCCCCCAACAGTGTATTCAC-3'	5'-TAGCCCGAAATCGCAAATCTT-3'	
Kitl	5'-GAATCTCCGAAGAGGCCAGAA-3'	5'-GCTGCAACAGGGGGGTAACAT-3'	
Fshr	5'-CCTTGCTCCTGGTCTCCTTG-3'	5'-CTCGGTCACCTTGCTATCTTG-3'	
Star	5'-ATGTTCCTCGCTACGTTCAAG-3'	5'-CCCAGTGCTCTCCAGTTGAG-3'	

Supplementary Table S2 Real-time PCR Primers