Raman spectroscopic signatures of carotenoids and polyenes enable label-free visualization of microbial distributions within pink biofilms

Hiroto Horiue¹, Mai Sasaki¹, Yuki Yoshikawa¹, Masanori Toyofuku^{2,3} & Shinsuke Shigeto¹

¹Department of Chemistry, Graduate School of Science and Technology, Kwansei Gakuin University, 2-1 Gakuen, Sanda, Hyogo 669-1337, Japan.

²Faculty of Life and Environmental Sciences, University of Tsukuba, Ibaraki 305-8572, Japan.
³Microbiology Research Center for Sustainability, University of Tsukuba, Ibaraki 305-8572, Japan.

Correspondence and requests for materials should be addressed to S.S. (email: shigeto@kwansei.ac.jp).

Table of contents

Supplementary Figure S1	 S-2
Supplementary Figure S2	 S-2
Supplementary Figure S3	 S-3
Supplementary Figure S4	 S-3



Supplementary Figure S1. Plot of singular values obtained from SVD of the Raman hyperspectral imaging data of the pink biofilm.



Supplementary Figure S2. MCR-ALS results assuming four components (k = 4). (A) Intrinsic Raman spectra of the four components (denoted 1'-4'). Each spectrum has been normalized so that the sum of intensities at all Raman shifts is equal to unity. The spectra are vertically offset for clarity of display. (**B**-**E**) Raman images of components 1' (**B**), 2' (**C**), 3' (**D**), and 4' (**E**). The Raman images are shown in pseudocolor.



Supplementary Figure S3. Molecular structure of flexirubin.



Supplementary Figure S4. Schematic diagram of the laboratory-built, 532 nm-excited confocal Raman microspectrometer. Abbreviations: ND, neutral density; HWP, half-wave plate; VBG, volume Bragg grating; BS, beam splitter; DM, dichroic mirror; CCD, charge-coupled device.