

Supplementary Figure S1:

In the study presented the inflammatory impact of CD4⁺CD127^{dim/-} enriched PBMC, CD4⁺CD25⁺CD127^{dim/-} regulatory T cells (Treg) and a population of Treg depleted CD4⁺CD25⁻CD127^{dim/-} enriched PBMC was analyzed in a chondrocyte co-culture experiment. In preliminary studies we confirmed successful Treg isolation process using the CD4⁺CD25⁺CD127^{dim/-} Regulatory T Cell Isolation Kit II human (Miltenyi Biotec, Germany) by flow cytometry.

A representative dot plot confirming successful Treg isolation (~95%) is shown here. A/B: Autofluorescence, C: FITC CD4 PE CD25 staining, D: PE CD25 PerCP-Cy5 CD127 staining