



S4 Fig. MES₅₅₀₀ increases protein level of IL-10, but does not have effect on PI3K downstream molecules. (A-B) Jurkat T cells were treated with MES₅₅₀₀ at 2V/cm for 10 min. (A) After treatment, the medium was changed and cells were stimulated with PMA/Io for 6 hr at 37 °C, and PathScan analysis was performed according to manufacturer's protocol (*Cell Signaling Technology*). (B) Western blotting was performed on whole cells lysates from cells stimulated with PMA/Io for 3 hr at 37 °C.