

Supplementary Table 1. Fractions obtained by column chromatography of hexanic extracts from each system.

System	Fraction	acetate:hexane	Volume (mL)
T1	1-12	90:10	100
	13-24	80:20	100
	25-36	70:30	100
	37-59	60:40	200
	60-72	50:50	100
T2	1-12	90:10	100
	13-24	80:20	100
	25-36	70:30	100
	37-59	60:40	200
	60-72	50:50	100
T3	1-12	90:10	100
	13-24	80:20	100
	25-36	70:30	100
	37-59	60:40	200
	60-72	50:50	100