Matter, Volume 3

Supplemental Information

Quantitative Method for Comparative Assessment

of Particle Removal Efficiency of Fabric Masks

as Alternatives to Standard Surgical Masks for PPE

Amy V. Mueller, Matthew J. Eden, Jessica M. Oakes, Chiara Bellini, and Loretta A. Fernandez

SUPPLEMENTAL INFORMATION

Table SI1. Mask details, mean filtration efficiency (\overline{x}), standard deviation of mean filtration efficiency between replicates* (s_r), and standard deviation of filtration efficiency over one minute runs (s_t).								
		worn as designed			worn with overlayer			
Sample	Description	x	S _r	St	x	Sr	S _t	
N95-1		99.2%	0.4%	0.8%	_	-	_	
N95-2	Makrite model 9500-N95	90.6%	5.9%	4.6%	95.2%	0.9%	4.7%	
S-1	3M surgical mask model 1826	74.6%	4.1%	9.5%	90.3%	1.5%	3.9%	
S-2	Keystone surgical mask model FM-EL-BLUE	59.3%	3.3%	13.0%	86.0%	3.2%	2.7%	
S-3	Hong Da Wei Cai surgical mask labeled for medical use	53.4%	4.4%	12.6%	90.0%	6.0%	5.6%	
0-1	surgical style 4 layer mask with black "charcoal" layer (no brand information available)	73.4%	4.1%	9.7%	86.8%	0.4%	5.2%	
CS-1	cloth surgical-style mask with earloops and wired nose bridge, layers (3): two cotton quilting fabric and one Pellon interfacing fabric	58.6%	5.0%	11.6%	77.5%	6.2%	0.8%	
CS-2	fabric surgical style mask with earloops, no wire at bridge of nose, layers: two cotton plain weave	28.2%	5.9%	24.3%	73.2%	1.4%	1.2%	
CS-3	fabric surgical style mask with ties, wired nose bridge, layers (6): two Smartfab nonwoven fabric, two disposable baby wipe (dry), one massage table non- woven fabric cover, one meltblown filter (BFE85)	85.0%	1.3%	5.4%	81.3%	3.4%	7.9%	
CS-4	fabric surgical style mask with ties, wired nose bridge, layers (2): two cotton duck	72.9%	8.8%	7.1%	78.5%	12.3%	6.7%	
CS-5	fabric surgical style mask with ties, no wire at bridge of nose, layers (2): two layers of cotton twill (sold by Reformation clothing company at thereformation.com)	56.0%	3.9%	13.1%	66.9%	1.7%	10.2%	

efficiency between replicates* (s_r), and standard deviation of filtration efficiency over one minute runs (s_t).								
		worn as designed			worn with overlayer			
Sample	Description	\overline{X}	Sr	St	\overline{X}	Sr	s _t	
CS-6	fabric surgical style mask with earloops, no wire at bridge of nose, layers (2): woven nylon	47.1%	2.3%	12.2%	56.8%	5.9%	8.7%	
CC-1	commercially produced nuisance dust mask modified with cloth liner, layers (4): two Smartfab nonwoven fabric, one disposable baby wipe (dry), one meltblown filter (BFE84)	85.9%	6.3%	4.7%	89.3%	1.5%	3.8%	
CC-2	commercially produced nuisance dust mask	60.3%	3.2%	10.4%	61.1%	2.8%	9.4%	
CC-3	fabric cone-shaped mask with elastic head band and wired nose bridge, layers (6): two cotton muslin fabric, two disposable baby wipe (dry), one massage table cover non-woven fabric, one meltblown filter (BFE85)	86.2%	1.0%	5.5%	88.5%	0.9%	3.8%	
CC-4	fabric cone-shaped mask with elastic head band, layers (6): two Smartfab nonwoven fabric, two disposable baby wipe (dry), one massage table non-woven fabric cover, one meltblown filter (BFE85)	89.1%	1.7%	3.4%	91.7%	2.8%	4.3%	
CC-5	fabric cone-shaped mask with elastic head band, wired nose bridge, PM2.5 filter insert, layers (4, including pocket): three cotton muslin, one massage table non-woven fabric cover	80.2%	2.5%	7.1%	84.3%	2.5%	5.9%	

Table SI1 (cont.). Mask details, mean filtration efficiency (\overline{x}), standard deviation of mean filtration

Table SI1 (cont.). Mask details, mean filtration efficiency (\overline{x}), standard deviation of mean filtration efficiency between replicates* (s_r), and standard deviation of filtration efficiency over one minute runs (s_t).

		worn as designed			worn with overlayer		
Sample	Description	x	Sr	s _t	x	Sr	s _t
CC-6	fabric cone-shaped mask with elastic head band, layers (5): two Smartfab nonwoven fabric, one massage table non-woven fabric cover, two meltblown filter (BFE85)	90.7%	0.8%	3.1%	91.5%	1.1%	3.1%
CC-7	fabric cone-shaped mask with elastic head band, wired nose bridge, layers (4): two Smartfab nonwoven, one massage table non-woven fabric cover, two meltblown filter (BFE85)	85.3%	2.2%	4.6%	87.2%	0.9%	4.4%
CC-8	fabric cone-shaped mask with two sets of ties, wired nose bridge, layers (3): two cotton fabric, one non-woven polypropylene (recycled grocery bag)	82.6%	1.2%	5.7%	81.3%	2.4%	8.5%
CD-1	duck-bill shaped mask with elastic head band, wired nose bridge, layers (6): 4 cotton fabric, 2 Pellon interfacing	64.2%	11.0%	9.5%	80.2%	1.8%	6.3%
N only	woven nylon stocking	7.0%	2.5%	18.0%	-	-	-

* n=4 replicates for mask CS-1, all other masks n=3 replicates



Figure S1. Two TSI PortaCount model 8028 used in this work. Sample tubes are of equal length and are connected to right-hand ports labeled "sample". Instruments were operated in count mode with "Mask"-labeled instrument sampling air from inside the mask and "Ref"-labeled instrument sampling ambient air just outside of the mask.



N95-1







S-1

S-2







CS-3



CS-4



CS-6

CC-1





Figure S2. Gallery of mask images. Masks ordered by sample ID. Descriptions included in Table S1.

Additional and updated results are available through a web portal at masktestingatNU.com.