

Suppl. Table 1. Oligonucleotide sequences of forward primer, probe (in case it was used) and reverse primer for the studied genes

Gene symbol	Method	Assay ID. TaqMan, Life Technologies	Forward Primer	Probe	Reverse Primer
PTSG2	TaqMan	Ss03394694_mL			
IL17A	TaqMan	Ss03391803_mL			
TGF β	TaqMan	Ss03382325_uL			
IL6	TaqMan		5'-tcctcgcaaaatctctgcaat	FAM-5'-tttcacacatctccttctc	5'-tccttgctgtttcacacttctcat
IL8	TaqMan		5'-tcctgctttctgcagctctct	FAM-5'-tgaggctgcagttct	5'-catcgaagtctgcacttactcttg
IL10	TaqMan		5'-gaggaggtgaagagtgcttta	FAM-5'-cctctcttgagcttgc	5'-ctcacatggctttagaca
TNF α	TaqMan		5'-aacctctggccaagga	FAM-5'-tcagatcatcgtctcaaac	5'-ggcgacggcttatctga
NF κ β	TaqMan		5'-ctagtgaaccgaaacctttctactat	FAM-5'-ctgaaatcaaagataaagagg	5'-gcttctgccgttctcttgt
β -actin	TaqMan		5'-tccagaggcgctctcca	VIC-5'-tcctggcatggagt	5'-cgcacttcatgatcgagttga
GAPDH	TaqMan		5'- gtcggagtgaacggattgg	5'-cgctggtcaccagggtgct	5'-caatgtccactttgccagagttaa