



Strains

A: UMB0248	H: UMB3444
B: UMB0567	I: UMB4683
C: UMB0595	J: UMB4708
D: UMB0619	K: UMB7052
E: UMB0622	L: UMB8390
F: UMB0633	M: UMB8616
G: UMB2128	N: UMB8710

Supplemental Figure 1. Agarose gels for Trial 2 PCR-based detection of predicted prophage sequences in viral lysate (yellow) and bacterial colony controls (light blue). Negative controls (no DNA) are indicated by “-”.

Supplemental Table 3. Primers for detecting lytic phages.

Phage Group #	Primer Pair	Expected Amplicon Size (bp)
2	CTGCAACCTCATCATTGC & CCGGCGCTGTCTTATATC	475
4	ATTCGCGCTAAGAAGTGC & TGCTCAGAATGTGCTTGG	349
3	ATTGCGCAAGGACAGC & CGAATTGGTGCGACTATG	460
14	TCGCTCAATCATCTCATCC & GATATGCCGGTCTTGGAG	447
0	TGAGCTTGCGTAGGTCAG & AGCGCAGACTCAGAGAGG	307