

Supplementary Table 2. Photopic electroretinography protocol.

Step	Trials	Inter sweep delay	Sweep length	Delay to next step	Description	cd.s/m²
1	25	0	500	12,500	Single 0.5Hz	0.3
2	25	0	500	12,500	Single 0.5Hz	1
3	25	0	500	12,500	Single 0.5Hz	3
4	25	0	500	12,500	Single 0.5Hz	10
5	25	0	500	12,500	Single 0.5Hz	25
6	25	0	500	12,500	Continuous 20Hz	10
7	25	0	500	12,500	Continuous 30Hz	10