

Table S1. Primers used for analysis.

	Forward Primers
F508del allele	5'-GGCACCATTAAAGAAAATATCATTGG-3'
"Other" allele (normal F508 codon)	5'-GGCACCATTAAAGAAAATATCATCTT-3'
	Reverse Primers
"Nml" (exon 11-12 junction)	5'-CAAACCTGGAGATGTCCTCTTC-3'
" +6nt" (retention of 6 intronic nt.)	5'-CAAACCTGGAGATGTCCCATTA-3'
"Dex12" (exon 11-13 junction)	5'-CAGCATCTTTGTATACTGCTCTCTTC-3'