

Analytical and Bioanalytical Chemistry

Electronic Supplementary Material

**Interlaboratory study of an optimised peptide mapping workflow
using automated trypsin digestion for monitoring monoclonal antibody
product quality attributes**

Silvia Millán-Martín, Craig Jakes, Sara Carillo, Tom Buchanan, Marc Guender,
Dan Bach Kristensen, Trine Meiborg Sloth, Martin Ørgaard, Ken Cook, Jonathan Bones

Additional files available under [10.1007/s00216-020-02809-z](https://doi.org/10.1007/s00216-020-02809-z)

Table S1	Analytical Instrumentation
Table S2	96-well plate scheme for trypsin digestion
Table S3	BPF software parameter settings for peptide mapping data analysis
Table S4	CDS software parameter settings for peptide mapping data processing and PTMs quantitation.
Table S5	MS Component table of CDS software for the studied PTMs (see separate .xlsx file)
Figure S1	UV Chromatogram traces (280nm) for trastuzumab RP-UV-MS Intact analysis during digestion time course study using buffer 1
Figure S2	Base Peak Chromatograms (BPC) for trastuzumab RP-UV-MS Intact Analysis using 5mM TCEP as reducing agent (a), and without TCEP addition (b)
Figure S3	Sequence coverage map of trastuzumab obtained using automated tryptic digestion with low pH buffer 1 after 5min.
Figure S4	Sequence coverage map of trastuzumab obtained using automated tryptic digestion with SMART buffer 2 after 5min
Table S6	Sequence coverage for the time course of digestion for trastuzumab
Table S7	Summary of PTMs identified and quantified for trastuzumab after digestion with trypsin beads, low pH buffer and 5mM TCEP addition on the KingFisher Duo Prime system
Table S8	Summary of PTMs identified and quantified for trastuzumab after digestion with trypsin beads, SMART buffer and 5mM TCEP addition on the KingFisher Duo Prime system
Table S9	Component list from BPF3.1 software for sequence coverage and Venn diagram (see separate .xlsx file)
Table S10	Analysis of Variance (ANOVA) to measure intra and inter-lab precision for NISTmAb complete compliant peptide mapping workflow

Table S1 Analytical Instrumentation used for the 4 Laboratories across Europe

	Site A (Ireland)	Site B (Denmark)	Site C (UK)	Site D (Switzerland)
UHPLC	Vanquish™ Flex Binary	Vanquish™ Duo Ternary	Vanquish™ Horizon Binary	Vanquish™ Flex Binary
HRMS	Q Exactive™ Plus	Q Exactive™ Plus	Q Exactive™ Plus	Q Exactive™ Plus
SMART Digest	Magnetic Beads KingFisher™ Duo	Magnetic Beads KingFisher™ Duo	Magnetic Beads KingFisher™ Duo	Manual Magnetic beads

Table S2 96-well plate scheme for trypsin digestion

Lane	Content	Volume (µL)
A	Low pH buffer	150
	Sample (100µg)	10
	TCEP (final concentration 5mM)	2
	Water	38
B	Tip Comb	
C	Empty	
D	Magnetic SMART™ Beads	15
	Bead Buffer (Low pH buffer)	100
E	Bead Wash Buffer (Bead Buffer 1:4 (v/v))	200
F	Waste Lane (Water)	250

Table S3 BPF software parameter settings for peptide mapping data analysis

Component Detection	Setting
Absolute MS Signal Threshold	1.0 x 10 ⁴ counts
Typical chromatographic peak width	0.3
Relative MS signal threshold (% base peak)	1
Relative Analog threshold (% of highest peak)	1
Width of Gaussian filter (represented as 1/n of chromatographic peak width)	3
Minimum valley to be considered as two chromatographic peaks	80.0%
Minimum MS peak width (Da)	1.2
Maximum MS peak width (Da)	4.2
Mass tolerance (ppm for high-res or Da for low-res)	4.00
Maximum retention time shift (min)	1.00
Maximum mass (Da)	30,000
Mass Centroiding Cutoff (% from base)	15
Identification	Setting
Maximum peptide mass	7,000
Mass Accuracy	5ppm
Minimum Confidence	0.8
Maximum Number of Modifications for a Peptide	1
Unspecified Modification	-58 to +162 Da
N-Glycosylation	CHO
Protease Specificity	High
Static Modifications	Setting
Side Chain	-
Variable Modifications	Setting
N Terminal	Gln - > Pyro Glu
C Terminal	Loss of lysine
Side Chain	Deamidation (N) Deamidation(Q) Glycation (K) Oxidation (MW) AsnAsu (N) AspAsu (D)

Table S4 Chromeleon™ Data System 7.2.9 parameter settings for peptide mapping data processing and PTMs quantitation

MS Detection	Setting
Extracted Ion Chromatogram	MS Default Detection Settings
Detection Algorithm	ICIS
Area noise factor	5
Peak noise factor	10
Baseline window	40
Noise method	INCOS
Min peak width	3
Multiplet resolution	10
Area tail extension	5
Area scan window	0
MS Component table	Setting
Supplementary Table S5 (.xlsx file)	.wbpf workbook from BPF3.1
MS Settings	MS Chromatogram Settings
Mass Precision	5 decimal places
Mass Tolerance (manually defined)	8.0 ppm
Smoothing	none
Peptide Table	Setting
Imported from BPF3.1 peptide workbook	-
Composite Scoring	Setting
Pass score if at least	2 criteria passed
Fail score if less than	1 criteria passed
MS Criteria	General MS
Isotopic dot product	≥ 0.9000
Mass accuracy	≤ 5.00 ppm
Peal apex alignment	≤ 0.50 min




Fig. S1 UV Chromatogram traces (280nm) for trastuzumab RP-UV-MS Intact analysis during digestion time course study using buffer 1




Fig. S2 Base Peak Chromatograms (BPC) for trastuzumab RP-UV-MS Intact Analysis using 5mM TCEP as reducing agent (a), and without TCEP addition (b). Digestion was performed for 35min and buffer 2




Fig. S3 Sequence coverage map of trastuzumab heavy (upper panel) and light (lower panel) chains, obtained using automated tryptic digestion with buffer 1 after 5 minutes. The coloured bars show the identified peptides, with the numbers in the bars reflecting the retention time. The different colors indicate the peptide recovery in the MS1 scan: red >50%, orange >20% and yellow >10% represent good recovery. Green, >5%, light blue >2% and cyan >1% represent fair recovery and grey-white scale shows poor recovered peptides




Fig. S4 Sequence coverage map of trastuzumab heavy (upper panel) and light (lower panel) chains, obtained using automated tryptic digestion with buffer 2 after 5 minutes. The coloured bars show the identified peptides, with the numbers in the bars reflecting the retention time. The different colors indicate the peptide recovery in the MS1 scan: red >50%, orange >20% and yellow >10% represent good recovery. Green, >5%, light blue >2% and cyan >1% represent fair recovery and grey-white scale shows poor recovered peptides

Table S6 Sequence coverage for the time course of digestion for trastuzumab

Proteins		Avg Number of MS Peaks (n=3)	Average Sequence Coverage (n=3)
Buffer 1, 5min	HC	353	100%
	LC	103	100%
Buffer 1, 10min	HC	449	100%
	LC	134	100%
Buffer 1, 15min	HC	510	100%
	LC	152	100%
Buffer 1, 20min	HC	578	100%
	LC	167	100%
Buffer 1, 25min	HC	570	100%
	LC	161	100%
Buffer 1, 30min	HC	553	100%
	LC	158	100%
Buffer 1, 35min	HC	536	100%
	LC	147	100%
Buffer 1, 40min	HC	594	100%
	LC	180	100%
Buffer 2, 5min	HC	353	100%
	LC	95	100%
Buffer 2, 10min	HC	548	100%
	LC	156	100%
Buffer 2, 15min	HC	593	100%
	LC	169	100%
Buffer 2, 20min	HC	617	100%
	LC	172	100%
Buffer 2, 25min	HC	621	100%
	LC	165	100%
Buffer 2, 30min	HC	599	100%
	LC	155	100%
Buffer 2, 35min	HC	555	100%
	LC	149	100%
Buffer 2, 40min	HC	649	100%
	LC	172	100%

Table S7 Summary of PTMs identified and quantified for trastuzumab after digestion with trypsin beads, buffer 1 and 5mM TCEP addition on the KingFisher™ Duo Prime system

Modification	5min (n=3)	10min (n=3)	15min (n=3)	20min (n=3)	25min (n=3)	30min (n=3)	35min (n=3)	40min (n=3)
HC N55+Deamidation	—	0.63	0.61	0.98	1.17	1.78	2.23	3.64
HC N77+Deamidation	0.07	0.04	0.06	0.12	0.16	0.24	0.25	0.26
HC N84+Deamidation	0.02	0.02	0.04	0.07	0.07	0.12	0.12	0.12
HC N289+Deamidation	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
HC N318+Deamidation	0.05	0.11	0.05	0.22	0.16	0.28	0.90	0.26
HC N364+Deamidation	0.01	0.04	0.02	0.12	0.09	0.09	0.33	0.44
LC N30+Deamidation	6.82	8.41	8.24	8.89	8.94	9.23	9.80	9.12
LC ~N137+Deamidation	0.88	0.30	0.21	0.72	1.26	3.73	1.30	1.53
LC N210+Deamidation	0.00	0.00	0.00	0.00	0.01	0.02	0.02	0.07
HC M255+Oxidation	1.36	1.06	1.11	1.02	1.13	1.24	1.10	1.25
HC M431+Oxidation	0.07	0.19	0.19	0.18	0.21	0.14	0.22	0.15
HC G449+Lys	1.08	1.44	1.60	1.62	1.66	1.60	1.83	1.69
HC N55+AsnAsu	1.00	1.21	1.45	2.14	2.36	2.93	3.39	3.60
HC N318+AsnAsu	0.48	1.05	1.74	2.77	3.38	3.93	3.41	3.90
LC N30+AsnAsu	0.19	1.48	1.53	1.08	0.87	0.50	0.50	0.31
LC ~N138+AsnAsu	9.67	1.49	1.71	3.06	3.87	7.43	2.04	1.92

Table S7 (continued) Summary of PTMs identified and quantified for trastuzumab after digestion with trypsin beads, buffer 1 and 5mM TCEP addition on the KingFisher™ Duo Prime system

Modification	5min (n=3)	10min (n=3)	15min (n=3)	20min (n=3)	25min (n=3)	30min (n=3)	35min (n=3)	40min (n=3)
HC ~D102+AspAsu	0.21	0.48	0.73	0.92	1.10	1.52	1.27	1.47
HC ~D283+AspAsu	0.53	1.49	2.57	2.68	2.57	4.02	3.64	4.54
HC ~D404+AspAsu	0.28	0.03	0.01	1.43	1.64	1.56	1.29	1.033
HC D224+isomerization	69.11	55.64	70.37	52.95	66.25	66.19	75.22	84.46
HC N300+A1G0	2.07	2.35	2.16	2.15	2.10	2.10	2.00	2.16
HC N300+A1G0F	11.43	8.76	7.43	7.28	6.40	5.48	7.27	6.32
HC N300+A1G1F	3.26	1.99	1.75	1.76	1.45	1.19	1.64	1.78
HC N300+A2G0	3.78	5.16	5.46	5.56	5.69	5.86	5.62	5.82
HC N300+A2G0F	48.99	50.91	52.15	52.18	53.23	54.56	52.07	51.93
HC N300+A2G1F	22.04	22.69	23.41	23.44	23.50	23.19	23.86	23.27
HC N300+A2G2F	2.00	1.84	1.91	1.75	1.80	1.79	2.16	2.17
HC N300+M5	2.17	2.11	2.08	2.15	2.16	2.29	2.05	2.26
HC N300+Unglycosylated	4.26	4.18	3.63	3.73	3.65	3.53	3.32	4.28

Table S8 Summary of PTMs identified and quantified for trastuzumab after digestion with trypsin beads, buffer 2 and 5mM TCEP addition on the KingFisher™ Duo Prime system

Modification	5min (n=3)	10min (n=3)	15min (n=3)	20min (n=3)	25min (n=3)	30min (n=3)	35min (n=3)	40min (n=3)
HC N55+Deamidation	0.64	0.48	0.85	1.04	1.69	1.95	3.14	3.58
HC N77+Deamidation	0.01	0.02	0.08	0.13	0.15	0.19	0.29	0.36
HC N84+Deamidation	0.14	0.02	0.05	0.07	0.09	0.09	0.08	0.15
HC N289+Deamidation	0.01	0.02	0.04	0.05	0.14	0.075	0.04	0.09
HC N318+Deamidation	0.05	0.43	0.28	1.17	0.62	1.43	1.023	3.26
HC N364+Deamidation	0.00	0.01	0.01	0.13	0.21	0.13	0.42	0.47
LC N30+Deamidation	6.92	8.90	8.88	9.65	9.34	9.58	10.06	9.40
LC ~N137+Deamidation	0.00	0.27	0.62	1.11	2.10	5.00	1.33	2.60
LC N210+Deamidation	0.00	0.00	0.00	0.01	0.02	0.03	0.04	0.07
HC M255+Oxidation	1.18	0.48	0.47	0.43	0.48	0.85	0.29	0.26
HC M431+Oxidation	0.09	0.17	0.17	0.15	0.15	0.14	0.17	0.17
HC G449+Lys	0.57	1.44	1.53	1.63	1.63	1.56	1.81	1.76
HC N55+AsnAsu	1.30	1.22	1.89	2.29	2.73	3.24	3.57	3.52
HC N318+AsnAsu	0.51	1.52	3.03	6.31	5.49	9.19	3.90	4.44
LC N30+AsnAsu	0.31	1.32	0.89	0.73	0.48	0.39	0.27	0.21
LC ~N137+AsnAsu	0.00	0.00	0.02	0.02	0.00	0.00	0.00	0.00

Table S8 (continued) Summary of PTMs identified and quantified for trastuzumab after digestion with trypsin beads, buffer 2 and 5mM TCEP addition on the KingFisher™ Duo Prime system.

Modification	5min (n=3)	10min (n=3)	15min (n=3)	20min (n=3)	25min (n=3)	30min (n=3)	35min (n=3)	40min (n=3)
HC ~D102+AspAsu	0.37	0.16	0.80	0.81	0.96	1.38	1.14	1.11
HC ~D283+AspAsu	0.90	0.46	4.02	6.45	5.28	6.92	3.57	4.54
HC ~D404+AspAsu	0.03	1.32	2.02	2.21	1.13	1.84	0.68	1.90
HC D224+isomerization	0.03	52.43	45.51	45.68	64.12	54.92	49.51	47.36
HC N300+A1G0	1.98	2.35	2.15	2.07	2.10	2.10	2.11	2.04
HC N300+A1G0F	7.54	10.03	7.75	7.12	6.38	5.43	7.56	6.73
HC N300+A1G1F	1.85	2.35	1.51	1.72	1.34	1.16	1.52	1.17
HC N300+A2G0	4.48	5.03	5.29	5.42	5.64	5.93	5.29	5.51
HC N300+A2G0F	53.65	50.11	52.68	52.89	53.19	54.71	52.33	52.80
HC N300+A2G1F	23.28	22.74	23.59	23.64	23.99	23.17	23.82	24.43
HC N300+A2G2F	1.77	2.05	2.15	2.18	2.17	1.98	2.22	2.27
HC N300+M5	1.99	2.01	2.03	2.06	2.16	2.26	2.02	2.14
HC N300+Unglycosylated	3.45	3.33	2.84	2.89	3.02	3.26	3.12	2.90

Table S10 Analysis of Variance (ANOVA) to measure intra and inter-lab precision for NISTmAb complete compliant peptide mapping workflow

HC N328+Deam						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.1001	3	0.0334	153.9872	2.05E-07	4.0662
Within Groups	0.0017	8	0.0002			
Intra-Lab precision (RSD)	10.22					
Inter-Lab precision (RSD)	44.43					

HC N364+Deam						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.1659	3	0.0553	49.1556	1.69E-05	4.0662
Within Groups	0.0090	8	0.0011			
Intra-Lab precision (RSD)	42.16					
Inter-Lab precision (RSD)	101.64					

HC ~N392/N387+ Deam						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.7405	3	0.2468	219.4049	5.09E-08	4.0662
Within Groups	0.0090	8	0.0011			
Intra-Lab precision (RSD)	20.22					
Inter-Lab precision (RSD)	105.25					

HC D283+Succ						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	8.3645	3	2.7882	420.8549	3.86E-09	4.0662
Within Groups	0.0530	8	0.0066			
Intra-Lab precision (RSD)	8.61					
Inter-Lab precision (RSD)	62.24					

HC N318+Succ						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.7734	3	0.2578	52.6100	1.31E-05	4.0662
Within Groups	0.0392	8	0.0049			
Intra-Lab precision (RSD)	8.28					
Inter-Lab precision (RSD)	20.68					

HC ~N392/N387+Succ						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.8415	3	0.2805	125.5970	4.56E-07	4.0662
Within Groups	0.0179	8	0.0022			

Intra-Lab precision (RSD)	5.54
Inter-Lab precision (RSD)	21.73

Table S10 (continued) Analysis of Variance (ANOVA) to measure intra and inter-lab precision for NISTmAb complete compliant peptide mapping workflow

HC M255+Oxid

<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.2813	3	0.0938	468.8194	2.51E-09	4.0662
Within Groups	0.0016	8	0.0002			
Intra-Lab precision (RSD)	3.54					
Inter-Lab precision (RSD)	27.04					

HC K450 Lys Loss

<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	21.0218	3	7.0073	225.5554	4.57E-08	4.0662
Within Groups	0.2485	8	0.0311			
Intra-Lab precision (RSD)	0.56					
Inter-Lab precision (RSD)	2.96					

HC Q1+Gln->Pyro-Glu

<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.4365	3	0.1455	831.3651	2.57E-10	4.0662
Within Groups	0.0014	8	0.0002			
Intra-Lab precision (RSD)	0.04					
Inter-Lab precision (RSD)	0.38					

HC N300+M5

<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.6293	3	0.2098	143.0152	2.74E-07	4.0662
Within Groups	0.0117	8	0.0015			
Intra-Lab precision (RSD)	8.14					
Inter-Lab precision (RSD)	34.11					

Table S10 (continued) Analysis of Variance (ANOVA) to measure intra and inter-lab precision for NISTmAb complete compliant peptide mapping workflow

HC N300+A1G0F						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	94.4941	3	31.4980	438.9969	3.26E-09	4.0662
Within Groups	0.5740	8	0.0718			
Intra-Lab precision (RSD)	9.10					
Inter-Lab precision (RSD)	67.20					

HC N300+A2G0F						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	36.8955	3	12.2985	55.2535	1.09E-05	4.0662
Within Groups	1.7807	8	0.2226			
Intra-Lab precision (RSD)	3.43					
Inter-Lab precision (RSD)	8.79					

HC N300+A1G1F						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	14.1954	3	4.7318	183.6998	1.03E-07	4.0662
Within Groups	0.2061	8	0.0258			
Intra-Lab precision (RSD)	9.18					
Inter-Lab precision (RSD)	43.65					

HC N300+A2G1F						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	38.3946	3	12.7982	196.0907	7.93E-08	4.0662
Within Groups	0.5221	8	0.0653			
Intra-Lab precision (RSD)	2.01					
Inter-Lab precision (RSD)	9.90					

HC N300+A2G2F						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	4.7004	3	1.5668	104.6852	9.28E-07	4.0662
Within Groups	0.1197	8	0.0150			
Intra-Lab precision (RSD)	4.65					
Inter-Lab precision (RSD)	16.62					

HC N300+ A2Ga1G1F						
<i>Source of Variation</i>	<i>SS</i>	<i>df</i>	<i>MS</i>	<i>F</i>	<i>P-value</i>	<i>F crit</i>
Between Groups	0.3505	3	0.1168	175.2333	1.23E-07	4.0662
Within Groups	0.0053	8	0.0007			
Intra-Lab precision (RSD)	6.24					
Inter-Lab precision (RSD)	28.99					