Supplementary Material

Simple and Fast Assay for Apolipoprotein E Phenotyping and Glycotyping: Discovering Isoform-Specific Glycosylation in Plasma and Cerebrospinal Fluid

Supplementary Table 1. Characteristics of the study population.

Variable	N=22		
ApoE, N (%)			
	5 (22.7%)		
	ε2/ε4	2 (9%)	
	ε3/ε3	5 (22.7%)	
	ε3/ε4	5 (22.7%)	
	ε4/ε4	5 (22.7%)	
Age, median (min, max)		66 (52, 81)	
Gender, N (%)			
	female	12 (54.5%)	
	male	10 (45.5%)	
Education, y, median (min,	17 (12, 20)		
Ethnicity, N (%)			
nc	18 (81.8%)		
	4 (18.2%)		
MoCA, median (min, max)	(missing = 6)	25 (20, 30)	
CDR, N (%)	(missing = 1)		
	0	16 (72.7%)	
	0.5	5 (22.7%)	

Supplementary Table 2. ApoE Assay reproducibility. Plasma samples from two individuals $(\epsilon 3/\epsilon 4$ and $\epsilon 4/\epsilon 4)$ were analyzed three times (intra-assay CVs) on three different days (inter-assay CVs).

		Intra-assay CVs			Inter-assay CVs		
Sample 1 (ε3/ε4)	E4/E3 Ratio		Day 1	Day 2	Day 3	·	
		STDEV	0.031	0.015	0.050	STDEV	0.090
		MEAN	1.38	1.41	1.56	MEAN	1.45
		CV (%)	2.25	1.06	3.21	CV (%)	6.21
	Glycosylation Level		Day 1	Day 2	Day 3		
		STDEV	0.015	0.011	0.012	STDEV	0.013
		MEAN	0.121	0.132	0.135	MEAN	0.130
		CV (%)	12.4	8.33	8.89	CV (%)	10.0
Sample 2 (ε4/ε4)	Glycosylation Level		Day 1	Day 2	Day 3		
		STDEV	0.008	0.002	0.003	STDEV	0.008
		MEAN	0.134	0.125	0.118	MEAN	0.126
		CV (%)	5 97	1 60	2 54	CV (%)	6 35

Supplementary Figure 1. Mass spectra resulting from the analysis of matched plasma and CSF from (a,b) heterozygous $\varepsilon 2/\varepsilon 3$, and (c,d) heterozygous $\varepsilon 2/\varepsilon 4$ individual. Matrix adduct peaks are labeled with *. For the glycan peak labeling, structure, and predicted masses, see Table 1.

