

S1 Table. Details of serologic test results from AFIAS COVID-19 Ab and/or EDI™ Novel Coronavirus COVID-19 IgG/IgM ELISA Kit for 17 symptomatic COVID-19 patients

COVID-19 Patients	d PSO (d)	IgM		IgG	
		AFIAS (COI)	ELISA (S/P ratio)	AFIAS (COI)	ELISA (S/P ratio)
A	13	0.13	0.89	20.94	1.74
	19	0.04	0.56	24.72	2.61
	23	0.02	0.57	23.77	2.54
	27	0.00	0.53	22.26	2.56
	42	0.00	0.52	23.74	2.42
B	8	0.00	0.47	0.00	0.35
	11	0.31	1.09	0.00	0.62
	15	0.38	1.10	23.20	2.26
	20	0.29	1.11	22.71	1.99
	24	0.40	0.89	22.90	1.94
	27	0.31	0.66	21.61	1.83
	31	0.49	0.63	14.09	1.98
	53	0.17	0.60	26.38	1.69
C	67	0.04	0.59	23.20	1.66
	5	0.00	0.61	0.00	0.35
	8	0.00	0.48	0.00	0.29
	9	0.00	0.46	0.00	0.28
	12	0.00	0.48	0.00	0.28
D	38	0.00	0.50	22.03	1.95
	16	1.15	1.46	14.82	3.02
	20	1.47	1.04	18.82	3.10
	22	1.63	1.00	20.46	3.15
E	30	2.10	0.86	23.34	3.13
	11	0.00	0.51	0.00	0.46
	12	0.02	0.59	0.00	0.44
	14	0.13	0.78	0.00	0.77
	16	0.11	0.52	0.00	0.48
	19	0.16	1.04	20.04	2.38
	24	0.06	0.92	22.72	2.73
	29	0.14	0.96	28.92	2.68
32	0.10	1.15	24.95	2.58	

(S1 Table continue)

COVID-19 Patients	d PSO (d)	IgM		IgG	
		AFIAS (COI)	ELISA (S/P ratio)	AFIAS (COI)	ELISA (S/P ratio)
F	9	0.05	0.48	0.00	0.72
	11	0.20	0.58	15.93	1.45
	15	2.21	0.81	26.03	2.92
	16	2.51	0.78	23.88	3.00
	17	2.02	0.70	22.78	2.85
	18	3.37	0.71	22.80	2.91
	23	3.27	0.68	19.48	2.95
	43	5.15	0.23	18.39	2.75
G	45	0.26	0.48	9.99	1.18
	47	0.33	1.10	23.73	2.69
	53	0.32	1.15	26.82	2.89
	60	0.28	1.11	25.30	2.91
H	8	0.07	0.41	0.00	0.57
	13	4.46	0.82	20.29	2.45
	14	3.21	0.78	21.20	2.41
	36	6.03	0.00	23.66	0.00
I	8	0.01	0.38	0.00	0.27
	13	6.59	1.50	0.00	0.76
	16	5.77	3.19	0.00	2.18
	17	5.66	2.73	0.00	2.08
	20	4.87	1.29	3.96	2.29
J	3	0.14	0.42	0.00	0.48
	4	0.17	0.44	0.00	0.46
	6	0.21	0.48	0.00	0.45
	9	0.15	0.47	0.00	0.52
	13	0.11	0.50	13.52	0.63
	42	0.11	0.17	14.14	0.82

(S1 Table continue)

COVID-19 Patients	d PSO (d)	IgM		IgG	
		AFIAS (COI)	ELISA (S/P ratio)	AFIAS (COI)	ELISA (S/P ratio)
K	4	0.29	0.50	0.00	0.36
	5	0.29	0.53	0.00	0.36
	7	0.30	0.50	0.00	0.33
	10	0.25	0.42	0.00	0.42
	14	0.15	0.50	28.11	0.89
	28	0.06	0.49	27.59	1.09
	43	0.00	0.31	31.46	1.50
L	3	1.59	0.42	0.00	0.28
	5	1.49	0.38	0.00	0.28
	12	1.29	0.00	0.00	0.00
	16	1.52	0.48	0.00	0.28
M	40	0.00	0.45	17.21	0.32
	41	0.00	0.62	23.64	0.74
	42	0.01	0.62	23.93	0.83
	44	0.09	0.63	21.59	1.07
	55	0.05	0.34	27.89	1.03
N	5	0.06	0.56	0.00	0.91
	6	0.04	0.55	0.00	1.45
	9	0.04	0.53	0.00	0.74
	12	0.37	0.58	11.88	0.63
	14	0.43	0.54	23.97	0.64
	15	0.30	0.66	22.54	1.09
	19	0.26	0.26	24.94	1.53
O	3	0.00	0.81	0.00	1.50
	4	0.00	0.89	0.00	1.94
	5	0.00	0.52	0.00	0.76
	9	0.01	0.57	0.00	0.79
	16	1.14	2.38	23.42	2.83

(S1 Table continue)

COVID-19 Patients	d PSO (d)	IgM		IgG	
		AFIAS (COI)	ELISA (S/P ratio)	AFIAS (COI)	ELISA (S/P ratio)
P	3	0.06	0.60	0.00	0.73
	4	0.04	0.62	0.00	0.89
	5	0.04	0.55	0.00	0.48
	9	0.24	0.54	13.63	0.51
	16	0.22	1.19	25.19	1.62
Q	2	0.13	0.53	9.85	0.47
	3	0.10	0.69	8.16	0.72
	4	0.10	0.51	10.20	0.69
	5	0.11	0.52	1.41	0.65
	9	0.87	0.57	26.49	0.69
	16	1.40	2.39	24.84	2.36

Abbreviation: COVID-19, coronavirus disease 2019; d PSO, days post symptom onset; AFIAS, AFIAS COVID-19 Ab; ELISA, EDITM Novel Coronavirus COVID-19 IgG/IgM ELISA Kit Ab; COI, cut-off index; S/P ratio, sample-to-positive control optical densities ratio.