



Supplemental Figure 1: Flow cytometry validation of negatively immunoselected splenocytes. Mouse splenocytes were incubated for 20 min at room temperature with anti-mouse fluorochrome-conjugated monoclonal antibodies against CD3 and CD4 (Miltenyi Biotec), fixed for 7 min with 4% paraformaldehyde and analyzed on a FACS CyAn ADP analyzer (Beckman Coulter). The upper diagram shows the distribution of double positive events. The lower diagram shows the distribution of CD3 positive events (T cells).