

Figure S1. in vitro does test of combined cell-therapy with MDSC and Treg on Treg/Th17 regulation

CD4+T cell(1 × 10⁶) isolated from normal C57BL/6 mice were cocultured with MDSC or Treg alone or combined MDSC and Treg cells for 3 days in the presence of anti-mouse CD3 and CD28 antibody and analyzed by flow cytometry. A plot from one representative experiment displays the proportions of IL-17+, CD25+Foxp3+ cells among CD4+ T cells. Numbers in the plots indicate percentages of gated cells.