

Supporting information

Antitubercular and Antiparasitic 2-Nitroimidazopyrazinones with Improved Potency and Solubility

Chee Wei Ang^{1,#}, Lendl Tan^{2,3}, Melissa L. Sykes⁴, Neda AbuGharbiyeh⁵, Anjan Debnath⁵, Janet
C. Reid¹, Nicholas P. West^{2,3}, Vicky M. Avery⁴, Matthew A. Cooper^{1,3}, Mark A.T.
Blaskovich^{1,3*}

¹ Centre for Superbug Solutions, Institute for Molecular Bioscience, The University of
Queensland, St Lucia, Queensland 4072, Australia

² School of Chemistry and Molecular Bioscience, The University of Queensland, St Lucia,
Queensland 4072, Australia

³ Australian Infectious Diseases Research Centre, St Lucia, Queensland 4067, Australia

⁴ Discovery Biology, Griffith Institute for Drug Discovery, Griffith University, Nathan,
Queensland 4111, Australia

⁵ Center for Discovery and Innovation in Parasitic Diseases, Skaggs School of Pharmacy and
Pharmaceutical Sciences, University of California, San Diego, La Jolla, California 92093 United
States

Table of Contents

LC/MS/MS detection and analysis parameters for plasma protein binding, microsomal and plasma stability	S2
Table S1. Mass spectrometer parameters of analogs tested for their plasma protein binding, microsomal and plasma stability.....	S3
HEK293 mammalian cytotoxicity assay (assay B) for <i>T.b.brucei</i> selectivity	S3
Table S2. Additional cytotoxicity data to determine selectivity indices of compounds 17a–k , 26a–k , and 34a–p	S4
¹ H and ¹³ C NMR spectra.....	S6
References.....	S82

LC/MS/MS detection and analysis parameters for plasma protein binding, microsomal and plasma stability

LC/MS/MS analysis was performed using AB Sciex 4000 QTRAP System and Shimadzu Nexara UPLC System equipped with Waters Atlantis T3, 2.1 × 50 mm, 5 μm, with guard column and column temperature was set to 40 °C. Injection volume was 3–5 μL and the autosampler cooler was set to 12 °C. The flow rate was 0.2 mL/min with gradient: 5% B for 1 minute then 5–100% B for 5 minutes (A: 0.1% (v/v) formic acid in water and B: 0.1% (v/v) formic acid in acetonitrile). Mass spectrometer parameters in Selected Reaction Monitoring (SRM) mode are tabulated below:

Table S1. Mass spectrometer parameters of analogs tested for their plasma protein binding, microsomal and plasma stability.

Compound	Ionization Mode	Declustering Potential	Collision Energy	<i>m/z</i>
pretomanid 2	positive	70	32	360.1→174.8 360.1→200.9
verapamil	positive	70	35	455.3→165.3 455.3→303.3
eucatropine	positive	70	30	292.3→109.4
carbutamide	positive	70	25	272.1→155.8
sulfamethoxazole	positive	60	17	254.3→155.9
17k	positive	70	35	329.2→105.1
26j	positive	70	35	322.1→143.2
34e	positive	70	35	348.2→168.1
34g	positive	70	40	432.0→252.2
34j	positive	70	35	432.1→252.2
34k	positive	70	40	432.2→252.2
34p	positive	70	40	354.2→174.2

HEK293 mammalian cytotoxicity assay (assay B) for *T. b. brucei* selectivity

Compounds were tested for activity against human embryonic kidney cells (HEK293), as described in Sykes et al (2012), with a slight modification, by using resazurin as a viability indicator. Briefly, 7.27×10^4 HEK cells/ml were added in 55 μ L of DMEM, high glucose medium to 384-well plates. Following 24 h incubation at 37°C / 5% CO₂, compounds were added following pre-dilution in DMEM on a Minitrak compound handling system (PerkinElmer). Cells were grown for 72 h before the addition of 10 μ L of a final concentration of 49 μ M of resazurin, prepared in DMEM. Following 4 h incubation, plates were read on an Envision plate reader (PerkinElmer).

IC₅₀ values were calculated from two independent experiments in GraphPad Prism 6 (GraphPad Software).

Table S2. Additional cytotoxicity data to determine selectivity indices of compounds 17a–k, 26a–k, and 34a–p.


Compound	HEK293 (Assay A) ^a CC ₅₀ (μM)	HEK293 (Assay B) ^b CC ₅₀ (μM)	3T3 ^c CC ₅₀ (μM)
17a	>100	>79.4	>73.3
17b	>100	>79.4	>73.3
17c	>100	>79.4	>73.3
17d	>100	>79.4	>73.3
17e	>100	>79.4	>73.3
17f	>100	>79.4	>73.3
17g	>100	>79.4	ND
17h	>100	ND	ND
17i	>100	>79.4	>73.3
17j	>100	>79.4	>73.3
17k	>100	>79.4	>73.3
26a	>100	>73.3	ND
26b	>100	>73.3	ND
26c	>100	>73.3	ND
26d	>100	>73.3	ND
26e	>100	>73.3	ND
26f	>25	>18.3	ND
26g	>118	ND	ND
26h	>94	>79.4	>36.6
26i	>100	>73.3	ND
26j	>100	>79.4	>73.3
26k	>100	>36.6	ND

Compound	HEK293 (Assay A)^a CC₅₀ (μM)	HEK293 (Assay B)^b CC₅₀ (μM)	3T3^c CC₅₀ (μM)
34a	>100	>73.3	>36.6
34b	>100	>79.4	>73.3
34c	>50	>36.6	>36.6
34d	>100	>73.3	>73.3
34e	>92	>79.4	>36.6
34f	>50	>36.6	>18.3
34g	>50	>39.7	>36.6
34h	>50	>36.6	>18.3
34i	>50	>36.6	>36.6
34j	>74	>79.4	>36.6
34k	>100	>79.4	>36.6
34l	>94	>79.4	>36.6
34m	>100	>79.4	>36.6
34n	>50	>36.6	>18.3
34o	>50	>39.7	>18.3
34p	>100	>73.3	>36.6

^a Assay A was performed based on 20 h of incubation. ^b Assay B was performed based on 72 h of incubation, following a procedure as described in S3. ^c Assay was performed based on 48 h of incubation, to determine the selectivity to *T. cruzi* amastigotes over 3T3 host cells. ND – not determined.


¹H and ¹³C NMR spectra

2-Nitro-7-(2-oxo-2-phenylethyl)imidazo[1,2-a]pyrazin-8(7H)-one (17a)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8692_034_new_ppt_ACN_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	32
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	2.2807
16 Acquisition Date	2017-10-08T13:27:40
17 Modification Date	2017-10-08T13:27:40
18 Spectrometer Frequency	600.13
19 Spectral Width	7183.9
20 Lowest Frequency	-591.3
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(2-oxo-2-phenylethyl)imidazo[1,2-a]pyrazin-8(7H)-one (17a)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_034_new_ppt_ACN_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	10000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-10-10T07:39:28
17 Modification Date	2017-10-10T07:39:28
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	13C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(2-oxo-2-(4-(trifluoromethoxy)phenyl)ethyl)imidazo[1,2-a]pyrazin-8(7H)-one (17b)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8692_088_ppt_ACN_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	155.2
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	2.2807
16 Acquisition Date	2017-10-04T12:02:56
17 Modification Date	2017-10-04T12:02:56
18 Spectrometer Frequency	600.13
19 Spectral Width	7183.9
20 Lowest Frequency	-591.3
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(2-oxo-2-(4-(trifluoromethoxy)phenyl)ethyl)imidazo[1,2-a]pyrazin-8(7H)-one (17b)




¹³C NMR (150 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8692_088_ppt_ACN_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	5800
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-10-04T17:17:55
17 Modification Date	2017-10-04T17:17:55
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(2-oxo-2-(*p*-tolyl)ethyl)imidazo[1,2-*a*]pyrazin-8(7*H*)-one (17c)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_072_ppt_ACN_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	172.8
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	2.2807
16 Acquisition Date	2017-09-28T17:11:06
17 Modification Date	2017-09-28T17:11:06
18 Spectrometer Frequency	600.13
19 Spectral Width	7183.9
20 Lowest Frequency	-591.3
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(2-oxo-2-(*p*-tolyl)ethyl)imidazo[1,2-*a*]pyrazin-8(7*H*)-one (17c)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_072_ppt_ACN_MeOH.3.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	10000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-09-28T22:32:10
17 Modification Date	2017-09-28T22:32:10
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-(2-Nitro-8-oxoimidazo[1,2-a]pyrazin-7(8H)-yl)-N-phenylacetamide (17d)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8692_039_combined_ppt _ACN_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2017-09-21T18:49:11
17 Modification Date	2017-09-21T18:49:11
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-(2-Nitro-8-oxoimidazo[1,2-a]pyrazin-7(8H)-yl)-N-phenylacetamide (17d)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_039_combined_ppt_ACN_MeOH.3.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	7500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-09-21T23:01:27
17 Modification Date	2017-09-21T23:01:27
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




N-Benzyl-2-(2-nitro-8-oxoimidazo[1,2-*a*]pyrazin-7(8*H*)-yl)acetamide (17e)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_062_ppt.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	32
11 Receiver Gain	140.4
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2017-09-17T13:21:29
17 Modification Date	2017-09-17T13:21:29
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




N-Benzyl-2-(2-nitro-8-oxoimidazo[1,2-*a*]pyrazin-7(8*H*)-yl)acetamide (17e)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_062_ppt.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	7000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-09-17T17:15:58
17 Modification Date	2017-09-17T17:15:58
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-(2-Nitro-8-oxoimidazo[1,2-a]pyrazin-7(8H)-yl)-N-(3-(trifluoromethoxy)benzyl)acetamide (17f)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_111_ppt_ACN.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	32
11 Receiver Gain	116.1
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2017-10-14T16:29:47
17 Modification Date	2017-10-14T16:29:47
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-(2-Nitro-8-oxoimidazo[1,2-a]pyrazin-7(8H)-yl)-N-(3-(trifluoromethoxy)benzyl)acetamide (17f)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_111_ppt_ACN.3.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	10000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-10-18T22:58:00
17 Modification Date	2017-10-18T22:58:00
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-(2-Nitro-8-oxoimidazo[1,2-a]pyrazin-7(8H)-yl)ethyl acetate (17g)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8412_082_crude.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	109.3
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	2.2807
16 Acquisition Date	2017-06-01T10:13:19
17 Modification Date	2017-06-01T10:13:19
18 Spectrometer Frequency	600.13
19 Spectral Width	7183.9
20 Lowest Frequency	-592.6
21 Nucleus	1H
22 Acquired Size	16384
23 Spectral Size	65536




2-(2-Nitro-8-oxoimidazo[1,2-a]pyrazin-7(8H)-yl)ethyl acetate (17g)




^{13}C NMR (150 MHz, $\text{DMSO}-d_6$)

Parameter	Value
1 Title	CWA8412_082_crude.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zgpg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	2560
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.9044
16 Acquisition Date	2017-06-01T11:38:16
17 Modification Date	2017-06-01T11:38:16
18 Spectrometer Frequency	150.92
19 Spectral Width	36231.9
20 Lowest Frequency	-3025.7
21 Nucleus	^{13}C
22 Acquired Size	32768
23 Spectral Size	65536




7-(2-Hydroxyethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (17h)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8412_083_purified.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	32
11 Receiver Gain	116.1
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2017-06-06T18:17:20
17 Modification Date	2017-06-06T18:17:20
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1211.0
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




7-(2-Hydroxyethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (17h)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8412_083_purified.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	8000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-06-06T22:34:17
17 Modification Date	2017-06-06T22:34:17
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




7-(2-(Benzyloxy)ethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (17i)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8412_112IIIA_f16-17+2f13-14.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse	zg
Sequence	
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	140.4
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	2.2807
16 Acquisition Date	2017-07-27T11:58:39
17 Modification Date	2017-07-27T11:58:39
18 Spectrometer Frequency	600.13
19 Spectral Width	7183.9
20 Lowest Frequency	-594.7
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




7-(2-(Benzyloxy)ethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (17i)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8412_112IIA_f16-17+2f13-14.3.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zgpg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	3000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.9044
16 Acquisition Date	2017-07-27T13:49:29
17 Modification Date	2017-07-27T13:49:29
18 Spectrometer Frequency	150.92
19 Spectral Width	36231.9
20 Lowest Frequency	-3188.5
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(2-((4-(trifluoromethoxy)benzyl)oxy)ethyl)imidazo[1,2-a]pyrazin-8(7H)-one (17j)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8692_042_f9-11_C18_co mbined1.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	32
11 Receiver Gain	140.4
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2017-09-01T16:44:02
17 Modification Date	2017-09-01T16:44:02
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1211.4
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(2-((4-(trifluoromethoxy)benzyl)oxy)ethyl)imidazo[1,2-*a*]pyrazin-8(7*H*)-one (17j)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_042_f9-11_C18_combined1.4.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	10000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-09-01T22:33:15
17 Modification Date	2017-09-01T22:33:15
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2313.9
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




7-(2-((4-Methylbenzyl)oxy)ethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (17k)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8692_045_ppt.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2017-09-04T16:42:44
17 Modification Date	2017-09-04T16:42:44
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1212.4
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




7-((4-Methylbenzyl)oxy)ethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (17k)




¹³C NMR (150 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8692_045_ppt.3.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	10000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2017-09-04T22:11:13
17 Modification Date	2017-09-04T22:11:13
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(thiazol-5-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26a)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8967_023_ppt_DCM_MeOH.6.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	10.0000
13 Pulse Width	7.0100
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-04-12T16:23:02
17 Modification Date	2018-04-12T16:23:02
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(thiazol-5-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26a)




^{13}C NMR (150 MHz, DMSO- d_6)

Parameter	Value
1 Title	CWA8967_023_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	9500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-04-10T23:00:22
17 Modification Date	2018-04-10T23:00:22
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	^{13}C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(thiazol-2-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26b)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8805/ CWA8805_066/ NMR/ CWA8805_066_ppt_DCM_MeOH/ 1/ fid
2 Title	CWA8805_066_ppt_DCM_MeOH.1.fid
3 Comment	¹ H NMR of CWA8805_066_ppt_DCM_MeOH 1.7mg in 500uL DMSO-d ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	zg
12 Experiment	1D
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	16
15 Receiver Gain	116.1
16 Relaxation Delay	5.0000
17 Pulse Width	7.0100
18 Presaturation Frequency	
19 Acquisition Time	1.9464
20 Acquisition Date	2018-01-25T16:46:34
21 Modification Date	2018-01-25T16:46:34
22 Class	
23 Spectrometer Frequency	600.13
24 Spectral Width	8417.5
25 Lowest Frequency	-1208.1
26 Nucleus	¹ H
27 Acquired Size	16384
28 Spectral Size	65536




2-Nitro-7-(thiazol-2-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26b)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8805/ CWA8805_066/ NMR/ CWA8805_066_ppt_DCM_MeOH/ 2/ fid
2 Title	CWA8805_066_ppt_DCM_MeOH.2.fid
3 Comment	¹³ C JMOD NMR of CWA8805_066_ppt_DCM_MeOH 1.7mg in 500uL DMSO- <i>d</i> ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	jmod
12 Experiment	JMOD
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	9000
15 Receiver Gain	191.6
16 Relaxation Delay	1.0000
17 Pulse Width	11.5000
18 Presaturation Frequency	
19 Acquisition Time	0.8738
20 Acquisition Date	2018-01-27T13:55:32
21 Modification Date	2018-01-27T13:55:32
22 Class	
23 Spectrometer Frequency	150.92
24 Spectral Width	37500.0
25 Lowest Frequency	-2151.2
26 Nucleus	¹³ C
27 Acquired Size	32768
28 Spectral Size	65536




7-((1-Methyl-1*H*-imidazol-2-yl)methyl)-2-nitroimidazo[1,2-*a*]pyrazin-8(7*H*)-one (26c)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_059_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-07-04T13:44:59
17 Modification Date	2018-07-04T13:44:59
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




7-((1-Methyl-1*H*-imidazol-2-yl)methyl)-2-nitroimidazo[1,2-*a*]pyrazin-8(7*H*)-one (26c)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_059_ppT_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	4096
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-07-04T16:01:59
17 Modification Date	2018-07-04T16:01:59
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




7-((1-Methyl-1H-pyrazol-5-yl)methyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (26d)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8967/ CWA8967_005/ NMR/ CWA8967_005_ppt_DCM_MeOH/ 1/ fid
2 Title	CWA8967_005_ppt_DCM_MeOH.1.fid
3 Comment	¹ H NMR of CWA8967_005_ppt_DCM_MeOH 1.5mg in 500uL DMSO-d ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	zg
12 Experiment	1D
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	32
15 Receiver Gain	96.4
16 Relaxation Delay	5.0000
17 Pulse Width	7.0100
18 Presaturation Frequency	
19 Acquisition Time	1.9464
20 Acquisition Date	2018-03-16T17:02:24
21 Modification Date	2018-03-16T17:02:24
22 Class	
23 Spectrometer Frequency	600.13
24 Spectral Width	8417.5
25 Lowest	-1208.1




7-((1-Methyl-1H-pyrazol-5-yl)methyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (26d)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8967/ CWA8967_005/ NMR/ CWA8967_005_ppt_DCM_MeOH/ 2/ fid
2 Title	CWA8967_005_ppt_DCM_MeOH.2.fid
3 Comment	¹³ C JMOD NMR of CWA8967_005_ppt_DCM_MeOH 1.5mg in 500uL DMSO- <i>d</i> ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	jmod
12 Experiment	JMOD
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	9500
15 Receiver Gain	191.6
16 Relaxation Delay	1.0000
17 Pulse Width	11.5000
18 Presaturation Frequency	
19 Acquisition Time	0.8738
20 Acquisition Date	2018-03-16T22:13:59
21 Modification Date	2018-03-16T22:13:59
22 Class	
23 Spectrometer Frequency	150.92
24 Spectral Width	37500.0
25 Lowest Frequency	-2151.2
26 Nucleus	¹³ C
27 Acquired Size	32768




7-((1*H*-benzo[*d*]imidazol-2-yl)methyl)-2-nitroimidazo[1,2-*a*]pyrazin-8(7*H*)-one (26e)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8967/ CWA8967_007/ NMR/ CWA8967_007_f7-9_ppt_DCM_MeOH/ 1/ fid
2 Title	CWA8967_007_f7-9_ppt_DCM_MeOH.1.fid
3 Comment	1H NMR of CWA8967_007_f7-9_ppt_DCM_MeOH 1.5mg in 500uL DMSO-d6; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	zg
12 Experiment	1D
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	16
15 Receiver Gain	59.7
16 Relaxation Delay	5.0000
17 Pulse Width	7.0100
18 Presaturation Frequency	
19 Acquisition Time	1.9464
20 Acquisition Date	2018-04-03T16:07:29
21 Modification Date	2018-04-03T16:07:29
22 Class	
23 Spectrometer Frequency	600.13
24 Spectral Width	8417.5
25 Lowest Frequency	-1208.1
26 Nucleus	1H
27 Acquired Size	16384
28 Spectral Size	65536




7-((1*H*-benzo[*d*]imidazol-2-yl)methyl)-2-nitroimidazo[1,2-*a*]pyrazin-8(7*H*)-one (26e)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_007_f7-9_ppt_DCM_MeOH.4.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	9000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	11.5000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-04-03T22:50:43
17 Modification Date	2018-04-03T22:50:43
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




7-((1-Methyl-1H-benzo[d]imidazol-2-yl)methyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (26f)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8967/ CWA8967_046/ NMR/ CWA8967_046_ppt_DCM_MeOH/ 1/ fid
2 Title	CWA8967_046_ppt_DCM_MeOH.1.fid
3 Comment	¹ H NMR of CWA8967_046_ppt_DCM_MeOH 1mg in 500uL DMSO-d ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	zg
12 Experiment	1D
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	64
15 Receiver Gain	109.3
16 Relaxation Delay	5.0000
17 Pulse Width	8.0000
18 Presaturation Frequency	
19 Acquisition Time	1.9464
20 Acquisition Date	2018-05-09T12:07:23
21 Modification Date	2018-05-09T12:07:24
22 Class	
23 Spectrometer Frequency	600.13
24 Spectral Width	8417.5
25 Lowest Frequency	-1208.1




7-((1-Methyl-1H-benzo[d]imidazol-2-yl)methyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (26f)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_046_ppt_DCM_MeOH.4.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	11000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-05-10T00:31:53
17 Modification Date	2018-05-10T00:31:53
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(pyridin-3-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26g)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8692_064_repurified_com bined1.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	291.7
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	32
11 Receiver Gain	155.2
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2019-06-19T17:14:50
17 Modification Date	2019-06-19T17:14:50
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	1H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(pyridin-3-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26g)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8692_064_repurified_com bined1.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	291.7
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	10000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2019-06-19T22:47:09
17 Modification Date	2019-06-19T22:47:09
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	13C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((6-(trifluoromethyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (26h)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9171_103_filtrate_f3_cry oprobe.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2019-05-01T17:17:49
17 Modification Date	2019-05-01T17:17:49
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((6-(trifluoromethyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (26h)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9171_103_filttrate_f3_cry oprobe.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	7500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2019-05-01T21:18:17
17 Modification Date	2019-05-01T21:18:17
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	13C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(pyrazin-2-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26i)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8805/ CWA8805_111/ NMR/ CWA8805_111_ppt_DCM_MeOH/ 1/ fid
2 Title	CWA8805_111_ppt_DCM_MeOH.1.fid
3 Comment	¹ H NMR of CWA8805_111_ppt_DCM_MeOH 1mg in 480uL DMSO-d ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	zg
12 Experiment	1D
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	32
15 Receiver Gain	116.1
16 Relaxation Delay	5.0000
17 Pulse Width	7.0100
18 Presaturation Frequency	
19 Acquisition Time	1.9464
20 Acquisition Date	2018-03-09T15:03:45
21 Modification Date	2018-03-09T15:03:45
22 Class	
23 Spectrometer Frequency	600.13
24 Spectral Width	8417.5
25 Lowest Frequency	-1208.1
26 Nucleus	¹ H
27 Acquired Size	16384




2-Nitro-7-(pyrazin-2-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26i)




¹³C NMR (150 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8805/ CWA8805_111/ NMR/ CWA8805_111_ppt_DCM_MeOH/ 2/ fid
2 Title	CWA8805_111_ppt_DCM_MeOH.2.fid
3 Comment	¹³ C JMOD NMR of CWA8805_111_ppt_DCM_MeOH 1mg in 480uL DMSO-d ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	jmod
12 Experiment	JMOD
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	9000
15 Receiver Gain	191.6
16 Relaxation Delay	1.0000
17 Pulse Width	11.5000
18 Presaturation Frequency	
19 Acquisition Time	0.8738
20 Acquisition Date	2018-03-09T21:58:02
21 Modification Date	2018-03-09T21:58:02
22 Class	
23 Spectrometer Frequency	150.92
24 Spectral Width	37500.0
25 Lowest Frequency	-2151.2
26 Nucleus	¹³ C
27 Acquired Size	32768
28 Spectral Size	65536




2-Nitro-7-(quinolin-6-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26j)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ Laboratory/NoteBooks/ CWA8967/ CWA8967_044/ NMR/ CWA8967_044_f16-40_ppt_D CM_MeOH/ 1/ fid
2 Title	CWA8967_044_f16-40_ppt_D CM_MeOH.1.fid
3 Comment	¹ H NMR of CWA8967_044_f16-40_ppt_D CM_MeOH 1.7mg in 500uL DMSO-d ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	zg
12 Experiment	1D
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	32
15 Receiver Gain	140.4
16 Relaxation Delay	5.0000
17 Pulse Width	8.0000
18 Presaturation Frequency	
19 Acquisition Time	1.9464
20 Acquisition Date	2018-05-11T11:21:13
21 Modification Date	2018-05-11T11:21:13
22 Class	
23 Spectrometer Frequency	600.13
24 Spectral Width	8417.5
25 Lowest Frequency	-1208.1
26 Nucleus	¹ H
27 Acquired Size	16384




2-Nitro-7-(quinolin-6-ylmethyl)imidazo[1,2-a]pyrazin-8(7H)-one (26j)




¹³C NMR (150 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8967/ CWA8967_044/ NMR/ CWA8967_044_f16-40_ppt_D CM_MeOH/ 2/ fid
2 Title	CWA8967_044_f16-40_ppt_D CM_MeOH.2.fid
3 Comment	¹³ C JMOD NMR of CWA8967_044_f16-40_ppt_D CM_MeOH 1.7mg in 500uL DMSO-d ₆ ; cryoprobe
4 Origin	Bruker BioSpin GmbH
5 Owner	nmr
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	jmod
12 Experiment	JMOD
13 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
14 Number of Scans	7500
15 Receiver Gain	191.6
16 Relaxation Delay	1.0000
17 Pulse Width	12.0000
18 Presaturation Frequency	
19 Acquisition Time	0.8738
20 Acquisition Date	2018-05-15T20:10:48
21 Modification Date	2018-05-15T20:10:49
22 Class	
23 Spectrometer Frequency	150.92
24 Spectral Width	37500.0
25 Lowest Frequency	-2151.2
26 Nucleus	¹³ C
27 Acquired Size	32768
28 Spectral Size	65536




7-(Benzo[d][1,3]dioxol-5-ylmethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (26k)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8805/ CWA8805_073/ NMR/ CWA8805_073_ppt_DCM_MeOH/ 1/ fid
2 Title	CWA8805_073_ppt_DCM_MeOH.1.fid
3 Comment	¹ H NMR of CWA8805_073_ppt_DCM_MeOH 2.2mg in 500uL DMSO-d ₆ ; LC-NMR AA_1H DMSO C:\ c.ang 10
4 Origin	Bruker BioSpin GmbH
5 Owner	biodiversity
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	zg
12 Experiment	1D
13 Probe	5 mm PASEI 1H/ D-13C Z-GRD Z866801/ 0003
14 Number of Scans	16
15 Receiver Gain	287.0
16 Relaxation Delay	5.0000
17 Pulse Width	9.8000
18 Presaturation Frequency	
19 Acquisition Time	1.9465
20 Acquisition Date	2018-01-31T17:24:00
21 Modification Date	2018-01-31T17:25:06
22 Class	
23 Spectrometer Frequency	600.08




7-(Benzo[d][1,3]dioxol-5-ylmethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (26k)




¹³C NMR (150 MHz, DMSO-d₆)

Parameter	Value
1 Data File Name	Z:/ Projects/ 029_Anaerobic/ 029E_Chemistry/ LaboratoryNoteBooks/ CWA8805/ CWA8805_073/ NMR/ CWA8805_073_ppt_DCM_MeOH/ 2/ fid
2 Title	CWA8805_073_ppt_DCM_MeOH.2.fid
3 Comment	¹³ C JMOD NMR of CWA8805_073_ppt_DCM_MeOH 2.2mg in 500uL DMSO-d ₆ ; LC-NMR A_13C_JMOD DMSO C:\ \ c.ang 10
4 Origin	Bruker BioSpin GmbH
5 Owner	biodiversity
6 Site	
7 Instrument	spect
8 Author	
9 Solvent	DMSO
10 Temperature	298.0
11 Pulse Sequence	jmod
12 Experiment	JMOD
13 Probe	5 mm PASEI 1H/ D-13C Z-GRD Z866801/ 0003
14 Number of Scans	10000
15 Receiver Gain	2050.0
16 Relaxation Delay	1.0000
17 Pulse Width	14.0000
18 Presaturation Frequency	
19 Acquisition Time	0.8739
20 Acquisition Date	2018-01-31T22:30:00
21 Modification Date	2018-01-31T22:47:51
22 Class	
23 Spectrometer Frequency	150.91
24 Spectral Width	37500.0
25 Lowest Frequency	-2227.5
26 Nucleus	¹³ C
27 Acquired Size	32768
28 Spectral Size	65536




7-([1,1'-Biphenyl]-4-ylmethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (34a)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8967_048_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	32
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-05-11T11:27:49
17 Modification Date	2018-05-11T11:27:49
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	1H
22 Acquired Size	16384
23 Spectral Size	65536




7-([1,1'-Biphenyl]-4-ylmethyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (34a)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_048_ppT_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	6000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-05-11T14:54:55
17 Modification Date	2018-05-11T14:54:55
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((4'-(trifluoromethoxy)-[1,1'-biphenyl]-4-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34b)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_089_ppt_MeOH_acetone.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9465
16 Acquisition Date	2018-07-27T11:20:00
17 Modification Date	2018-07-28T12:13:16
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((4'-(trifluoromethoxy)-[1,1'-biphenyl]-4-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34b)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_089_ppt_MeOH_acetone.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse	jmod
Sequence	
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	10240
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-07-28T17:41:24
17 Modification Date	2018-07-28T17:41:25
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




7-((4'-Methyl-[1,1'-biphenyl]-4-yl)methyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (34c)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9117_033_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-09-07T15:10:38
17 Modification Date	2018-09-07T15:10:38
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	1H
22 Acquired Size	16384
23 Spectral Size	65536




7-((4'-Methyl-[1,1'-biphenyl]-4-yl)methyl)-2-nitroimidazo[1,2-a]pyrazin-8(7H)-one (34c)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9117_033_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	5120
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-09-07T19:54:41
17 Modification Date	2018-09-07T19:54:41
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((3'-(trifluoromethoxy)-[1,1'-biphenyl]-4-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34d)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_106_ppt_MeOH_DCM.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	109.3
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-08-16T18:02:27
17 Modification Date	2018-08-16T18:02:27
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((3'-(trifluoromethoxy)-[1,1'-biphenyl]-4-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34d)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_106_ppt_MeOH_DCM.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	9000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-08-16T22:51:47
17 Modification Date	2018-08-16T22:51:47
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(4-(pyridin-3-yl)benzyl)imidazo[1,2-a]pyrazin-8(7H)-one (34e)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9171_097_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2019-04-24T17:07:32
17 Modification Date	2019-04-24T17:07:32
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	1H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(4-(pyridin-3-yl)benzyl)imidazo[1,2-a]pyrazin-8(7H)-one (34e)




¹³C NMR (150 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9171_097_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	7500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2019-04-24T21:08:24
17 Modification Date	2019-04-24T21:08:24
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((6-phenylpyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34f)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9117_042_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-09-10T17:16:22
17 Modification Date	2018-09-10T17:16:22
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((6-phenylpyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34f)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9117_042_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	6500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-09-10T21:01:32
17 Modification Date	2018-09-10T21:01:32
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((6-(4-(trifluoromethoxy)phenyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34g)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9117_058_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	155.2
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-09-12T13:08:02
17 Modification Date	2018-09-12T13:08:02
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((6-(4-(trifluoromethoxy)phenyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34g)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9117_058_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	6500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-09-14T20:27:30
17 Modification Date	2018-09-14T20:27:30
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((6-(*p*-tolyl)pyridin-3-yl)methyl)imidazo[1,2-*a*]pyrazin-8(7*H*)-one (34h)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9117_023_ppt_MeOH_DCM.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	155.2
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-09-03T11:12:58
17 Modification Date	2018-09-03T11:12:58
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((6-(*p*-tolyl)pyridin-3-yl)methyl)imidazo[1,2-*a*]pyrazin-8(7*H*)-one (34h)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9117_023_ppt_MeOH_DCM.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	4600
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-09-03T14:00:53
17 Modification Date	2018-09-03T14:00:53
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((6-(3-(trifluoromethoxy)phenyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34i)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8967_114_ppt_MeOH_DCM.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	116.1
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-08-16T17:57:51
17 Modification Date	2018-08-16T17:57:51
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((6-(3-(trifluoromethoxy)phenyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34i)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_114_ppT_MeOH_DCM.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	9500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-08-17T05:55:30
17 Modification Date	2018-08-17T05:55:30
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((5-(4-(trifluoromethoxy)phenyl)pyridin-2-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34j)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9171_091_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	9000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2019-04-18T03:42:49
17 Modification Date	2019-04-18T03:42:49
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((5-(4-(trifluoromethoxy)phenyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34k)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9117_099_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	140.4
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-10-07T11:53:36
17 Modification Date	2018-10-07T11:53:36
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((5-(4-(trifluoromethoxy)phenyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34k)




^{13}C NMR (150 MHz, DMSO- d_6)

Parameter	Value
1 Title	CWA9117_099_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	9500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-10-07T17:18:52
17 Modification Date	2018-10-07T17:18:52
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	^{13}C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((5-(*p*-tolyl)pyridin-3-yl)methyl)imidazo[1,2-*a*]pyrazin-8(7*H*)-one (34l)




¹H NMR (600 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9171_090_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2019-04-17T17:05:39
17 Modification Date	2019-04-17T17:05:39
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((5-(*p*-tolyl)pyridin-3-yl)methyl)imidazo[1,2-*a*]pyrazin-8(7*H*)-one (34l)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9171_090_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	6500
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2019-04-17T20:34:03
17 Modification Date	2019-04-17T20:34:04
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((5-(3-(trifluoromethoxy)phenyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34m)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA91117_115_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	155.2
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-10-16T17:18:51
17 Modification Date	2018-10-16T17:18:51
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	1H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((5-(3-(trifluoromethoxy)phenyl)pyridin-3-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34m)




¹³C NMR (150 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9117_115_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	10000
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-10-16T22:39:45
17 Modification Date	2018-10-16T22:39:45
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	¹³ C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(4-(thiazol-2-yl)benzyl)imidazo[1,2-a]pyrazin-8(7H)-one (34n)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9117_072_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	biodiversity
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	5 mm PASEI 1H/ D-13C Z-GRD Z866801/0003
10 Number of Scans	16
11 Receiver Gain	287.0
12 Relaxation Delay	5.0000
13 Pulse Width	9.8000
14 Presaturation Frequency	
15 Acquisition Time	1.9465
16 Acquisition Date	2018-09-18T14:53:00
17 Modification Date	2018-09-18T14:53:58
18 Spectrometer Frequency	600.08
19 Spectral Width	8417.5
20 Lowest Frequency	-308.1
21 Nucleus	1H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(4-(thiazol-2-yl)benzyl)imidazo[1,2-a]pyrazin-8(7H)-one (34n)




¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA9117_072_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	biodiversity
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	297.9
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	5 mm PASEI 1H/ D-13C Z-GRD Z866801/0003
10 Number of Scans	7500
11 Receiver Gain	2050.0
12 Relaxation Delay	1.0000
13 Pulse Width	14.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8739
16 Acquisition Date	2018-09-18T21:45:00
17 Modification Date	2018-09-18T21:45:36
18 Spectrometer Frequency	150.91
19 Spectral Width	37500.0
20 Lowest Frequency	-2227.5
21 Nucleus	13C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-((2-phenylthiazol-5-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34o)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA9117_100_ppt_DCM_MeOH.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	16
11 Receiver Gain	109.3
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-10-23T12:59:20
17 Modification Date	2018-10-23T12:59:20
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	1H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-((2-phenylthiazol-5-yl)methyl)imidazo[1,2-a]pyrazin-8(7H)-one (34o)




^{13}C NMR (150 MHz, DMSO- d_6)

Parameter	Value
1 Title	CWA9117_100_ppt_DCM_MeOH.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	jmod
8 Experiment	JMOD
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	5200
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.8738
16 Acquisition Date	2018-10-23T15:46:56
17 Modification Date	2018-10-23T15:46:56
18 Spectrometer Frequency	150.92
19 Spectral Width	37500.0
20 Lowest Frequency	-2151.2
21 Nucleus	^{13}C
22 Acquired Size	32768
23 Spectral Size	65536




2-Nitro-7-(4-(piperidin-1-yl)benzyl)imidazo[1,2-a]pyrazin-8(7H)-one (34p)




¹H NMR (600 MHz, DMSO-d₆)

Parameter	Value
1 Title	CWA8967_055_ppt_DCM_Me OH_cryoprobe.1.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	32
11 Receiver Gain	59.7
12 Relaxation Delay	5.0000
13 Pulse Width	8.0000
14 Presaturation Frequency	
15 Acquisition Time	1.9464
16 Acquisition Date	2018-05-16T02:12:37
17 Modification Date	2018-05-16T02:12:37
18 Spectrometer Frequency	600.13
19 Spectral Width	8417.5
20 Lowest Frequency	-1208.1
21 Nucleus	¹ H
22 Acquired Size	16384
23 Spectral Size	65536




2-Nitro-7-(4-(piperidin-1-yl)benzyl)imidazo[1,2-a]pyrazin-8(7H)-one (34p)



¹³C NMR (150 MHz, DMSO-*d*₆)

Parameter	Value
1 Title	CWA8967_055_ppt_DCM_MeO H_cryoprobe.2.fid
2 Origin	Bruker BioSpin GmbH
3 Owner	nmr
4 Spectrometer	spect
5 Solvent	DMSO
6 Temperature	298.0
7 Pulse Sequence	zgpg
8 Experiment	1D
9 Probe	Z129649_0009 (CP TCI 600S3 H&F-C/ N-D-05 Z)
10 Number of Scans	2560
11 Receiver Gain	191.6
12 Relaxation Delay	1.0000
13 Pulse Width	12.0000
14 Presaturation Frequency	
15 Acquisition Time	0.9044
16 Acquisition Date	2018-05-16T03:36:59
17 Modification Date	2018-05-16T03:36:59
18 Spectrometer Frequency	150.92
19 Spectral Width	36231.9
20 Lowest Frequency	-3025.7
21 Nucleus	13C
22 Acquired Size	32768
23 Spectral Size	65536



References

1. Sykes, M. L.; Baell, J. B.; Kaiser, M.; Chatelain, E.; Moawad, S. R.; Ganame, D.; Ioset, J.-R.; Avery, V. M. Identification of compounds with anti-proliferative activity against *Trypanosoma brucei* strain 427 by a whole cell viability based HTS campaign. *PLoS Negl. Trop. Dis.* **2012**, 6, e1896-e1896.