

**TABLE S2: Pulse Sequences and Parameters Used for T2-Weighted MRI of Seminal Vesicle Evaluation**

Imaging Plane	FOV (mm)	Resolution (mm)	Section Thickness (mm)	No. of Sections	TR/TE	Echo-Train Length	Flip Angle (°)	Duration
Axial	140	0.3	3	26	4400/120	26	90	2 min 48 s
Sagittal	140	0.3	3	24	2900/120	26	90	1 min 50 s
Coronal	140	0.3	3	20	2400/120	26	90	2 min 20 s