Table SI. Colorimetric result of the semen sample interpreted by the smartphone and the corresponding total motile sperm count

Category: Low TMSC																						
Sample No.	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
Test strip color																						
TMSC	0	0	0	0	1	0	0	0	0	0	0	0	0	0	0	0	3.2	15	14	13	8.7	3.3
Result	\checkmark	\checkmark	\checkmark	\checkmark	Х	\checkmark																

Category: Normal TMSC																				
Sample No.	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42
Test strip color																				
TMSC	67	56	248	86	775	67	46	92	122	86	266	140	96	74	41	36	80	58	27	194
Result	\checkmark	Х	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	\checkmark	Х	\checkmark	Х	Х	Х	Х	\checkmark	\checkmark	\checkmark	\checkmark	Х	\checkmark

TMSC: total motile sperm count

*All the sample colors are changed to the gray scale and subtracted by the blank value