



A**B**

9 microdevices

C4 arrays of 9 microdevices glued
in a 35mm culture plate

VIDEOMICROSCOPY

Supplementary Figure S5. Microdevice to study clustering at the single-cell scale

a Mask used for the fabrication of the silicon wafer. **b** One array of 9 PDMS micro-wells (outer diameter: 650 μ m, inner diameter: 450 μ m, and height: 200 μ m) that are **(c)** glued on the bottom of the compartments of CELLview™ cell culture dishes for monitoring by time-lapse video-microscopy.