

Supplemental Figure 1

Main

Confirm Setup

Patient ID Operator Location/Site ID MRI Info

Treatment ID Notes

NOTE: Select Modify Setup to change values below

Transducer ID XDCR Connection Frequency (KHz)

Pr.0 (MPa)	Ispta.0(mW/cm ²)	Pulse Width (ms)	PRF (Hz)
<input type="text" value="0.636"/>	<input type="text" value="1017"/>	<input type="text" value="5"/>	<input type="text" value="10"/>

Pr.3(MPa)	Ispta.3(mW/cm ²)	Duty Factor	XMT Volts
<input type="text" value="0.535"/>	<input type="text" value="719.73"/>	<input type="text" value="0.05"/>	<input type="text" value="14.2"/>

Pr.a (MPa)	Ispta.a(mW/cm ²)	Skull Thickness (mm)
<input type="text" value="0.301"/>	<input type="text" value="227.7"/>	<input type="text" value="5"/>

Sonication Duration (s)

Time between Sonications (s)

Number of Sonications

MRI Sync Enable

MRI Sync Delay(ms)

Sonication Indicator

Sonication File

BC Pulsar 1002 Tabletop
Monday, October 12, 2020 01:47:31 PM
Version 3.0.0