

On-line Table: Summary of MR imaging sequence parameters for included patients

Patient	Center	Manufacturer	Model	TR (ms)	TE (ms)	No. of Averages ^a	Field Strength	Spacing between Sections (mm)	Flip Angle
1	1	Siemens ^b	Avanto	24	7	1	1.5T	0.7	25°
2	1	Siemens	Avanto	25	7.15	1	1.5T	0.7	25°
3	2	GE Healthcare ^c	Signa HDxt	26	3.6	1	1.5T	0.5	17°
4	3	Philips Healthcare	Achieva	25	3.42	1	3T	0.6	20°
5	4	Philips Healthcare ^d	Achieva	25	3.42	1	3T	0.6	20°
6	4	GE Healthcare	Unspecified	11.3	4.84	0.63	1.5T	0.7	10°
7	5	GE Healthcare	Signa HDxt	Min	Min	1	1.5T	0.5	17°
8	6	GE Healthcare	Signa HDxt	Min	Min	1	1.5T	0.5	17°
9	6	GE Healthcare	Signa HDxt	Min	Min	1	1.5T	0.5	17°
10	6	Philips Healthcare	Intera	23	2.34	1	1.5T	0.5	20°
11	7	Philips Healthcare	Achieva	25	3.42	1	3T	0.6	20°
12	8	Philips Healthcare	Intera	23	3.09	1	1.5T	0.5	20°
13	9	Siemens	Avanto	25	7.15	1	1.5T	0.7	25°
14	9	Philips Healthcare	Intera	23	3.09	1	1.5T	0.5	20°
15	9	Philips Healthcare	Intera	23	3.09	1	1.5T	0.5	20°
16	9	Siemens	Avanto	25	7.15	1	1.5T	0.7	25°
17	9	Siemens	Avanto	25	7.15	1	1.5T	0.7	25°
18	9	Philips Healthcare	Intera	23	3.09	1	1.5T	0.5	20°
19	9	Siemens	Skyra	20	3.57	1	3T	0.98	18°
20	9	Siemens	Avanto	25	7.15	1	1.5T	0.7	25°
21	10	Siemens	Sonata	37	4.63	1	1.5T	0.26	25°
22	10	Siemens	Sonata	37	4.63	1	1.5T	0.26	25°
23	11	GE Healthcare	Signa HDxt	26	4	1	3T	0.3	17°
24	12	GE Healthcare	Signa HDxt	30	3	1	3T	0.29	15°
25	6	GE Healthcare	Signa HDxt	Min	Min	1	1.5T	0.5	17°
26	6	GE Healthcare	Signa HDxt	Min	Min	1	1.5T	0.5	17°
27	6	GE Healthcare	Discovery MR750	Min	Min	1	3T	0.5	15°
28	6	GE Healthcare	Discovery MR750	Min	Min	1	3T	0.5	15°
29	6	GE Healthcare	Discovery MR750	Min	Min	1	3T	0.5	15°
30	6	GE Healthcare	Discovery MR750	Min	Min	1	3T	0.5	15°
31	6	GE Healthcare	Discovery MR750	Min	Min	1	3T	0.5	15°

Note:—Min indicates minimum.

^a Number of averages is an absolute number.

^b Erlangen, Germany.

^c Milwaukee, Wisconsin.

^d Best, the Netherlands.