

Supplementary Table S2. Antibodies used for multiparametric immunofluorescence stainings.

Antigen	Clone	Dilution	Source
CD3	polyclonal	1:150	Dako
CD4	EPR6855	1:50	Abcam
CD4	SP35	prediluted	Ventana Medical Systems
CD8	C8/144B	1:100	Dako
CD8	SP57	prediluted	Ventana Medical Systems
CD20	L26	1:500/1:100	Thermo Fisher Scientific
BDCA-2	polyclonal	1:50/1:75	R&D Systems
FOXP3	236A/E7	1:100	Abcam
pan-cytokeratin	AE1/AE3	1:100	Thermo Fisher Scientific
pan-cytokeratin	AE1/AE3/PCK26	prediluted	Ventana Medical Systems
GrzB	Grb-7	1:50	Dako
IRF7	G-8	1:75	Santa Cruz Biotechnology