

Supplementary Table 2. Stability of symmetric dimethylarginine (SDMA), measured in duplicate by immunoassay, in 10 equine serum samples stored at -10, 4, or 25°C for 2 and 7 days.

-20°C					
Sample	Day 0	Day 2	Bias	Day 7	Bias
#1	14.9	13.2	-1.7	13.8	-1.0
#2	16.0	15.0	-1.0	16.2	0.2
#3	15.0	15.2	0.2	14.5	-0.4
#4	16.0	15.2	-0.8	14.5	-1.5
#5	14.4	13.9	-0.5	15.0	0.6
#6	17.0	16.3	-0.7	17.7	0.7
#7	15.5	15.1	-0.4	15.0	-0.4
#8	15.9	16.1	0.2	15.9	0.0
#9	14.2	13.7	-0.5	14.1	-0.1
#10	18.1	16.5	-1.6	16.9	-1.2
4°C					
Sample	Day 0	Day 2	Bias	Day 7	Bias
#1	14.9	13.3	-1.6	14.4	-0.5
#2	16.0	14.8	-1.3	15.8	-0.2
#3	15.0	14.0	-1.0	14.0	-1.0
#4	16.0	14.0	-2.0	15.6	-0.4
#5	14.4	13.6	-0.8	13.3	-1.1
#6	17.0	17.1	0.1	16.5	-0.5
#7	15.5	14.8	-0.6	16.0	0.5
#8	15.9	15.0	-0.9	16.2	0.3
#9	14.2	14.7	0.5	13.1	-1.1
#10	18.1	17.1	-1.0	17.2	-0.9
25°C					
Sample	Day 0	Day 2	Bias	Day 7	Bias
#1	14.9	14.5	-0.4	13.1	-1.8
#2	16.0	15.9	-0.2	14.8	-1.2
#3	15.0	15.6	0.7	15.1	0.1
#4	16.0	15.0	-1.0	15.4	-0.6
#5	14.4	13.6	-0.8	13.6	-0.8
#6	17.0	17.4	0.4	17.1	0.1
#7	15.5	15.4	-0.1	14.9	-0.6
#8	15.9	15.5	-0.4	15.5	-0.4
#9	14.2	14.6	0.4	14.8	0.6
#10	18.1	16.9	-1.2	17.0	-1.1