Draft instructions for self collection


Mouth rinse sample collection for COVID-19 Testing

Please read all instructions before taking your sample

Make sure you are wearing a mask especially if you are doing indoor sample collection



Wash your hands with soap and water or hand sanitizer



Open the pink saline (salt water) vial by twisting the top holding it with your thumb & finger.

Pull your face mask down below your mouth, if you are wearing one.



Squeeze all of the saline into your mouth, but do not swallow it.

Put your face mask back on.



Rinse your mouth for 5 seconds by swishing the saline in your mouth. Then tilt head back and gargle for 5 seconds.

Repeat swishing and gargling 3 times for a total of 30 seconds.



Twist the cap off the container.

Pull your face mask down and spit the liquid from your mouth into the container. Put your face mask back on



Put the lid back on the container and screw on tightly.



Put container in the plastic bag.



Wash your hands with soap and water or



Give the bag back to your health hand sanitizer. care provider.







Saliva Sample




Wash your hands with soap and water.



Spread a paper towel down on a flat surface.





Put the lid with the inside surface facing up on the paper towel.

Don't eat, drink, smoke, or chew gum for <u>I hour</u> before providing your saliva sample.



Swallow and bring up saliva from the back of your throat



Spit into the tube until your saliva reaches the marked line at 10ml.



Put the lid back on the container.



Put container in the plastic bag.

If you're having trouble creating saliva:

- Think about your favorite food
- Gently rub your cheeks from the outside (pressing them against your teeth and gums)



Wash your hands with soap and water.



Give the bag back to the nurse.

