



Supplemental Material

Supplemental Figure S1: *CHEK2* c.475T>C chromatograms



Supplemental Figure S2: Cell sorting strategy



Abbreviations: APC, allophycocyanin; BV785, Brilliant Violet 785; FITC, fluorescein isothiocyanate; FSC, forward scatter; MM, myelomonocytic; PB, peripheral blood; SSC, side scatter.

Supplementary Table S1: QuikChange oligonucleotides used to generate the CHK2 variants

Variant	Oligonucleotide pair
Y159H	F: CCACTGTGATCTTCTATGTGTGCAATGTAAGAGTTTTTAGGAC R: GTCCTAAAAACTCTTACATTGCACACATAGAAGATCACAGTGG
I157T	F: CTTCTATGTATGCAGTGTAAGAGTTTTTAGGACCCACTTCCC R: GGGAAAGTGGGTCCCTAAAAACTCTTACACTGCATACATAGAAG
H143Y	F: TCCCTGAAAATCCGAAAATATTTCTTGCTGTATGTTCCGGTATTTATCTGTT R: AACAGATAAATACCGAACATACAGCAAGAAATATTTTCGGATTTTCAGGGA