



Supplementary Figure 2: Representative electron microscopy negative staining image showing EQ reagent contamination. Electron micrograph shows aggregates of vesicles (red arrowheads). Other unidentified structures seen in the background (yellow arrowheads) could be protein aggregates and/or dense aggregates of EVs. Fibrous-like shapes (blue arrowheads) could result from carry-over of EQ reagent. Scale bar reflects the magnification at the camera. 100k magnification, scale bar 200 nm. EQ: ExoQuick-TC.



© The Author(s) 2021. Open Access This article is licensed under a Creative Commons Attribution 4.0 International License

(<https://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, sharing, adaptation, distribution and reproduction in any medium or

format, for any purpose, even commercially, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made.