

Supplementary table 1. Schedule of IntelliCage analysis

Session name	Days
Free adaptation	2
Nose poke adaptation	3
Drinking session adaptation	3
Place preference test	7
Place preference reversal test	5
Impulsivity test	4
Attention test	5
Place avoidance test	9
Delay discounting test	10

Supplementary table 2. Primer sets in this study.

A. Primer sets for sgRNA (5'->3')

	top	bottom
<i>Macf1</i>	CACCGGTAGGTCTCGTGCCTGCAA	AAACTTGCAGGACACGAGACCTACC
<i>Ehd1</i>	CACCGGGACTTCACCATCTCCTTCT	AAACAGAAGGAGATGGTGAAGTCC

B. ssODN for knock-in of mutations (5'->3')

<i>Macf1</i>	TGATGCCACTGACAGATAAATCACTCTGAGAAGCAGGCATTATGGGAGCCTTGCAGG ACGACGAGACCTACATGGGTAGGCTCTGTGTGAATAGGGTCTCCTGACTGTTTTGGTA TTGAT
<i>Ehd1</i>	TTCTACACACTGTCTCCTGTCAACGGCAAGATCACAGGTGCTAATGCTAAGAA GGAGATGTGAAGTCCAAGCTGCCAACACAGTGCTGGGGAAGATCTGGAAGT TGGCAGATGTGGACA

C. Genotyping primers

	Fw	Rv	PCR conditions
<i>Macf1</i>	CAGGTAAAGGAGAGCTTCCAAGCAG	CACACTCCACTTCGATATGCCTGAG	94°C, 1 min → (98°C, 10 sec → 68°C, 1 min) ×30
<i>Ehd1</i>	CAGACCCAGGACTTCAGCAAGTTCC	GGATGAGTGTGGGGCTGTCAGACAG	

D. Sequencing primers for genotyping

	5'->3'
<i>Macf1</i>	CCTCTGAAACAGTAGATCAAGACCAACC
<i>Ehd1</i>	CCTCCGTTTGGAGGGTGGGAATGAG