

Additional file 1: Table S1. Antibodies used for surface and intracellular cytokine staining

Antibody	Clone	Assay concentration ($\mu\text{g/mL}$)	Company
CD3-BUV395	SK7	0.8	BD Biosciences
CD14-BUV395	M ϕ P9	2.5	BD Biosciences
CD19-BV421	HIB19	1.3	BioLegend
CD24-BV605	ML5	3.0	BioLegend
CD27-BV480	L128	5.0	BD Biosciences
CD38-BV785	HIT2	5.3	BioLegend
IL-10-PE-Dazzle594	JES3-19F1	1.5	BioLegend
EBI3 (IL-27 subunit)- APC (IL-35)	ebic6	0.3	Invitrogen
LAP (TGF- β)-FITC	TW4-2F8	10.0	BioLegend