

Table S1. Primers and probes used. NA = not applicable.

Primer name	Sequence	Catalogue number IDT
E Sarbeco_F1 Forward	ACAGGTACGTTAATAGTTAATAGCGT	10006889
E Sarbeco_R2 Reverse	ATATTGCAGCAGTACGCACACA	10006891
E Sarbeco_P1 Probe	ACACTAGCCATCCTTACTGCGCTTCG	10006893
sgLeadSARSCoV2-F	CGATCTCTTGTAGATCTGTTCTC	NA