

Table 3. Sequence information for ChIP-PCR primers used in described studies:

ZKSCAN3 DP1	F-	TGGTGTTTAAGGCACAGTGGGA
ZKSCAN3 DP1	R-	TGTCATGTGTTTACTGTCCTCCA
ZKSCAN3 DP2	F-	AGGACAGTAAACACATGACATAGCA
ZKSCAN3 DP2	R-	AGGTTCAAGCTCGGGACCTTC
ZKSCAN3 DP3	F-	GGTCCCGAGCTTGAACCTCA
ZKSCAN3 DP3	R-	CTTGAGGGACTGCAGCCTCA
GAPDH		Comes from the ChIP kit (Thermo Fisher. Catalog No.26156)