



Supplementary Figure 7. Presence of ExoS, T and Y or Exo S and ExoT in bacterial culture supernatants of PAO1F and the DY mutant, respectively. Concentrated supernatants were analyzed by PAGE on 4-12 % NuPage gradient gels run in MES buffer and stained with Bio-Safe Coomassie G-250 (BIO - RAD). The identity of the proteins was determined by processing the bands cut from the gel by tryptic digests and analysis by tandem mass spectrometry (LC-MS/MS) revealed. A minimum of 30 % sequence coverage was achieved for protein identification.