

Supplemental Table 4. Cell markers used for flow cytometry cell subset identification for *in vitro* tracer association data in Suppl. Fig. 2.

Cell subset name	Cell markers used for identification
T Cell	CD45 ⁺ CD3 ⁺
B Cell	CD45 ⁺ CD3 ⁻ CD11b ⁻ CD11c ⁻
Dendritic Cell	CD45 ⁺ CD3 ⁻ CD11c ⁺
Macrophage	CD45 ⁺ CD3 ⁻ CD11b ⁺ F4/80 ⁻ SSC ^{lo}