

Supplemental Table 5. Cell markers used for flow cytometry cell subset identification in NO treatment tracer association data in Fig. 3.

Cell subset name	Cell markers used for identification
T Cell	CD45 ⁺ CD3 ⁺ B220 ⁻
B Cell	CD45 ⁺ B220 ⁺ CD3 ⁻ CD11b ⁻
cDC	CD45 ⁺ CD11c ⁺ B220 ⁻
pDC	CD45 ⁺ CD11c ⁺ B220 ⁺
SSM	CD45 ⁺ B220 ⁻ CD3 ⁻ CD11b ⁺ CD11c ⁺ CD169 ⁺ F4/80 ⁻
MCM	CD45 ⁺ B220 ⁻ CD3 ⁻ CD11b ⁺ CD11c ⁺ CD169 ⁺ F4/80 ⁺
MSM	CD45 ⁺ B220 ⁻ CD3 ⁻ CD11b ⁺ CD11c ⁺ CD169 ⁺ F4/80 ⁺
dDC	CD45 ⁺ CD11c ^{hi} DEC205 ⁺
LC	CD45 ⁺ CD11c ^{lo} DEC205 ⁺