
Suppl Table 1: PCR Primers

Gene	Forward	Reverse	
CHGA	GATCCTTTCCATTCTGAGACATCA	GAACCTCTGAGTTCATCTTCA	
SYP	ACAAGACCGAGAGTGACCT	AGTCCCCAACTAAGAAGACCT	
NKX6.3	GGAAGAAGCAGACACTCCTG	TTCCGCCACTTGGTCCT	
NGN3	CTCATCGCTCTCTATTCTTTGC	GGCAGGTCACTTCGTCTTC	
NEUROD	GACTATCACTGCTCAGGACCTA	CACTCGTCTGTCCAGCTTG	
HPRT1	TTGTTGTAGGATATGCCCTTGA	GCGATGTCAATAGGACTCCAG	

Suppl Table 2: Primary & Secondary Antibodies

Antibodies	Vendor	Catalog #	Dilution/Application
Gastrin	Dako Agilent	A0568	1:1000 (IF)
NFкB-p65	Cell Signaling Technologies	8242S	1:100 (IHC/IF) 1:1000 (WB)
Nkx6.3	Novus Biologicals	NBP1-92185	1:50 (IHC) 1:1000 (WB)
Synaptophysin	Cell Signaling Technologies	9020S	1:200 (IF)
α-SMA	Thermo Scientific	MS-113-P0	1:200 (IF)
S100b	Abcam	Ab52642	1:100 (IF)
IL-17b	R&D Systems	AF1248	1:100 (IF)
TNF-α	R&D Systems	AF-410-NA	1:50 (IF)
TNFAIP2	Novus Biologicals	NBP1-33480	1:100 (IF)
β-catenin	R&D Systems	AF1329	1:100 (IF)
pSTAT3 (Y705)	Cell Signaling Technologies	9145S	1:100 (IF)
STAT3	Cell Signaling Technologies	9139	1:1000 (WB)
Histone H3	Cell Signaling Technologies	12648S	1:10000 (WB)
β-Tubulin	Cell Signaling Technologies	5346S	1:2000 (WB)
GAPDH	Cell Signaling Technologies	2118S	1:2000 (WB)
Alexa Fluor Donkey anti-Rabbit IgG 488	Invitrogen	A21206	1:500 (IF)
Alexa Fluor Donkey anti-Mouse IgG 647	Invitrogen	A31571	1:500 (IF)
Alexa Fluor Donkey anti-Goat IgG 594	Invitrogen	A11058	1:500 (IF)


Supplementary Figure S1





Α


В

Supplementary Figure S2













Supplementary Figure S3









Rico K, et al. BMJ Open Gastro 2021; 8:e000765. doi: 10.1136/bmjgast-2021-000765

