

Supplemental Table 2. List of antibodies use for flow cytometry.

Location	Antibody	Clone	Source
	Live/dead	65-0865-18	eBio
	CD45	30-F11	BioLegend
	CD4	RM4-5	BioLegend
	CD8	53-6.7	BioLegend
	NKp46	29A1.4	BioLegend
	ICOS	C398.4A	BioLegend
	PD-1	29F.1A12	BioLegend
	Tim-3	RMT3-23	BioLegend
	Lag-3	eBioC9B7W	eBio
	OX40	OX86	BioLegend
	CD122	TM-B1	BD Biosciences
	CD25	PC61.5	BioLegend
	CD62L	MEL14	eBio
Surface	AH1-A5 tetramer		NIH Tetramer Core
	MHCII	M5/114.15.2	eBio
	CD24	M1/69	BD Biosciences
	F4/80	BM8	BioLegend
	CD11c	N418	BioLegend
	CD11b	M1/70	BD Biosciences
	CD103	2E7	BioLegend
	CD274 (PD-L1)	B7-H1	BioLegend
	Arg1		R&D Systems
	Ly6C	HK1.4	BioLegend
	Ly6G	1A8	BD Biosciences
	4-1BB	17B5	eBio
	CD27	LG.3A10	BioLegend
	CTLA-4	14D3	eBio
	iNos	CXNFT	eBio
Intracellular	Granzyme A	GzA-3G8.5	Invitrogen
	FoxP3	FKJ16a	eBio
	Ki-67	B56	BD Biosciences
	IFN-g	XMG1.2	BioLegend
	TNF-a	TN3-19	BioLegend