
**Supplementary Figure 1** 



## Supplementary Figure 1 legend

## The limit of detection of *polA* and *Tpp47* in plasma and urine by nPCR assay.

Different number of *T. pallidum* (Nichols strains) were diluted in plasma (A) and urine (B) and *polA* gene was detected by nPCR assay; Different number of *T. pallidum* (Nichols strains) were diluted in plasma (C) and urine (D) and *Tpp47* gene was detected by nPCR assay. M: marker; line1, 3\*10<sup>4</sup> *T. pallidum*/mL; line 2, 3\*10<sup>3</sup> *T. pallidum*/mL; line 3, 3\*10<sup>2</sup> *T. pallidum*/mL; line 4, 3\*50 *T. pallidum*/mL; line 5, 3\*10 *T. pallidum*/mL; line 6, 3\*5 *T. pallidum*/mL; line 7, 3 *T. pallidum*/mL; line 8, 0.3 *T. pallidum*/mL; line 9, negative control (plasma or urine without *T. pallidum*).