

Guide for FiO₂ estimation among participants with oxygen supplementation in whom a temporary interruption of the supplementation is judged inappropriate.

Administration device	FiO₂ (Liters)
Nasal canulae	0.24-0.28 (Max 5 lts)
Venturi Mask	0.24 (3 lts)
	0.35 (8 lts)
	0.40 (10 lts)
	0.50 (12 lts)
Reservoir mask	0.55 (6 lts)
	0.60 (7 lts)
	0.70 (8 lts)
	0.80 (9 lts)
	0.90-0.99 (10-15 lts)