

Supplementary Table 1 – Primers and probe sequences.

Target	Sequence (5'-3')
CAR	F: TTTTGGACAAGAGACGTGGC R: CCTTGGTGGCTGTACTGAGA
SSTR2_737	F: CCTAAGACTGAAAAGCAGCCAT
SSTR2_723	R: CCTAAGACTGAAAAGCAGCCAT
SSTR5_1600	F: TTCTTCTCTTGCAGAGCCTGA
SSTR5_1616	R: GACTGTCAGGCAGAAGACACTG
sgRNA_SSTR2	F: CATGAGCATCGACCGATAACC R: TTGCTTGTCAGGTCATAGTA
sgRNA_SSTR5	F: GAACACGCTGGTCATCTACG R: CGCAGAACGCCGCGTCCTTC
CAR-FAM probe	Bio-Rad Id: dCNS496782297*
MKL2-HEX probe	Bio-Rad Id: dHsaCNS879832622**
*MIQE [46] context provided at https://www.bio-rad.com/digital-assays/assay-detail/dCNS496782297	
**MIQE [46] context provided at https://www.bio-rad.com/digital-assays/assay-detail/dHsaCNS879832622	

References

1. Bustin SA, Benes V, Garson JA, Hellemans J, Huggett J, Kubista M *et al.* The MIQE Guidelines: Minimum Information for Publication of Quantitative Real-Time PCR Experiments. *Clin Chem.* **2009**; 55(4):611-22.