

Supplementary Table SI Preimplantation genetic testing for aneuploidies (PGT-As) methods for each clinic.

IVF clinic	Name of testing laboratory	Genetic testing method	Genetic testing platform	Thresholds		
				Euploid	Aneuploid	Mosaic
Ovation—Austin San Antonio IVF Midwest Fertility Specialists California Fertility Partners Ovation—Baton Rouge	Ovation Fertility Genetics Laboratories	NGS	VeriSeq (Illumina, CA, USA)	≤30% abnormal cells	>70% abnormal cells	31–70% abnormal cells (mosaic-low = 31–50%, mosaic-high = 51–70%)
Alpha IVF Women's Specialists	Alpha IVF Women's Specialists	NGS	Ion Torrent (Thermo Fisher Scientific, MA USA)	≤20% abnormal cells	>80% abnormal cells	21–80% abnormal cells (mosaic-low = 21–50%, mosaic-high = 51–80%)
Wings IVF Women's Hospital	Genexplore Diagnostics and Research Centre	NGS	Ion Torrent (Thermo Fisher Scientific, MA USA)	≤20% abnormal cells	>80% abnormal cells	Not provided (no mosaic embryos in dataset)
IVF-Spain—Alicante IVF-Spain—Madrid IVF-Spain—San Sebastián	Integrated Genetic Lab Services	NGS	VeriSeq (Illumina, CA, USA)	≤30% abnormal cells	>70% abnormal cells	31–70% abnormal cells (information for mosaic-low/mosaic-high not provided)

NGS, next generation sequencing.