

Supplemental Table 1. Contamination measured in WEINAS robotic system. — Raw data

WEINAS	Cyclophosphamide			Fluorouracil			Gemcitabine		
	(ng/cm <sup>2</sup> )			(ng/cm <sup>2</sup> )			(ng/cm <sup>2</sup> )		
	Day1	Day2	Day3	Day1	Day2	Day3	Day1	Day2	Day3
Clamping arm*	35400	31700	28700	765	648	627	670	633	523
Disc	40.76	35.67	34.24	0.23	0.20	0.18	10.83	8.04	7.73
Waste inlet	44.70	35.79	29.74	10.96	10.26	9.57	13.29	12.26	11.70
A(inside)	0.71	0.59	0.51	0.058	0.045	0.044	0.11	0.09	0.06
B(inside)	0.34	0.32	0.30	0.01	0.009	0.008	0.011	0.01	0.008
C(inside)	0.39	0.26	0.02	0.03	0.029	0.026	0.054	0.049	0.044
D(inside)	0.45	0.38	0.37	0.027	0.024	0.02	0.034	0.029	0.025
E(inside)	0.53	0.44	0.34	0.091	0.087	0.084	0.037	0.033	0.027
F(outside)	0.16	0.11	0.11	0.028	0.023	0.021	0.015	0.01	0.011
Gloves	14.67	14.40	13.00	2.25	2.15	2.33	0.64	0.61	0.65
Infusion bags	7.20	13.20	11.45	3.94	3.19	3.28	1.92	2.31	3.60