

Supplementary Materials for
Conformationally responsive dyes
enable protein-adaptive differential scanning fluorimetry

Taiasean Wu, Joshua C. Yu, Arundhati Suresh, Zachary J. Gale-Day, Matthew G. Alteen,
Amanda S. Woo, Zoe Millbern, Oleta T. Johnson, Emma C. Carroll, Carrie L. Partch, Denis
Fourches, Nelson R. Vinueza, David J. Vocadlo, Jason E. Gestwicki*

Correspondence to: jason.gestwicki@ucsf.edu

This PDF file includes:

Materials and Methods
Figs. S1 to S18
Tables S1 to S3
Caption for Movie S1
Captions for Data S1 to S4

Other Supplementary Materials for this manuscript include the following:

Movie S1
Data S1 to S4 (Dye library sources and extended details, Full dye screen results for all proteins, Protein panel sources and extended details, Dye screen results for all proteins.)

Materials and Methods

Making and storing dye stock solutions

All commercial dyes were used without further purification. The purity of all Max Weaver dyes (19) was confirmed >80% by LCMS at North Carolina State University prior to addition to the dye library. All LCMS analyses of the Max Weaver library dyes were performed on an Agilent Technologies 1260 High-Performance Liquid Chromatography (HPLC) system coupled with an Agilent 6520 Q-TOF mass spectrometer. (Agilent Technologies, CA, USA). Ionization was performed via electrospray ionization with the following parameters: gas temperature 350 °C, drying gas 6 liters per minute, nebulizer 35 psi, V_{cap} voltage 3500 V and fragmentor voltage at 160 V. Powders were stored in RT, sealed and parafilm, in the dark. From each dye, a 5 mM DMSO was made by transferring ~1 mg of solid dye into a pre-massed 1.5 mL microcentrifuge tube using a fresh microspatula (disposable smartSpatula(R),L anti-static, Cat# Z561762-300EA). The tube was massed again to obtain the exact mass, and the solid dye was resuspended in enough DMSO to create a 5 mM stock solution. To dissolve, the solution was vortexed for 30 seconds, left in the dark at room temperature for 16 hours, and vortexed again. Final dissolved solution was distributed into 15 μ L aliquots in PCR strip tubes (0.2mL Split-Strip 8-Strip Tubes w/Individually Attached Flat Caps, Clear) using an E100 ClipTip p125 Matrix Pipette (Yellow). 5 mM stocks were stored in 1.5 mL microcentrifuge tubes, as well as in 15 μ L aliquots in PCR strip tubes in the dark in -80 °C. Fresh dye stocks were prepared for the entire library every 12-18 months. In addition, 40 μ L of each dye was placed in a LabCyte Echo Qualified 384PP Source Microplate (Cat# PP0200).

Dye library DMSO stock handling

To a LabCyte Echo Qualified 384PP Source Microplate (Cat# PP0200), 40 μ L of each dye library dye was added to create a “mother plate” from which single-use dye screening plates “daughter plates” could be dispensed. At regular intervals, the mother plate was thawed, and 250 nL of each dye was transferred into a single well of a 384-well V-bottom PP Greiner BioOne Microplate (Cat #781280) using a LabCyte Echo. Transfer instructions were created using an in-house web app. The mother and daughter plates were sealed in foil (Excel Scientific PolarSeal Aluminum Microplate Seals, FPS-EK-100), and thawed prior to use.

General protein handling.

All proteins were stored and handled as recommended by each protein source. All proteins were screened in the preferred biochemical buffer for that specific protein, including any additives. Proteins with no preferred buffer were screened in 10 mM HEPES pH 7.20, 200 mM NaCl. See table S2 for information on specific protein sources, buffers, and screening concentrations.

Calculation of structural and biochemical protein properties.

Primary sequences were analyzed using the Bio.SeqUtils.ProtParam python module. The specific functions used for each parameter were: `.molecular_weight()` for **Molecular weight**, `.gravy()` for Grand average of **hydropathicity** index (GRAVY), `.get_amino_acids_percent()` for **primary sequence composition**, and `.secondary_structure_fraction()` for fraction of amino acids found primarily in **helix, turn, and sheet secondary structures**. Code for this analysis was sourced with minor adaptations from “Visualizing and Analyzing Proteins in Python” by Aren Carpenter published online (<https://towardsdatascience.com/visualizing-and-analyzing-proteins-in-python-bd99521ccd>).

The presence of **disordered regions** was calculated using the IUPred2A(36) online server (<https://iupred.elte.hu/>). A multi-FASTA file containing all protein primary sequences was uploaded, and two separate queries were submitted, (i) searching for long regions of disorder, using the settings “Analysis type: IUPred2A long disorder” (>30 residues) with “Context-dependent predictions (default ANCHOR2)”, and (ii) searching for short regions of disorder, using the settings “Analysis type: IUPred2A long disorder” (<30 residues) with “Context-dependent predictions (default ANCHOR2)”. A single “percent predicted disorder”, presented in table S3 was calculated for each protein by dividing the number of residues with a IUPred long disorder probability ≥ 0.5 by the total number of residues in that protein.

Calculation of protein sequence similarity

Primary sequences for all tested proteins (see data S2) were converted to FASTA format using basic string handling functions in R, and then read as sequences using the function `msa::readAAStringSet()`, and a multiple sequence alignment was calculated between all sequences using `msa::msa()`. The resulting multiple sequence alignment was converted to a `seqinr`-compatible class using `msa::msConvert(type = “seqinr::alignment”)`. The resulting object was used to calculate pairwise distance matrices for both sequence similarity and sequence identity using `seqinr::dist.alignment(matrix = “similarity”)` and `seqinr::dist.alignment(matrix = “identity”)`, respectively. The resulting root distance matrix was converted to a similarity matrix (S) by $S = 1 - (\text{root distance})^2$.

Dye screens

Each protein was screened for DSF compatibility with library dyes in the following manner. (A generalizable protocol for dye screening is provided in fig. S4 and movie S1.) For individual protein sources, concentrations, and screening buffers, see table S2.

Resuspension of library dyes A single “daughter plate”, containing 250 nL of 5 mM DMSO stock of each library dye was thawed from -80 °C. After thawing, the daughter plate was visually inspected to ensure no wells were missing dyes. 10 mL of the desired buffer was placed in a buffer reservoir (VWR Catalog #89094-670), and 20 uL of buffer was transferred into each well of the daughter plate using a E100 ClipTip p125 Matrix Pipette (Thermo Fisher #4671040BT, “Yellow”) and associated tips (ClipTip 384 125 Pipette Tips, #94410153) . The resuspended daughter plate was spun-down in a salad spinner for 30 seconds to remove bubbles, and placed on the deck of an Opentrons OT-2 liquid handling robot at room temperature.

Protein dispensing and dye addition. Protein was diluted to 5X screening concentration, and 2 μL of this solution was dispensed into each well of a white, low-volume 384-well qPCR plate (Axygen PCR-284-LC480WNFBC, lot # 23517000) using an E100 ClipTip matrix pipette, and the plate was spun down for 30 seconds to settle droplets to the bottoms of the wells. The plate was placed on the deck of the Opentrons OT-2 liquid handling robot, and the OT-2 was used to transfer 8 μL of resuspended dye from the resuspended dye daughter plate into the protein-containing qPCR plate. Any dyes which remained adhered to the daughter plate were not further mixed or disrupted prior to transferring. Following transfer, each well contained 10 μL of volume, with a typical protein concentration of 1 μM and typical dye concentration of 50 μM . Exact conditions for each screen can be found in table S2.

Heating and measurement. The qPCR plate was sealed with optically clear sealing film (Applied Biosystems MicroAmp Optical Adhesive Film #4311971) and spun for 30 seconds to settle droplets into the bottoms of wells and remove bubbles. The qPCR plate was placed in an Analytik Jena 384G qTower qPCR and heated from 25 °C to 94 °C in 69 increments of 1 °C every 30 seconds with lid heating to 98 °C, and fluorescence measurements were taken at each degree in six channels (“FAM” ex/em = 470/520, “JOE” ex/em = 515 nm/545 nm, “TAMRA” ex/em = 535 nm/580 nm, “ROX” ex/em = 565 nm/605 nm, “Cy5” ex/em = 630 nm/670 nm, and “Cy5.5” ex/em = 660 nm/705 nm). This heating and measuring procedure results in a heating rate of 1 °C per 1 minute, and a total runtime of approximately 70 minutes. Raw fluorescence readings were exported and analyzed using an in-house script. Tmas were extracted using DSFworld (15).

Hit validation. Following hit selection, a custom dye daughter plate was created, including the initially-tested dye concentration, as well as three two-fold serial dilutions for each chosen hit dye. Validation daughter plates were dispensed as described above in “Dye library DMSO stock handling”; serial dilutions of dye stocks were performed during the LabCyte Echo transfer by dispensing the appropriate amounts of either dye stock or DMSO into the well. The above procedure for “Dye screens” was repeated using this validation plate, yielding a four-point dose-response for each dye. Raw fluorescence readings were exported and analyzed using an in-house script.

Dye screening of SARS-CoV2 proteins

Dye screens for all SARS-CoV2 proteins were conducted following the procedure described in Fig. 1G and fig. S3, and outlined above “Dye screens”. See table S2 for protein sources. Proteins were screened in the following buffer, termed “CoV2 DSF buffer”: 50 mM Tris HCl pH 7.5, 150 mM NaCl, 1 mM TCEP, 0.01% Triton X-100. For all experiments, TCEP and Triton were added immediately prior to use. All Tmas were calculated using DSFworld (15).

paDSF of PLPro with GRL0617

GRL0617 was obtained from TOCRIS (Cat # 7280) and resuspended in DMSO to a stock concentration of 50 mM. In 10 mL of CoV2 DSF buffer (50 mM Tris HCl pH 7.5, 150 mM NaCl, 1 mM TCEP, 0.01% Triton X-100) a 2.5X GRL0617 stock was prepared by diluting 5 µL of 50 mM GRL0617 stock in 995 µL buffer. Three 2-fold serial dilutions were prepared from this stock. 80 µL of 5X PLPro (5 µM) was prepared in CoV2 DSF buffer from a 32 µM stock solution. 1 mL of 2.5X dye was prepared for each dye. These components were combined and mixed well to yield the final tested conditions: 5 µM dye; 100, 10, and 0 µM GRL0617; 1 µM PLPro, or for no-protein controls, buffer alone in place of protein. 10 µL of solution was transferred to each well of the 384-well plate, and DSF was performed as described in “*Heating and measurement*” above.

paDSF of nsp3-mac1 with ADP ribose

A daughter plate containing 250 nL of each dye was dispensed using the LabCyte Echo as described above in “*Dye library DMSO stock handling*”. ADP ribose was obtained from Sigma Aldrich (A0752) and resuspended in DMSO to a stock concentration of 10 mM. To the CoV2 DSF buffer, EDTA was added to a final concentration of 1 mM. In this buffer, seven serial dilutions of ADP ribose were prepared at 1.25X final concentration. The top concentration was prepared by

adding 74 μL of 10 mM ADP ribose to CoV2 buffer. To the dye daughter plate, 20 μL of the appropriate concentration of ADP ribose was added to each well. A 5X protein solution was prepared by adding 13.3 μL of 300 μM purified nsp3-mac1 to 386.7 μL of buffer for a final concentration of 10 μM , and 2 μL of this protein solution was added to each well of a white qPCR plate. To the white qPCR plate, 8 μL of dye-ADPribose solution was added. DSF was performed as described above in “Heating and Measurement.” All dyes were tested at a final concentration of 5 μM with the exception of L094, which was tested at 25 μM . ADP ribose was tested at 300, 150, 75, 37.5, 18.8, 9.3, 4.6 and 0 μM . Protein was tested at a final concentration of 2 μM , or 0 μM for the no-protein controls.

paDSF of nsp3-mac1 with screening compounds

To the CoV2 DSF buffer, EDTA was added to a final concentration of 1 mM, DMSO was added to a final concentration of 3%, and SYPRO Orange was added to a final concentration of 5X. In this buffer, a 2X protein solution of nsp3-mac1 was created by dilution of 270 μL of 37 μL nsp3-mac1 stock solution into a 2 mL final volume. Buffer and stock solution of a purine small molecule library was added to the DSF plate, for a final compound concentration of 20 μM . DSF was performed as described above in “heating and measurement.” For head-to-head comparison between dye T004 and SYPRO, the above procedure was repeated, using either 5X SYPRO Orange or 5 μM T004, and 192 compounds from the ZINC catalog, as described in (23).

Mass photometry of OGT

Slide cleaning and gasket preparation. Prior to experimental set up, glass slides (High Precision Deckgläser Microscope cover glass, 24 x 50 mm, 170 +/- 5 μm , No 1.5H) were cleaned by three repetitions of: submersion and swirling in a beaker of isopropyl alcohol, air dry, submersion and swirling in a beaker of miliQ water, air dry, followed by one final round of drying by pressurized air passed through a pasteur pipette tip. Mass photometry gaskets (Grace BioLabs Ref 103250 Culture Well Reuseable Gaskets, CW-50R-1.0, 50-3mm DIA x 1 mm depth, 3-10 μL) were prepared by cutting into 2 x 4 well strips.

Sample preparation. 50 mL of filtered OGT buffer was prepared by combining 2.5 mL Tissue-Culture Grade 1M HEPES (Gibco 1 M HEPES Buffer Solution, Ref 15630-106) with 3 mL 5M NaCl and 44.5 μL MiliQ water, both of which had been spin filtered immediately prior to use (Ultrafree PVDF 0.1 μm , Ref R7KA45465). To a 4 mL aliquot of this buffer, 4 μL of DTT was added, and the buffer was spin-filtered twice. In this final buffer, each of full-length OGT, OGT TPR domain, and OGT catalytic domain were diluted to 2 μM , and the ligand L4 was diluted from a 10 mM DMSO stock to a final concentration of 20 μM . All protein and ligand solutions were once again spin filtered. From the 20 μM L4 solution, six two-fold serial dilutions were made in the filtered buffer. To each of three L4 concentrations or L4-free buffer, a 1:1 mixture was made between ligand and each of the protein solutions, or protein-free buffer, for final L4 concentrations of 0, 0.6, and 2.5 μM , and final protein concentrations of 0 or 1 μM . The solutions were left to incubate for 2 hours at room temperature to equilibrate.

Data collection. Mass photometry data was collected on a Refeyn OneMP Mass Photometer as follows. Instrument was focused and calibrated to BSA (2 mg/mL BSA solution, Bradford Assay Kit) and thyroglobulin (Cytiva MW Native Marker Kit, Thyroglobulin standard Mr 669000) mass standards. Each sample was measured well as follows: into a fresh gasket well, 10 μL of filtered buffer was pipetted, followed by 2 μL of 2 μM protein-ligand mixture, followed by gentle mixing by pipetting up and down, for a final volume of 12 μL . Final conditions measured

were 0, 0.6, and 2.5 μM L4 in the presence of full length OGT, catalytic domain alone, TPR domain alone, or buffer (as a control for protein-independent ligand aggregation).

paDSF of OGT with dye TW408

OGT buffer: To 10 mL of OGT buffer, 10 μL of 1 M DTT was added. Following thorough mixing, 100 μL of a 5 mM DMSO stock of TW408 was added to this buffer, to yield the final OGT TW408 buffer (final buffer: 50 mM HEPES pH 7.4, 300 mM NaCl, 1 mM DTT, 50 μM TW408, 1% DMSO). All following solutions were prepared in this buffer.

2X protein solutions: 140 μL of 4X solutions of purified TPR and catalytic domains were prepared by dilution of concentrated protein stocks to a final 4X concentration of 4 μM . These 4X protein solutions were used to prepare 140 μL of 2X solutions of either (1) TPR alone (70 μL 4 μM TPR + 70 μL buffer), (2) catalytic alone (70 μL 4 μM catalytic + 70 μL buffer), or (3) 1:1 mix of catalytic and TPR domains (70 μL 4 μM TPR + 70 μL 4 μM catalytic). A 2X solution of OGT was then prepared by first creating a 1.3X concentrated TW408 buffer, to offset the volume displaced by the OGT stock solution, and then combining 210 μL of this buffer with 70 μL of 8 μM OGT stock.

Compound solutions: A concentrated solution of the highest tested concentration for each compound was created in OGT buffer. From this solution, six 2-fold serial dilutions were made in OGT buffer.

UDP GlcNAc: Ultra Pure UDP-GlcNAc (Promega Cat # V7071) was prepared in buffer in seven 2-fold serial dilutions, from an highest concentration of 100 μM .

L4: 1 μL of 50 mM DMSO stock of L4 was diluted in 9 μL of DMSO and 990 μL of TW408 OGT buffer, to yield a final 5X concentrated L4 solution of 50 μM . From this 50 μM solution, seven 2-fold serial dilutions were created in OGT buffer.

Combination of protein and ligand, and DSF: In PCR strip tubes, appropriate volumes of concentrated protein, compound, and buffer were combined to a final 1X volume of 36 μL , and the resulting solution was mixed well by pipetting up and down. 10 μL of each sample was then transferred in triplicate into a white, low-volume 384-well qPCR plate (Axygen PCR-284-LC480WNFBC, lot # 23517000) using an E100 ClipTip matrix pipette. The plate was spun down for 30 seconds to settle droplets to the bottoms of the wells. and DSF measurement was performed as described in ‘*Heating and measurement*’ (above).

Data analysis.

Transition ratios. Ratios were calculated from paDSF experiments in which all constructs were present at 1 μM , and TW408 was present at 50 μM . paDSF detection ratios between isolated catalytic and TPR domains were calculated between the mean maximum RFU intensities reached, and correspond to the peak of the unfolding transition. paDSF detection ratios between transitions 1 and 2 in OGT (corresponding to $T_{\text{ma_OGT-1}}$ and $T_{\text{ma_OGT-2}}$) were calculated between the mean maximum RFU intensities reached by the two respective sigmoidal components of the fitted model from which $T_{\text{ma_OGT-1}}$ and $T_{\text{ma_OGT-2}}$ were calculated.

Statistical tests. Tests for significance in concentration-dependent thermal shifts with UDPGlcNAc or L4 were calculated using a Spearman correlation between ligand concentration and T_{ma} , using `stats::cor.test` function in R, with `method = "spearman"`. Tests for significance in differences between two groups (e.g. $T_{\text{ma_OGT-1}}$ vs $T_{\text{ma_cat}}$, $T_{\text{ma_OGT-1}}$ vs $T_{\text{ma_TPR}}$) were performed using Two-tailed Students T-Test, using `stats::t.test()` function in R. In all cases, correlations were considered significant if $p < 0.01$.

Variable temperature UDP-Glo Enzymatic activity assay of OGT

Material preparation. Buffer was prepared fresh by adding 4 μL of 1 M DTT to 4 mL of OGT buffer (final buffer: 50 mM HEPES pH 7.4, 300 mM NaCl, 1 mM DTT). Three 4X-concentrated L4 solutions were prepared by adding 1 μL 10 mM L4 in DMSO to 999 μL of buffer. Following thorough mixing, two 150 μL two-fold serial dilutions were prepared, yielding at least 150 μL each of 20, 10, 5, and 0 μM L4 solutions. A 2X-concentrated protein solution was prepared by adding 50 μL of 8 μM purified OGT to 150 μL buffer, yielding 200 μL of 2 μM OGT. The 4X-concentrated sugar-substrate mix was prepared by adding 2 μL of 100 mM Ultra Pure UDP-GlcNAc (Promega Cat # V7071) to 998 μL buffer, followed by the addition of 80 μL of 1 mM substrate peptide (Anaspec Cat. # 63726, Lot#1958507, >95%, sequence: KKKYPPGGSTPVSSANMM) to 720 μL of this UDP-GlcNAc buffer, yielding 800 μL of buffer containing 100 μM substrate and 200 μM UDP-GlcNAc.

Enzymatic activity at variable temperatures. Enzymatic reaction mixtures were prepared at each of the four L4 concentrations by combining 48 μL 4X OGT, 96 μL 2X sugar-substrate, and 48 μL 4X L4. Final conditions were: 500 nM OGT, 100 μM UDP-GlcNAc, 50 μM substrate peptide, and either 5, 2.5, 1.25, or 0 μM L4. As a negative control, the identical reactions were set up using buffer in place of OGT, yielding final conditions of 0 nM OGT, 100 μM UDP-GlcNAc, 50 μM substrate peptide, and either 5, 2.5, 1.25, or 0 μM L4. To monitor catalytic activity across a range of temperatures, all reactions were immediately divided into 22 μL aliquots across 8 PCR strip tubes, and these tubes were placed in a PCR thermocycler, and incubated across a gradient of 30 - 50 $^{\circ}\text{C}$ with a heated lid for 140 minutes.

Quenching, detection, and data collection of enzymatic activity. During incubation, the Promega UDP-Glo reagent Promega UDP-Glycosyltransferase assay (Cat # V6961) was prepared. The “UDP-Glo working solution” was prepared by addition of 1 μL “UDP-Glo enzyme” to 74 “Enzyme Dilution Buffer”. The “UDP detection reagent” was prepared by adding 10 μL of this “UDP-Glo working solution” to 990 μL “Nucleotide Detection Reagent”. Following the 140 min incubation, all samples were rapidly cooled to 12 $^{\circ}\text{C}$ in the thermocycler, and 5 μL of each solution was immediately transferred into a well of a white, low-volume, hard-sided 384-well plate (Greiner BioOne #784075, Lot 1209043), and 5 μL of the prepared “UDP detection reagent” was added to all wells, simultaneously quenching OGT activity and beginning the UDP detection reaction. The plate was shaken gently by hand to mix, centrifuged for 1 min at 1400 rcf to remove bubbles, and incubated at room temperature in the dark for 65 minutes. Following incubation, resulting luminescence was recorded on a Tecan Spark plate reader.

Calculation of dye physical chemical properties

In python, the SMILES string for each dye was converted to a mol object using Chem.MolFromSmiles. For each mol object, the number of hydrogen bond donors, hydrogen bond acceptors, and heteroatoms were calculated using the rdkit.Chem.Lipinski module; cLogP and Molecular Refractivity were calculated using rdkit.Chem.Crippen module, and molecular weight and Topological Polar Surface Area were calculated using the MolWt and TPSA functions from the rdkit.Chem.Descriptors module respectively. In R, the Pearson correlation was calculated between the dye hit rate and each calculated property using cor.test.

Calculation of ECFP4 Tanimoto coefficients and clustering

SMILES strings of all dyes were canonicalized using Chem.CannonSmiles from the 'RDKit' package in python 3. From canonicalized SMILES, pairwise ECFP4 Tanimoto Coefficient was calculated using the Morgan algorithm implemented in the 'RDKit' package in python 3. The resulting pairwise similarity matrix was converted to dissimilarities by subtracting each value from 1. This dissimilarity matrix was passed to hierarchical cluster analysis using Lance-Williams dissimilarity and 'complete' linkage clustering method via the 'hclust' function implemented in R package in the 'stats' package. Similarities were calculated within four Chembridge Screening Libraries, a 50K compound "diversity" library (DiverSET), and three "focused" libraries "GPCR", "Ioncore", and "NHRCore", in the same manner.

Scaffold extraction. To define chemical scaffolds, the same pairwise matrix was passed to Butina clustering, implemented using rdkit.ML.Cluster.Butina.ClusterData with cutoff = 0.7.

Calculation of dye "hit rate". Dye "hit rate" was calculated by dividing the number of proteins for which a dye was called a "hit" from the number of proteins against which it was tested. Pairwise dye hit rate similarity was calculated as the absolute value of the hit rate between two dyes, and a complete pairwise similarity matrix was calculated for the library. For example, if a dye was screened against 70 proteins, and was a "hit" for 45 of them, the hit rate for that dye would be 45/70, or 64.3 percent. Dye "sensitivity hit rate" was calculated in the same way, but using the number of proteins for which a dye was called "sensitive" instead of "hit."

Calculation of dye hit profile similarities

The pairwise similarity between paDSF dye screen results for different proteins was calculated using the Jaccard coefficient, $J(A,B) = |A \cap B| / |A \cup B|$, where A and B are each the set of outcomes (hit, sensitive, or none) for proteins for which a different dye was tested. Specifically, the union was defined as the number of dyes which received the same assignment between the two proteins (both hit, both sensitive, or both none). The intersection was defined as the number of dyes which were tested against both proteins, regardless of the assignment. Higher jaccard coefficients indicate more similar dye screen results: proteins which received different assignments for every dye would have $J = 0$, while proteins which received the same assignment for every dyes would have $J = 1$. Dyes which never received an assignment of "hit" or "sensitive" for any tested protein (103 dyes) were not included in this calculation. For each protein pair, the Jaccard coefficient was calculated using only the dyes which were tested against both proteins. The final dataset includes 2,211 unique comparisons made between the 67 screened proteins.

Calculation of correlation between protein sequence and dye fingerprint

Proteins which contained multiple distinct chains (eIF2B, nsp10/nsp16, nsp10/nsp14, proteasome, and nucleosome) were removed prior to calculation, as it was not possible to distinguish the contribution of each chain to the dye screen results. For the remaining set of all unique protein pairs (67 proteins, 2,211 unique pairs), the Spearman correlation between dye screen results (see above) and primary sequence similarity (see above) using the `base::cor.test(method = "spearman")` function in R.

References.

36. B. Mészáros, G. Erdos, Z. Dosztányi, IUPred2A: context-dependent prediction of protein disorder as a function of redox state and protein binding. *Nucleic Acids Res.* **46**, W329–W337 (2018).
37. P. Rajagopal, E. Tse, A. J. Borst, S. P. Delbecq, L. Shi, D. R. Southworth, R. E. Klevit, A conserved histidine modulates HSPB5 structure to trigger chaperone activity in response to stress-related acidosis. *Elife.* **4** (2015), doi:10.7554/eLife.07304.
38. C. L. Gustafson, N. C. Parsley, H. Asimgil, H.-W. Lee, C. Ahlback, A. K. Michael, H. Xu, O. L. Williams, T. L. Davis, A. C. Liu, C. L. Partch, A Slow Conformational Switch in the BMAL1 Transactivation Domain Modulates Circadian Rhythms. *Mol. Cell.* **66**, 447–457.e7 (2017).
39. M. Ravalin, P. Theofilas, K. Basu, K. A. Opoku-Nsiah, V. A. Assimon, D. Medina-Cleghorn, Y.-F. Chen, M. F. Bohn, M. Arkin, L. T. Grinberg, C. S. Craik, J. E. Gestwicki, Specificity for latent C termini links the E3 ubiquitin ligase CHIP to caspases. *Nat. Chem. Biol.* **15**, 786–794 (2019).
40. A. K. Michael, J. L. Fribourgh, Y. Chelliah, C. R. Sandate, G. L. Hura, D. Schneidman-Duhovny, S. M. Tripathi, J. S. Takahashi, C. L. Partch, Formation of a repressive complex in the mammalian circadian clock is mediated by the secondary pocket of CRY1. *Proc. Natl. Acad. Sci. U. S. A.* **114**, 1560–1565 (2017).
41. J. C. Tsai, L. E. Miller-Vedam, A. A. Anand, P. Jaishankar, H. C. Nguyen, A. R. Renslo, A. Frost, P. Walter, Structure of the nucleotide exchange factor eIF2B reveals mechanism of memory-enhancing molecule. *Science.* **359** (2018), doi:10.1126/science.aag0939.
42. D. R. Wassarman, K. Bankapalli, L. J. Pallanck, K. M. Shokat, Tissue-restricted inhibition of mTOR using chemical genetics. *Proc. Natl. Acad. Sci. U. S. A.* **119**, e2204083119 (2022).
43. E. Thompson, thesis, University of California, San Francisco, Ann Arbor, United States (2020).
44. S. Militi, E. S. Maywood, C. R. Sandate, J. E. Chesham, A. R. Barnard, M. J. Parsons, J. L. Vibert, G. M. Joynson, C. L. Partch, M. H. Hastings, P. M. Nolan, Early doors (Edo) mutant mouse reveals the importance of period 2 (PER2) PAS domain structure for circadian

Caption for Movie S1

Movie S1. Visual protocol for dye screen. Annotated video of the primary liquid handling steps of a dye screen.

Captions for Data S1 to S4

Data S1. Chemical structures and common names of Aurora dyes. See Data S3 for more information, includes specific dye commercial sources, CAS numbers, and MSDS where available.


Data S2. Full dye screen results for all tested proteins. These results can also be explored interactively at https://padsfdyes.shinyapps.io/Exp1243_heatmap_cache/.

Data S3. Primary sequences of protein panel constructs. See Table S1 for more information on sources and biochemical properties of proteins.


Data S4. Extended information on Aurora dyes. A spreadsheet containing, for each dye in the library: Aurora-library name (e.g. L095), common name (where applicable), CAS, source, catalog number, MSDS, and SMILES. the 312 dyes in the Aurora library.

a. Multiple improved dye-protein pairs provides chemical and optical flexibility

Protein: Nuclear-hormone receptor 23



b. Improved dye-protein pairs resolve multiple archetypes of SYPRO incompatibility



c. Additional representative examples of paDSF dye-protein hits




Fig. S1. Representative examples of paDSF screening results. For full screening results for all proteins, see Data S2. **(A)** Results of a dye screen against the SYPRO-incompatible protein nuclear-hormone receptor 23. Far left panel: representative unsuccessful DSF conducted with SYPRO Orange. Right panels: examples of successful paDSF data, collected with six different paDSF dyes. Plots show dyes at four different concentrations, demonstrating a dose-responsive relationship between RFU and dye concentration. **(B)** Representative examples of resolution of four archetypes of SYPRO Orange incompatibility by paDSF dyes. **(C)** Examples of successful paDSF data collected for four additional proteins. As in panel A, dyes are plots show dyes at four different concentrations, demonstrating a dose-responsive relationship between RFU and dye concentration.

***T_m* apparent of tested proteins with paDSF dyes**





Fig. S2. Apparent melting temperatures for paDSF dyes versus SYPRO Orange. T_{mas} reported by compatible dyes for all tested proteins (single point: single compatible dye), compared to SYPRO Orange values. Error bars represent +/- standard error.

See end of document for **Table S1. Chemical structures and common names for all Aurora library dyes.** Chemical structures, Aurora library name, and common name (where applicable) for all 312 tested dyes in the Aurora library. Unless otherwise noted, dyes are available either commercially, or through the North Carolina State Max Weaver Dye Library. See Data S1 for extended information on dyes and sources.

See end of document for **Fig. S3. Apparent melting temperatures for paDSF dyes versus SYPRO Orange.** T_{mas} reported by compatible dyes for all tested proteins (single point: single compatible dye), compared to SYPRO Orange values. Error bars represent +/- standard error.

a. nsp3 mac1 with ADP ribose



b. PLPro with GRL0617





Fig. S4. Apparent thermal shifts with paDSF dyes for SARS-CoV2 proteins. (A) Thermal upshift in the presence of ADP ribose with nsp3 mac1, with 6 paDSF dyes and SYPRO orange. Error bars represent +/- standard deviation. (B) Thermal upshift in the presence of GRL0617 with PLPro, with 3 paDSF dyes.

a. Compound incompatibility of SYPRO Orange with 27 of 320 compounds in pyridine pilot screen



b. Restoration of compound compatibility via use of dye T004 in place of SYPRO Orange




Fig. S5. Reduction in artifactual dye activation by screening compounds for nsp3 mac1. (A) Protein-independent activation of SYPRO Orange observed for compounds which produced uninterpretable raw DSF data from a 320-compound pyridine pilot screen using SYPRO Orange. (B) Comparison of DSF data collected using SYPRO Orange and dye T004, demonstrating improved compound compatibility of T004 over SYPRO Orange.




Fig. S6. Raw and analyzed paDSF with OGT, catalytic domains and TPR domain. paDSF results from the 10 hit dyes for OGT, with raw RFU data displayed alongside fitted models for (A) OGT alone, fitted to DSFworld Model 4 (two sigmoids, with initial RFU (B) catalytic domain alone, fitted to DSFworld Model 2 (one sigmoids, with initial RFU (C) TPR domain alone, fitted to DSFworld Model 2 (one sigmoids, with initial RFU (D) overlay of raw data and final fits for all three constructs. (E) T_{ma} s for the three constructs, extracted using the fits displayed in panels A-D. Data are presented as the mean of three technical replicates +/- standard deviation. No significant (ns) difference was found between T_{ma} OGT-1 and T_{ma} cat, nor T_{ma} OGT-2 and T_{ma} TPR. Significance determined using two-tailed t-test, and a significance threshold of $p < 0.01$.






Fig. S7. Diagram of OGT constructs used in this study. (A) Structural model of full-length OGT, composed from crystal structures of TPR (PDB 1W3) and catalytic (PDB 5C1D) domains. Cartoons of the (B) full-length OGT construct, which includes the full TPR and catalytic domains, (C) catalytic domain construct, which includes the full catalytic domain and TPR repeats 9 - 13.5, and (D) TPR domain construct, which includes TPR repeats 1-10.

a. Scheme of UDP-GlcNAc interaction with tested constructs



b. Raw paDSF with dye TW408, for UDP-GlcNAc titration



c. Comparison of thermal shifts in OGTs transitions and isolated domains






Fig. S8. Thermal shifts with UDP-GlcNAc from TW408 paDSF. paDSF results using dye TW408 to monitor titration of UDP-GlcNAc with OGT, catalytic domain, or TPR domain. (A) Schematic of interaction of UDP-GlcNAc with full length and isolated catalytic domain, but not TPR domain, of OGT. (B) Raw RFU for each construct in the presence of 0 - 50 μM UltraPure UDP-GlcNAc. (C) Comparison of resulting T_{mas} . Left panel: $T_{\text{ma OGT1}}$ and $T_{\text{ma cat}}$; right panel: $T_{\text{ma OGT2}}$ and $T_{\text{ma TPR}}$. Both $T_{\text{ma OGT1}}$ and $T_{\text{ma cat}}$ show significant thermal upshift with UDP-GlcNAc ($T_{\text{ma OGT1}}$ spearman estimate = 0.92, $p = 1.13\text{e-}10$; $T_{\text{ma cat}}$ spearman estimate = 0.94, $p = 1.48\text{e-}12$); while neither $T_{\text{ma OGT2}}$ and $T_{\text{ma TPR}}$ do ($T_{\text{ma OGT2}}$ spearman estimate = 0.095, $p = 0.66$; $T_{\text{ma TPR}}$ spearman estimate = 0.005, $p = 0.98$). Data are presented as the mean of three technical replicates +/- standard deviation.

a. Scheme of L4 interaction with tested constructs



b. Raw paDSF with dye TW408, for L4 titration



c. Comparison of thermal shifts in OGTs transitions and isolated domains





Fig. S9. Thermal shifts with L4 from TW408 paDSF. paDSF results using dye TW408 to monitor titration of L4 with OGT, catalytic domain, or TPR domain. (A) Schematic of L4 with the four tested conditions: full OGT, catalytic domain, TPR domain, and a 1:1 mix of catalytic and TPR domain. All constructs were tested at 1 μM final; 1:1 mix had 1 μM catalytic domain and 1 mM TPR domain. (B) Raw RFU for each construct in the presence of 0 - 10 μM L4. (C) Comparison of resulting T_{ma} s. Left panel: $T_{\text{ma OGT1}}$, $T_{\text{ma cat:TPR 1}}$, and $T_{\text{ma cat}}$; right panel: $T_{\text{ma OGT2}}$, $T_{\text{ma cat:TPR 2}}$, and $T_{\text{ma TPR}}$. Significant thermal upshifts with L4 treatment observed with $T_{\text{ma OGT2}}$, $T_{\text{ma cat:TPR 2}}$, and $T_{\text{ma TPR}}$ ($T_{\text{ma OGT2}}$ spearman estimate = 0.95, $p = 2.3\text{e-}12$; $T_{\text{ma cat:TPR 2}}$ spearman estimate = 0.97, $p = 4.7\text{e-}15$; $T_{\text{ma TPR}}$ spearman estimate = 0.96, $p = 3.6\text{e-}14$). Significant thermal downshift with L4 observed in $T_{\text{ma OGT1}}$ (data in Fig. 3D, included here for reference, $T_{\text{ma OGT1}}$ spearman estimate = -0.92, $p = 6.2\text{e-}11$). No dose-responsive thermal shift with L4 observed in $T_{\text{ma cat}}$ or $T_{\text{ma cat:TPR 1}}$ ($T_{\text{ma cat}}$ spearman estimate = -0.27, $p = 0.2$; $T_{\text{ma cat:TPR 1}}$ spearman estimate = -0.34, $p = 0.09$). Data are presented as the mean of three technical replicates \pm standard deviation.

a. Extracted molecular weights



b. Histograms of measured molecular weights





Fig. S10. OGT oligomerization state and molecular mass not are affected by L4. (A) Molecular weights of full length OGT, catalytic domain, and TPR domain, extracted from mass photometry data. Dashed lines: expected molecular weights for each species. OGT: 115 kDa, Catalytic domain; 82 kDa, TPR domain: 43, 86, 129 kDa (monomer, dimer, trimer). Error bars: +/- standard error in mass from mass photometry. (B) Histograms of counts of MW from mass photometry.

a L4- and Temperature-dependent OGT activity



b L4-dependent OGT activity at 41 °C




Fig. S11. L4 treatment increases temperature sensitivity of OGT catalytic activity. OGT catalytic activity with a model substrate peptide (KKKYPGGSTPVSSANMM), with enzymatic activity allowed to proceed at 29, 32, 35, 38, 41, 44, 47, or 50 °C. Activity was measured as RLU using Promega UDP-Glo Assay, with higher RLU values corresponding to higher activity. **(A)** Comparison of temperature-dependent loss of enzymatic activity (RLU) in the presence of increasing concentrations of L4. Solid lines and dots: OGT-containing samples; Dashed lines and Xs: protein-free control. **(B)** L4-dependent increase in temperature sensitivity of catalytic activity at 41 °C. In all panels, data are presented as the mean of three technical replicates +/- standard deviation.




Fig. S12. Chemical diversity of paDSF dyes versus ChemBridge DIVERSet library.



a dye cLogP vs DSF hit-rate**b** Distribution of cLogP values in Aurora

Fig. S13. cLogP of dyes and relationship to paDSF performance. (A) Comparison of cLogP, a calculated predictor of water solubility, to the paDSF performance of that dye. Filled circles: percent of tested proteins for which a dye was called a “hit”. Xs: percent of tested proteins for which a dye was called “sensitive”. No significant correlation was found between hit rate and cLogP (Pearson p value > 0.05, see panel for exact values) (B) Distribution of cLogP values for the full Aurora library (312 dyes), for reference.




Fig. S14. Presence and absence of known fluorogenic cores in paDSF dyes. Individual points represent a single dye, and dyes are arranged along the axis in order of increasing hit rate. Data are shown for the full Aurora library (312 dyes). Colored points indicate dyes which contain a given substructure, grey points indicate dyes lacking that substructure. A single substructure is presented in each panel, and panels are shown only for substructures present in at least one dye. All points left of the grey vertical line, in the area labeled “Never hit” were never a hit with any tested protein. Some dyes contain multiple substructures.













Fluorogenic motif		Dyes with motif			Assignments from dyes with this motif		
Name	Structure	Active	Inactive	Total	Hit	None	Percent of all hits
<i>Summary</i>							
No queried motif		44	158	202	207	9,330	23
Any queried motif	See below	63	48	111	704	4,581	77
<i>By substructure</i>							
BODIPY		2	7	9	2	412	0.06
6DMN		0	0	0	0		0
Prodan		0	0	0	0		0
Anthradan		0	0	0	0		0
Dansyl		0	0	0	0		0
Coumarin		6	22	28	7	1,333	0.2
Vinyl julolidine		3	4	7	41	268	1.4
Cyanine		25	8	33	1,463	8,121	49
Hemicyanine		21	2	23	1,368	2,694	46
Streptocyanine		3	2	5	47	429	47
Merocyanine		7	1	8	73	848	2.4
Oxonol		0	0	0	0		0

Table S3. Presence and activity of queried fluorogenic scaffolds in the full Aurora dye library.

Distribution of protein hit rates




Fig. S15. Distribution of protein hit rates of individual dyes. Density plot: the distribution of hit rates across 107 dyes which were assigned “hit” with at least one tested protein (Aurora-active). Points: single points represent the protein hit rates for individual dyes. For dyes which were assigned “hit” with at least one tested protein, the protein selectivities ranged from detection of 1 protein out of 63 tested at the most selective, to three dyes which detected over 50 percent of tested proteins at the most promiscuous. There were 206 dyes which showed no activity with any tested protein (maximum of 64 different proteins). These 206 dyes were not included in the plot or density calculation.

Distribution of fingerprint jaccard




Fig. S16. Dye fingerprints are unique to each protein. Density plot of the pairwise Jaccard indices, which quantify the extent of overlaps in the assignments given to each dye across two different proteins. Identical dye fingerprints would have a jaccard value of 1. Box plot contains standard summary information (lower whisker = smallest observation greater than or equal to lower hinge - 1.5 * IQR, lower hinge = 25% quantile, median = 50% quantile, upper hinge = 75% quantile, upper whisker = largest observation less than or equal to upper hinge + 1.5 * IQR; outlier points not displayed). Jittered points: individual jaccard indices (2,211 points representing the 2,211 unique pairs between 67 proteins).

Distribution of promiscuous and selective dyes as hits for each protein

Assignment for given screen ■ hit ■ sensitive



Dye hit rate across all proteins (percent)

Fig. S17. Distribution of promiscuous and selective dyes among the hit and sensitive dyes assigned for each protein. Individual histograms displaying the number of hit dyes for a given protein, binned by the global hit rates of the specific hit dyes. Total hit rates are divided into ten evenly-spaced bins, containing hit rates of up to 10, 20, 30, 40, 50, 60, 70, 80, 90, or 100% of all tested proteins. Green bars: for “hit” assignment. Grey bars: for “sensitive” assignment. The presence of both selective and promiscuous dye hits in an individual dye screen appears as the presence of hits in both lower hit rate (e.g. < 20%) and higher hit rate (e.g. >40%) bins. Results for individual proteins are presented in individual panels, with the protein name indicated in the panel title.



























































Fig. S18. Weak positive correlation between protein sequence similarity and dye fingerprint. Comparison of the relationship between primary sequence similarity and dye screen outcomes shows a weak positive correlation: 0.19, p-value $3.9e-20$. Data shown includes 2,211 unique pairings between 67 proteins (72 proteins tested total - 5 proteins tested only as hetero-complexes; see Supplementary Methods). To provide context, the number of dyes screened in common between two proteins is indicated in point size and color. Smaller, darker points: fewer dyes screened in common. Larger, lighter points: more dyes screened in common. See figure key for exact mappings. Calculations of the Jaccard Index (y axis values) and pairwise sequence similarity (x axis values) are described in Supplementary Methods.












See caption for complete list of dyes in Aurora concise























Class Amyloid Dye A001 Common name BF-188	
	<i>In Aurora concise</i>
Class Amyloid Dye A002 Common name PBB3	
	<i>In Aurora concise</i>
Class Amyloid Dye A003 Common name FSB	
	<i>In Aurora concise</i>
Class Amyloid Dye A004 Common name BSB	
Class Amyloid Dye A005 Common name K114	
Class Amyloid Dye A006 Common name Resorufin B.E.	
Class Amyloid Dye A007 Common name Resorufin E.E.	
Class Amyloid Dye A008 Common name Resorufin	
Class Amyloid Dye A009 Common name Congo Red	
	<i>In Aurora concise</i>
Class Amyloid Dye A010 Common name NIAD4	
	<i>In Aurora concise</i>
Class Amyloid Dye A011 Common name Thiazine Red	























Class Amyloid Dye A012 Common name Nile Red	
Class Amyloid Dye A013 Common name LDS 698	
	<i>In Aurora concise</i>
Class Amyloid Dye A014 Common name LDS 722	
	<i>In Aurora concise</i>
Class Amyloid Dye A015 Common name LDS 730	
	<i>In Aurora concise</i>
Class Amyloid Dye A016 Common name PBB5 (aka Amb615160)	
	<i>In Aurora concise</i>
Class Amyloid Dye A017 Common name CRANAD-2	
Class Amyloid Dye A018 Common name THK 265	
Class Amyloid Dye A019 Common name BAP1	
Class Amyloid Dye A020 Common name ThT	
	<i>In Aurora concise</i>
Class Amyloid Dye A021 Common name DCM 6490	
	<i>In Aurora concise</i>
Class Amyloid Dye A022 Common name LDS 798	
	<i>In Aurora concise</i>

Class Amyloid Dye A023 Common name Amb 1861079	 <i>In Aurora concise</i>	Class Common organic Dye C009 Common name Rose Bengal	
Class Amyloid Dye A024 Common name Amb 2442356	 <i>In Aurora concise</i>	Class Common organic Dye C010 Common name Victoria Blue B	 <i>In Aurora concise</i>
Class Amyloid Dye A025 Common name Amb 1484600		Class Common organic Dye C011 Common name Azophloxine	
Class Common organic Dye C001 Common name Meldola's Blue	 <i>In Aurora concise</i>	Class Common organic Dye C012 Common name Orange G	
Class Common organic Dye C002 Common name Xylene Cyanol FF		Class Common organic Dye C013 Common name Chrysoidine G	
Class Common organic Dye C003 Common name Victoria Blue R	 <i>In Aurora concise</i>	Class Common organic Dye C014 Common name Bismark Brown Y	
Class Common organic Dye C004 Common name Acid Yellow 17		Class Common organic Dye C015 Common name Janus Green B	
Class Common organic Dye C005 Common name Sudan II		Class Common organic Dye C016 Common name Martius Yellow	
Class Common organic Dye C006 Common name Naphthol Yellow S		Class Common organic Dye C017 Common name Tropaeolin O	
Class Common organic Dye C007 Common name Naphthol Green B		Class Common organic Dye C018 Common name DCVJ	
Class Common organic Dye C008 Common name Fuchsin Acid		Class Common organic Dye C019 Common name Titan Yellow	























Class Common organic Dye C020 Common name 1,8,-ANS	
Class Common organic Dye C021 Common name Amaranth Red	
Class Common organic Dye C022 Common name CPM	
Class Laser Dye L001 Common name LD390	
Class Laser Dye L002 Common name POPOP	
Class LaserL003Bis-MSB	
Class Laser Dye L004 Common name LD423	
Class Laser Dye L005 Common name LD425	
Class Laser Dye L006 Common name Coumarin440	
Class Laser Dye L007 Common name Coumarin445	
Class Laser Dye L008 Common name Coumarin450	























Class Laser Dye L009 Common name Coumarin456	
Class Laser Dye L010 Common name Coumarin460	
Class Laser Dye L011 Common name Coumarin461	
Class Laser Dye L012 Common name LD466	
Class Laser Dye L013 Common name LD473	
Class Laser Dye L014 Common name Coumarin478	
Class Laser Dye L015 Common name Coumarin480	
Class Laser Dye L016 Common name Coumarin481	
Class Laser Dye L017 Common name Coumarin485	
Class Laser Dye L018 Common name LD489	
Class Laser Dye L019 Common name Coumarin490	












Class Laser Dye L020 Common name LD490		Class Laser Dye L031 Common name Coumarin535	
Class Laser Dye L021 Common name Coumarin498		Class Laser Dye L032 Common name Coumarin540	
Class Laser Dye L022 Common name Coumarin500		Class Laser Dye L033 Common name Coumarin540A	
Class Laser Dye L023 Common name Coumarin503		Class Laser Dye L034 Common name Coumarin545	
Class Laser Dye L024 Common name Coumarin504		Class Laser Dye L035 Common name Pyromethene546	
Class Laser Dye L025 Common name Coumarin510		Class Laser Dye L036 Common name Pyromethene556	
Class Laser Dye L026 Common name Coumarin515		Class Laser Dye L037 Common name Pyromethene567	
Class Laser Dye L027 Common name Coumarin519		Class Laser Dye L038 Common name Pyromethene580	
Class Laser Dye L028 Common name Coumarin521		Class Laser Dye L039 Common name Pyromethene597	
Class Laser Dye L029 Common name Coumarin523		Class Laser Dye L040 Common name Pyromethene597-8C9	
Class Laser Dye L030 Common name Coumarin525		Class Laser Dye L041 Common name Pyromethene605	












<p>Class Laser Dye L042 Common name Pyromethene650</p>		<p>Class Laser Dye L053 Common name Rhodamine590Perchlorate</p>	
<p>Class Laser Dye L043 Common name Fluorescein548</p>		<p>Class Laser Dye L054 Common name Rhodamine610Chloride</p>	
<p>Class Laser Dye L044 Common name DisodiumFluorescein</p>		<p>Class Laser Dye L055 Common name Rhodamine610Perchlorate</p>	
<p>Class Laser Dye L045 Common name Fluorol555</p>		<p>Class Laser Dye L056 Common name Rhodamine610Tetrafluoro borate</p>	
<p>Class Laser Dye L046 Common name Rhodamine3BPerchlorate</p>		<p>Class Laser Dye L057 Common name KitonRed620</p>	
<p>Class Laser Dye L047 Common name Rhodamine560Chloride</p>		<p>Class Laser Dye L058 Common name Rhodamine640Perchlorate</p>	
<p>Class Laser Dye L048 Common name Rhodamine560Perchlorate</p>		<p>Class Laser Dye L059 Common name Sulforhodamine640</p>	
<p>Class Laser Dye L049 Common name Rhodamine575</p>		<p>Class Laser Dye L060 Common name DODChloride</p>	 <p><i>In Aurora concise</i></p>
<p>Class Laser Dye L050 Common name Rhodamine19Perchlorate</p>		<p>Class Laser Dye L061 Common name LD688</p>	
<p>Class Laser Dye L051 Common name Rhodamine590Chloride</p>		<p>Class Laser Dye L062 Common name LDS750</p>	 <p><i>In Aurora concise</i></p>
<p>Class Laser Dye L052 Common name Rhodamine590Tetrafluoro borate</p>		<p>Class Laser Dye L063 Common name LDS759</p>	 <p><i>In Aurora concise</i></p>










Dyes in Aurora concise:,,,












<p>Class Laser Dye L064 Common name LDS821</p>  <p><i>In Aurora concise</i></p>	<p>Class Laser Dye L075 Common name Oxazine725Chloride</p> 
<p>Class Laser Dye L065 Common name LDS925</p> 	<p>Class Laser Dye L076 Common name Oxazine750Perchlorate</p> 
<p>Class Laser Dye L066 Common name LDS950</p> 	<p>Class Laser Dye L077 Common name HIDClodide</p> 
<p>Class Laser Dye L067 Common name LD800</p> 	<p>Class Laser Dye L078 Common name DOTClodide</p> 
<p>Class Laser Dye L068 Common name Phenoxazine660</p> 	<p>Class Laser Dye L079 Common name DOTPerchlorate</p> 
<p>Class Laser Dye L069 Common name CresylViolet670Perchlorate</p> 	<p>Class Laser Dye L080 Common name HITCPerchlorate</p> 
<p>Class Laser Dye L070 Common name NileBlue690Perchlorate</p> 	<p>Class Laser Dye L081 Common name HITClodide</p> 
<p>Class Laser Dye L071 Common name LD690Perchlorate</p> 	<p>Class Laser Dye L082 Common name HITCSpecial</p> 
<p>Class Laser Dye L072 Common name LD700Perchlorate</p> 	<p>Class Laser Dye L083 Common name DTTClodide</p> 
<p>Class Laser Dye L073 Common name Oxazine720Perchlorate</p> 	<p>Class Laser Dye L084 Common name IR-26</p> 
<p>Class Laser Dye L074 Common name Oxazine725Perchlorate</p> 	<p>Class Laser Dye L085 Common name IR-125</p> 












Class Laser Dye L097 Common name DQTCI	 <i>In Aurora concise</i>	Class Laser Dye L086 Common name IR-140	
Class Laser Dye L098 Common name DQTCI(*)	 <i>In Aurora concise</i>	Class Laser Dye L087 Common name IR-143	
Class Laser Dye L099 Common name DTCI	 <i>In Aurora concise</i>	Class Laser Dye L088 Common name IR-144	
Class Laser Dye L100 Common name Cryptocyanine	 <i>In Aurora concise</i>	Class Laser Dye L089 Common name DASPI	 <i>In Aurora concise</i>
Class Laser Dye L101 Common name DTDCI	 <i>In Aurora concise</i>	Class Laser Dye L090 Common name PICI	 <i>In Aurora concise</i>
Class Laser Dye L102 Common name NCI	 <i>In Aurora concise</i>	Class Laser Dye L091 Common name DMETCI	 <i>In Aurora concise</i>
Class Laser Dye L103 Common name HDITCP	 <i>In Aurora concise</i>	Class Laser Dye L092 Common name DOCI	 <i>In Aurora concise</i>
Class Laser Dye L104 Common name DQTrCI	 <i>In Aurora concise</i>	Class Laser Dye L093 Common name HICI	 <i>In Aurora concise</i>
Class Laser Dye L105 Common name Q-SwitchI(Kodak9740)	 <i>In Aurora concise</i>	Class Laser Dye L094 Common name DQOCI	 <i>In Aurora concise</i>
Class Laser Dye L106 Common name Q-Switch5	 <i>In Aurora concise</i>	Class Laser Dye L095 Common name Pinacyanol	 <i>In Aurora concise</i>
Class Laser Dye L107 Common name MalachiteGreen	 <i>In Aurora concise</i>	Class Laser Dye L096 Common name DDBCi	 <i>In Aurora concise</i>












Class Laser Dye L108 Common name BDN	
Class Max Weaver Set 1 Dye MWA001 Common name MWA001	
Class Max Weaver Set 1 Dye MWA002 Common name MWA002	
Class Max Weaver Set 1 Dye MWA003 Common name MWA003	
Class Max Weaver Set 1 Dye MWA004 Common name MWA004	
Class Max Weaver Set 1 Dye MWA005 Common name MWA005	
Class Max Weaver Set 1 Dye MWA006 Common name MWA006	
Class Max Weaver Set 1 Dye MWA007 Common name MWA007	
Class Max Weaver Set 1 Dye MWA008 Common name MWA008	
Class Max Weaver Set 1 Dye MWA010 Common name MWA010	
Class Max Weaver Set 1 Dye MWA011 Common name MWA011	























Class Max Weaver Set 1 Dye MWA013 Common name MWA013	
Class Max Weaver Set 1 Dye MWA014 Common name MWA014	
Class Max Weaver Set 1 Dye MWA015 Common name MWA015	
Class Max Weaver Set 1 Dye MWA016 Common name MWA016	
Class Max Weaver Set 1 Dye MWA017 Common name MWA017	
Class Max Weaver Set 1 Dye MWA018 Common name MWA018	
Class Max Weaver Set 1 Dye MWA020 Common name MWA020	
Class Max Weaver Set 1 Dye MWA021 Common name MWA021	
Class Max Weaver Set 2 Dye MWB001 Common name MWB001	
Class Max Weaver Set 2 Dye MWB002 Common name MWB002	
Class Max Weaver Set 2 Dye MWB003 Common name MWB003	












Class Max Weaver Set 2 Dye MWB004 Common name MWB004	
Class Max Weaver Set 2 Dye MWB005 Common name MWB005	
Class Max Weaver Set 2 Dye MWB006 Common name MWB006	
Class Max Weaver Set 2 Dye MWB007 Common name MWB007	
<i>Dye MWB008 not received/tested</i>	
Class Max Weaver Set 2 Dye MWB009 Common name MWB009	 <i>In Aurora concise</i>
<i>Dye MWB011 not received/tested</i>	
Class Max Weaver Set 2 Dye MWB011 Common name MWB011	
Class Max Weaver Set 2 Dye MWB012 Common name MWB012	
Class Max Weaver Set 2 Dye MWB013 Common name MWB013	
Class Max Weaver Set 2 Dye MWB014 Common name MWB014	












Class Max Weaver Set 2 Dye MWB015 Common name MWB015	
Class Max Weaver Set 2 Dye MWB016 Common name MWB016	
Class Max Weaver Set 2 Dye MWB017 Common name MWB017	
Class Max Weaver Set 2 Dye MWB019 Common name MWB019	
Class Max Weaver Set 3 Dye MWC001 Common name MWC001	
Class Max Weaver Set 3 Dye MWC002 Common name MWC002	
Class Max Weaver Set 3 Dye MWC003 Common name MWC003	
Class Max Weaver Set 3 Dye MWC004 Common name MWC004	
Class Max Weaver Set 3 Dye MWC005 Common name MWC005	
Class Max Weaver Set 3 Dye MWC006 Common name MWC006	
Class Max Weaver Set 3 Dye MWC007 Common name MWC007	 <i>In Aurora concise</i>












Class Max Weaver Set 3 Dye MWC008 Common name MWC008	 <i>In Aurora concise</i>
Class Max Weaver Set 3 Dye MWC009 Common name MWC009	
Class Max Weaver Set 3 Dye MWC010 Common name MWC010	 <i>In Aurora concise</i>
Class Max Weaver Set 3 Dye MWC011 Common name MWC011	
Class Max Weaver Set 3 Dye MWC012 Common name MWC012	
Class Max Weaver Set 3 Dye MWC013 Common name MWC013	
Class Max Weaver Set 3 Dye MWC014 Common name MWC014	
Class Max Weaver Set 3 Dye MWC015 Common name MWC015	
Class Max Weaver Set 3 Dye MWC016 Common name MWC016	
Class Max Weaver Set 3 Dye MWC017 Common name MWC017	
Class Max Weaver Set 3 Dye MWC018 Common name MWC018	












Class Max Weaver Set 3 Dye MWC019 Common name MWC019	
Class Max Weaver Set 3 Dye MWC020 Common name MWC020	
Class Max Weaver Set 3 Dye MWC021 Common name MWC021	 <i>In Aurora concise</i>
Class Max Weaver Set 3 Dye MWC022 Common name MWC022	
Class Max Weaver Set 3 Dye MWC023 Common name MWC023	
Class Max Weaver Set 3 Dye MWC024 Common name MWC024	
Class Max Weaver Set 3 Dye MWC025 Common name MWC025	
Class Max Weaver Set 3 Dye MWC026 Common name MWC026	
Class Max Weaver Set 3 Dye MWC027 Common name MWC027	
Class Max Weaver Set 3 Dye MWC028 Common name MWC028	
Class Max Weaver Set 3 Dye MWC029 Common name MWC029	




























Class Max Weaver Set 3 Dye MWC052 Common name MWC052	
Class Max Weaver Set 3 Dye MWC053 Common name MWC053	 <i>In Aurora concise</i>
Class Max Weaver Set 3 Dye MWC054 Common name MWC054	
Class Max Weaver Set 3 Dye MWC055 Common name MWC055	
Class Max Weaver Set 3 Dye MWC056 Common name MWC056	
Class Max Weaver Set 3 Dye MWC057 Common name MWC057	
Class Max Weaver Set 3 Dye MWC058 Common name MWC058	
Class Max Weaver Set 3 Dye MWC059 Common name MWC059	
Class Max Weaver Set 3 Dye MWC060 Common name MWC060	
Class Max Weaver Set 3 Dye MWC061 Common name MWC061	 <i>In Aurora concise</i>
Class Max Weaver Set 3 Dye MWC062 Common name MWC062	

Class Max Weaver Set 4 Dye MWD01 Common name MWD01	
Class Max Weaver Set 4 Dye MWD02 Common name MWD02	
Class Max Weaver Set 4 Dye MWD03 Common name MWD03	
Class Max Weaver Set 4 Dye MWD04 Common name MWD04	
Class Max Weaver Set 4 Dye MWD05 Common name MWD05	
Class Max Weaver Set 4 Dye MWD06 Common name MWD06	
Class Max Weaver Set 4 Dye MWD07 Common name MWD07	
Class Max Weaver Set 4 Dye MWD08 Common name MWD08	
Class Max Weaver Set 4 Dye MWD09 Common name MWD09	
Class Max Weaver Set 4 Dye MWD10 Common name MWD10	
Class Max Weaver Set 4 Dye MWD11 Common name MWD11	

Class Max Weaver Set 4 Dye MWD12 Common name MWD12	
Class Max Weaver Set 4 Dye MWD13 Common name MWD13	
Class Max Weaver Set 4 Dye MWD14 Common name MWD14	
Class Max Weaver Set 4 Dye MWD15 Common name MWD15	
Class Max Weaver Set 4 Dye MWD16 Common name MWD16	
Class Max Weaver Set 4 Dye MWD17 Common name MWD17	
Class Max Weaver Set 5 Dye MWE01 Common name MWE01	
Class Max Weaver Set 5 Dye MWE02 Common name MWE02	
Class Max Weaver Set 5 Dye MWE03 Common name MWE03	
Class Max Weaver Set 5 Dye MWE04 Common name MWE04	
Class Max Weaver Set 5 Dye MWE05 Common name MWE05	

Class Max Weaver Set 5 Dye MWE06 Common name MWE06	
Class Max Weaver Set 5 Dye MWE07 Common name MWE07	
Class Max Weaver Set 5 Dye MWE08 Common name MWE08	
Class Max Weaver Set 5 Dye MWE09 Common name MWE09	
Class Max Weaver Set 5 Dye MWE10 Common name MWE10	
Class Max Weaver Set 5 Dye MWE11 Common name MWE11	
Class Max Weaver Set 5 Dye MWE12 Common name MWE12	
Class Max Weaver Set 5 Dye MWE13 Common name MWE13	
Class Max Weaver Set 5 Dye MWE14 Common name MWE14	
Class Max Weaver Set 6 Dye MWF01 Common name MWF01	
Class Max Weaver Set 6 Dye MWF02 Common name MWF02	

<p>Class Max Weaver Set 6 Dye MWF03 Common name MWF03</p>		<p>Class Commercial textile, order 1 Dye T009 Common name Direct Red 224</p>	
<p>Class Max Weaver Set 6 Dye MWF04 Common name MWF04</p>		<p>Class Commercial textile, order 1 Dye T010 Common name Direct Red 23</p>	
<p>Class Max Weaver Set 6 Dye MWF05 Common name MWF05</p>		<p>Class Commercial textile, order 1 Dye T011 Common name Direct Red 81</p>	
<p>Class Max Weaver Set 6 Dye MWF06 Common name MWF06</p>		<p>Class Commercial textile, order 1 Dye T014 Common name Direct Violet 51</p>	
<p>Class Max Weaver Set 6 Dye MWF07 Common name MWF07</p>		<p>Class Commercial textile, order 1 Dye T016 Common name Direct Blue 199</p>	
<p>Class Max Weaver Set 6 Dye MWF08 Common name MWF08</p>	 <i>In Aurora concise</i>	<p>Class Commercial textile, order 1 Dye T027 Common name Direct Black 80</p>	
<p>Class SYPRO Dye SYPRO Common name SYPRO Orange</p>	 <i>In Aurora concise</i>	<p>Class SYPROlogs Dye TW383 Common name SYPROlog1</p>	
<p>Class Commercial textile, order 1 Dye T001 Common name Direct Brown 210</p>	 <i>In Aurora concise</i>	<p>Class SYPROlogs Dye TW384* Common name SYPROlog2</p> <p><i>*synthesized in-house (not commercially available)</i></p>	 <i>In Aurora concise</i>
<p>Class Commercial textile, order 1 Dye T002 Common name Direct Yellow 96</p>		<p>Class SYPROlogs Dye TW385* Common name SYPROlog3</p> <p><i>*synthesized in-house (not commercially available)</i></p>	
<p>Class Commercial textile, order 1 Dye T004 Common name Quinoline Yellow WS</p>	 <i>In Aurora concise</i>	<p>Class SYPROlogs Dye TW408* Common name SYPROlog4</p> <p><i>*synthesized in-house (not commercially available)</i></p>	 <i>In Aurora concise</i>
<p>Class Commercial textile, order 1 Dye T008 Common name Direct Red 227</p>		<p>Class SYPROlogs Dye TW418* Common name SYPROlog5</p> <p><i>*synthesized in-house (not commercially available)</i></p>	 <i>In Aurora concise</i>

Class SYPROlogs Dye TW419* Common name SYPROlog6 <i>*synthesized in-house (not commercially available)</i>		<i>In Aurora concise</i>
Class SYPROlogs Dye TW495* Common name SYPROlog7 <i>*synthesized in-house (not commercially available)</i>		<i>In Aurora concise</i>
Class SYPROlogs Dye AW131* Common name SYPROlog8 <i>*synthesized in-house (not commercially available)</i>		<i>In Aurora concise</i>



Class SYPROlogs Dye AW134* Common name SYPROlog9 <i>*synthesized in-house (not commercially available)</i>		<i>In Aurora concise</i>
Class SYPROlogs Dye AW132* Common name SYPROlog10 <i>*synthesized in-house (not commercially available)</i>		

Table S1. Chemical structures and common names for all Aurora library dyes. Chemical structures, Aurora library name, and common name (where applicable) for all 312 tested dyes in the Aurora library. Unless otherwise noted, dyes are available either commercially, or through the North Carolina State Max Weaver Dye Library. Detailed information on dye sources is available in Data S1.

Complete list of dyes in Aurora concise: TW408, MWF08, TW495, MWC061, L063, L096, A022, A016, T004, L091, L097, L093, A014, L095, L098, L094, SYPRO Orange, A003, L062, A024, MWC007, A002, A020, L060, A001, L089, L099, TW384, A013, MWC053, TW419, A009, A015, C003, L064, A021, MWC034, AW134, MWC045, C010, L090, A023, MWC021, MWC008, A010, AW131, MWB009, C001, TW418

Protein				Dye screen conditions		Biochemical / structural properties							Source / references			
Abbreviated protein name	Full protein name	Construct details	Annotation	Screening conc. (µM)	Screening buffer	Sequence Length	Molecular Weight	GRAVY hydrophobicity	Percent disordered	Fraction Helix	Fraction Turn	Fraction Sheet	Protein source	Uniprot ID	PDB Code	Ref.
ACD_domain_protein_6	Alpha-crystallin B chain	Alpha crystallin domain of HSPB5	Eye lens refractivity, chaperone-like properties	10	25 mM MOPS, 150 mM NaCl, pH 7.50	89	10180.3614	-0.652809	10	0.30337079	0.2247191	0.21348315	Maria Janowska, Kleivit Lab, UW	P02511		Rajagopal, Ponn, et al. 2015
B5_pre_inc_phospho_protein_4	Alpha-crystallin B chain	phosphomimetic construct, incubated at 37C to activate	Oligomeric/conformational modulation of HSPB5	10	25 mM MOPS, 150 mM NaCl, pH 7.50	175	20242.6725	-0.5777143	21	0.30285714	0.23428571	0.21142857	Maria Janowska, Kleivit Lab, UW	P02511		Rajagopal, Ponn, et al. 2015
B5_with_zinc_protein_3	Alpha-crystallin B chain	HSPB5 in the presence of Zinc	Oligomeric/conformational modulation of HSPB5	10	25 mM MOPS, 150 mM NaCl, pH 7.50	175	20158.6422	-0.5314286	21	0.30285714	0.25142857	0.21142857	Maria Janowska, Kleivit Lab, UW	P02511		Rajagopal, Ponn, et al. 2015
Bag2	BAG family molecular chaperone regulator 2		Co-chaperone, nucleotide-exchange factor	5	10 mM HEPES, pH 7.40, 200 mM NaCl	211	23771.6663	-0.5919431	22	0.2464455	0.20853081	0.32701422	Rebecca Freilich, Gestwicki Lab, UCSF	O95816		
BSA	Bovine Serum Albumin		Blood osmotic pressure regulation; ion and small molecule binding	10	20 mM HEPES, pH 7.5, 150 mM KCl, 20 mM NaCl, 5 mM MgCl2	607	69292.6835	-0.4289951	0	0.28336079	0.14991763	0.2915902	Commercial - Sigma Aldrich Cat A2153-100G >96%	P02769	3V03	
caspase1_zymogen	Caspase-1	Zymogen form	Protease, inflammatory processes	1	50 mM HEPES, pH 7.40, 250 mM NaCl	404	45158.155	-0.3334158	11	0.26980198	0.22524752	0.25247525	Andy Ambrose, Arkin Lab, UCSF	P29466	18MQ	
caspase2_active	Caspase-2	Active form	Protease, apoptotic cascade	1	50 mM HEPES, pH 7.40, 250 mM NaCl	452	50684.7558	-0.3011062	6	0.28318584	0.20132743	0.29867257	Andy Ambrose, Arkin Lab, UCSF	P42575	1PYO	
caspase3_active	Caspase-3	Active form	Protease, apoptotic cascade	1	50 mM HEPES, pH 7.40, 250 mM NaCl	277	31607.5495	-0.4570397	12	0.28519856	0.23104693	0.22382671	Andy Ambrose, Arkin Lab, UCSF	P42574	1CP3	
caspase6_active	Caspase-6	Active form	Protease, programmed cell death, development, innate immunity, inflammation	1	50 mM HEPES, pH 7.40, 250 mM NaCl	293	33309.597	-0.3552901	6	0.29692833	0.19795222	0.24573379	Andy Ambrose, Arkin Lab, UCSF	P55212	2WDP	
caspase6_zymogen	Caspase-6	Zymogen form	Protease, programmed cell death, development, innate immunity, inflammation	1	50 mM HEPES, pH 7.40, 250 mM NaCl	293	33309.597	-0.3552901	6	0.29692833	0.19795222	0.24573379	Andy Ambrose, Arkin Lab, UCSF	P55212	2WDP	
caspase7_active	Caspase-7	Active form	Protease, apoptotic cascade	1	50 mM HEPES, pH 7.40, 250 mM NaCl	303	34276.4229	-0.5089109	9	0.27392739	0.21452145	0.21122112	Andy Ambrose, Arkin Lab, UCSF	P55210	1F1J	
caspase9_active	Caspase-9	Active form	Protease, apoptotic cascade	1	50 mM HEPES, pH 7.40, 250 mM NaCl	416	46280.3047	-0.2370192	10	0.29807692	0.24519231	0.25721154	Andy Ambrose, Arkin Lab, UCSF	P55212	1JXQ	
CBP_KIX	Histone lysine acetyltransferase CREBBP		Histone acetylation, transcriptional activation	10	10 mM HEPES, pH 7.20, 200 mM NaCl	81	9534.7203	-0.8037037	0	0.2962963	0.13580247	0.32098765	Partch Lab, UCSC	P45481	1KDX_1	Gustafson, Chelsea L, et al. 2017
CHIP	E3 ubiquitin-protein ligase CHIP		E3 ubiquitin ligase	1	10 mM HEPES, pH 7.20, 200 mM NaCl, 1 mM TCEP	303	34855.8215	-0.8108911	20	0.24422442	0.18151815	0.31353135	Matt Ravalin, Gestwicki Lab, UCSF	Q9UNE7	4KBQ	Ravalin M, et al. 2019
CLOCK_PASB_W362A	Circadian locomotor output cycles protein kaput	point mutant reduces formation of soluble aggregates at room temp	Transcriptional activator, circadian rhythm	10	10 mM HEPES, pH 7.20, 200 mM NaCl	361	42025.9864	-0.5260388	8	0.31301939	0.19944598	0.23822715	Partch Lab, UCSC	O08785		Michael, Alicia K, et al. 2017
CoV2Orf9bMark2	SARS-CoV2 complex Orf9b-Mark2	screened as pre-assembled complex	SARS-CoV2 host-virus protein complex	10	50 mM Tris HCL, pH 7.5, 150 mM NaCl, 0.1 mM TCEP								Amy Diallo, Verba Lab, UCSF	POD1D1		
OP_BMAL1_PASB_W427A	Aryl hydrocarbon receptor nuclear translocator-like protein 1	point mutant reduces formation of soluble aggregates at room temp	Transcriptional activator, circadian rhythm	10	10 mM HEPES, pH 7.20, 200 mM NaCl	387	44123.1476	-0.4434109	7	0.28423773	0.20671835	0.23514212	Partch Lab, UCSC	O00327	4F3I	Michael, Alicia K, et al. 2017
CP581_CBPKIX	CREB-binding protein		Histone acetylation, transcriptional activation	10	10 mM HEPES, pH 7.20, 200 mM NaCl	80	9405.6063	-0.77	0	0.3	0.1375	0.3125	Partch Lab, UCSC	Q92793	1S80	Gustafson, Chelsea L, et al. 2017
diUb	K48 linked di-Ubiquitin		Post-translational protein tag, protein regulation and degradation	5	50 mM Tris HCL, 150 mM NaCl, 1% DMSO, 0.01% Triton X-100, 1 mM TCEP	76	8548.6933	-0.5605263	26	0.28947368	0.19736842	0.22368421	Commercial - R&D Systems UCN-200	POCG48	1UBQ	
eIF2B_with_ISRI_B_all	Eukaryotic translation initiation factor 2 subunit 2	Screened as pre-assembled complex	Initiation factor in eukaryotic protein synthesis	0.25	50 mM TRIS, pH 7.4, 150 mM NaCl, 2 mM MgCl2, 1 mM DTT	2369	262753.142	-0.2064584		0.30434783	0.22119038	0.28366399	Michael Schoof, Walter Lab, UCSF	6CAJ		Tsai JC, et al. 2018
eIF2B_with_ISRI_B_alpha	Eukaryotic translation initiation factor 2 subunit 2, alpha subunit	Screened as pre-assembled complex	Initiation factor in eukaryotic protein synthesis	0.25	50 mM TRIS, pH 7.4, 150 mM NaCl, 2 mM MgCl2, 1 mM DTT	305	33711.8019	0.07409836	0	0.3442623	0.17704918	0.29508197	Michael Schoof, Walter Lab, UCSF	6CAJ_3		Tsai JC, et al. 2018
eIF2B_with_ISRI_B_beta	Eukaryotic translation initiation factor 2 subunit 2, beta subunit	Screened as pre-assembled complex	Initiation factor in eukaryotic protein synthesis	0.25	50 mM TRIS, pH 7.4, 150 mM NaCl, 2 mM MgCl2, 1 mM DTT	368	40952.1788	-0.2054348	22	0.28804348	0.22282609	0.30706522	Michael Schoof, Walter Lab, UCSF	6CAJ_4		Tsai JC, et al. 2018
eIF2B_with_ISRI_B_delta	Eukaryotic translation initiation factor 2 subunit 2, delta subunit	Screened as pre-assembled complex	Initiation factor in eukaryotic protein synthesis	0.25	50 mM TRIS, pH 7.4, 150 mM NaCl, 2 mM MgCl2, 1 mM DTT	523	57556.561	-0.3544933	32	0.27724665	0.23709369	0.29063098	Michael Schoof, Walter Lab, UCSF	6CAJ_5		Tsai JC, et al. 2018
eIF2B_with_ISRI_B_epsilon	Eukaryotic translation initiation factor 2 subunit 2, epsilon subunit	Screened as pre-assembled complex	Initiation factor in eukaryotic protein synthesis	0.25	50 mM TRIS, pH 7.4, 150 mM NaCl, 2 mM MgCl2, 1 mM DTT	721	80364.9077	-0.2768377	17	0.30651872	0.23162275	0.27461859	Michael Schoof, Walter Lab, UCSF	6CAJ_1		Tsai JC, et al. 2018
eIF2B_with_ISRI_B_gamma	Eukaryotic translation initiation factor 2 subunit 2, gamma subunit	Screened as pre-assembled complex	Initiation factor in eukaryotic protein synthesis	0.25	50 mM TRIS, pH 7.4, 150 mM NaCl, 2 mM MgCl2, 1 mM DTT	452	50239.7541	-0.1130531	4	0.31858407	0.21460177	0.26327434	Michael Schoof, Walter Lab, UCSF	6CAJ_2		Tsai JC, et al. 2018
FKBP12	Peptidyl-prolyl cis-trans isomerase FKBP1A		Kinase receptor regulation	10	20 mM Tris, pH 8.0	108	11950.569	-0.4111111	15	0.27777778	0.24074074	0.21296296	Douglas Wassarman, Shokat Lab, UCSF	P62942		Wassarman, Douglas R., et al. 2022
FL_B5_no_inc_protein_1	Alpha-crystallin B chain	Phosphomimetic mutant, no pre-incubation	Oligomeric/conformational modulation of HSPB5	10	25 mM MOPS, 150 mM NaCl, pH 7.50	175	20158.6422	-0.5314286	21	0.30285714	0.25142857	0.21142857	Maria Janowska, Kleivit Lab, UW	P02511		Rajagopal, Ponn, et al. 2015
FL_B5_pre_inc_protein_2	Alpha-crystallin B chain	pre-incubated at 37C to activate	Oligomeric/conformational modulation of HSPB5	10	25 mM MOPS, 150 mM NaCl, pH 7.50	175	20158.6422	-0.5314286	21	0.30285714	0.25142857	0.21142857	Maria Janowska, Kleivit Lab, UW	P02511		Rajagopal, Ponn, et al. 2015
GB1	Immunoglobulin G-binding protein G		Binds to immunoglobulin Fc region	10	10 mM HEPES, pH 7.20, 200 mM NaCl	56	6151.7266	-0.5392857	4	0.25	0.125	0.28571429	Partch Lab, UCSC	P06654	2J52	
GS_R324C	Glutamine synthetase	Pathogenic mutant of glutamine synthetase	Implicated in Glutamine Deficiency, a rare genetic disease	10	20 mM HEPES, pH 7.5, 150 mM KCl, 20 mM NaCl, 5 mM MgCl2	384	43252.2819	-0.5648437	5	0.25520833	0.26302083	0.22916667	Erin Thompson, Fraser Lab, UCSF	P15104		Thompson, Erin et al. 2020
GS_R324S	Glutamine synthetase	Pathogenic mutant of glutamine synthetase	Implicated in Glutamine Deficiency, a rare genetic disease	10	20 mM HEPES, pH 7.5, 150 mM KCl, 20 mM NaCl, 5 mM MgCl2	384	43236.2163	-0.5734375	7	0.25520833	0.265625	0.22916667	Erin Thompson, Fraser Lab, UCSF	P15104		Thompson, Erin et al. 2020
GS_R341C	Glutamine synthetase	Pathogenic mutant of glutamine synthetase	Implicated in Glutamine Deficiency, a rare genetic disease	10	20 mM HEPES, pH 7.5, 150 mM KCl, 20 mM NaCl, 5 mM MgCl2	384	43252.2819	-0.5648437	4	0.25520833	0.26302083	0.22916667	Erin Thompson, Fraser Lab, UCSF	P15104		Thompson, Erin et al. 2020
GS_WT	Glutamine synthetase	Glutamine synthetase	Metabolic enzyme, ATP-dependent conversion of glutamate to glutamine	10	20 mM HEPES, pH 7.5, 150 mM KCl, 20 mM NaCl, 5 mM MgCl2	384	43305.3247	-0.5830729	7	0.25520833	0.26302083	0.22916667	Erin Thompson, Fraser Lab, UCSF	P15104	2OJW	Thompson, Erin et al. 2020
GST	Glutathione S-transferase A3		Enzymatic conjugation of reduced glutathione to various electrophiles	1	10 mM HEPES, pH 7.20, 200 mM NaCl	222	25740.7871	-0.3189189	3	0.32882883	0.16216216	0.32882883	Partch Lab, UCSC	Q16772	5LCZ	
H104_FL_B5_pre_inc_protein_5	Alpha-crystallin B chain	Phosphomimetic mutant of HSPB5, incubated at 37C to activate	Oligomeric/conformational modulation of HSPB5	10	25 mM MOPS, 150 mM NaCl, pH 7.50	175	20149.6752	-0.5354286	20	0.30285714	0.25142857	0.21142857	Maria Janowska, Kleivit Lab, UW	P02511		Rajagopal, Ponn, et al. 2020

hACE2	Angiotensin-converting enzyme 2	Hormone regulation in cardiovascular homeostasis	Carboxypeptidase essential for cardiovascular homeostasis	0.5	50 mM Tris HCl, 150 mM NaCl, 1% DMSO, 0.01% Triton X-100, 1 mM TCEP	723	83594.7348	-0.4614108	2	0.31673582	0.22821577	0.26970954	Commercial - R&D Systems R&D 933-ZN-101	Q9BYF1	1R42	
HIP	Hsc70-interacting protein		State-dependent protein-interacting partner for nucleotide domain of Hsp70	1	10 mM HEPES, pH 7.20, 200 mM NaCl	369	41331.2623	-0.8574526	69	0.19783198	0.21680217	0.34146341	Rebecca Freilich, Gestwicki Lab, UCSF	P50502		
His_NusA	Transcription termination/antitermination protein NusA		rRNA transcription and antitermination complex component	10	10 mM HEPES, pH 7.20, 200 mM NaCl	329	36252.7006	-0.2179331	2	0.2887538	0.16717325	0.32826748	Partch Lab, UCSC	P0AFF6	5LM9	
HisPER2_AB_I324 N	Period circadian protein homolog 2	324N "Edo" mutant in interdomain linker	Transcriptional repressor complex component in circadian rhythm	10	10 mM HEPES, pH 7.20, 200 mM NaCl	304	34638.3344	-0.2555921	9	0.31907895	0.23355263	0.22697368	Partch Lab, UCSC	O54943	4F3L	Militti, Stefania, et al. 2016
HisSUMO	Small ubiquitin-related modifier 4		Post-translational protein tag	1	10 mM HEPES, pH 7.20, 200 mM NaCl	80	9237.4484	-0.945	30	0.2375	0.2	0.25	Partch Lab, UCSC	Q6EEV6	2ASQ	
Hsc70	Heat shock cognate 71 kDa protein		Molecular chaperone	10	10 mM HEPES, pH 7.20, 200 mM NaCl	646	70897.2343	-0.4560372	20	0.27089783	0.22755418	0.24148607	Rebecca Freilich, Gestwicki Lab, UCSF	P11142	3AGY	
Hsp10	10 kDa heat shock protein, mitochondrial		Mitochondrial molecular chaperonin	10	10 mM HEPES, pH 7.20, 200 mM NaCl	80	8343.6766	-0.02	4	0.3	0.2625	0.25	Dr. Hao Shao, Gestwicki Lab, UCSF	P61604	4JP1	
Hsp60	60 kDa heat shock protein, mitochondrial		Mitochondrial molecular chaperonin	10	DPBS	558	59777.8702	-0.168638	10	0.28853047	0.20430108	0.27598566	Dr. Hao Shao, Gestwicki Lab, UCSF	P10809	4JP1	
Hsp72_SBD	72kDa Heat shock protein, substrate binding domain	Substrate binding domain only	Molecular chaperone	10	25 mM Tween, 100 mM NaCl, 1mM TCEP, 1 mM EDTA	235	26241.3454	-0.5902128	23	0.26382979	0.20425532	0.28510638	Rebecca Freilich, Gestwicki Lab, UCSF	P54652		
ISG15	Ubiquitin-like protein ISGylation 15		Innate immune signalling	5	50 mM Tris HCl, 150 mM NaCl, 1% DMSO, 0.01% Triton X-100, 1 mM TCEP	165	17887.294	-0.1387879	5	0.31515152	0.27272727	0.26666667	Commercial - R&D Systems UL-601	P05161	122M	
KaiC_AE	Circadian clock protein kinase KaiC	Phosphomimetic of the S431/pT432 daytime state of KaiC	Cyanobacterial circadian rhythm oscillator, protein-protein interaction hub	10	25 mM Tris, pH 7.40, 130 mM NaCl, 1 mM MgCl2, 1 mM TCEP, 1 mM ATP	519	58014.2318	-0.2231214	6	0.30828516	0.22928709	0.24662813	Jeffery Swan, Partch Lab, UCSC	Q79V60	2GBL	
KaiC_EA	Circadian clock protein kinase KaiC	Phosphomimetic of the pS431/T432 nighttime state of KaiC	Cyanobacterial circadian rhythm oscillator, protein-protein interaction hub	10	25 mM Tris, pH 7.40, 130 mM NaCl, 1 mM MgCl2, 1 mM TCEP, 1 mM ATP	519	58014.2318	-0.2231214	6	0.30828516	0.22928709	0.24662813	Jeffery Swan, Partch Lab, UCSC	Q79V60	2GBL	Tseng, Roger, et al. 2017
klv_light_chain	Kappa-1 variable domain of immunoglobulin light chain		Aggregates in genetic disease Light Chain Amyloidosis	17.5	50 mM TRIS, pH 7.4, 150 mM NaCl, 2 mM MgCl2	108	11932.0587	-0.5222222	5	0.27777778	0.2962963	0.15740741	Marina Ramirez-Alvarado			
lysozyme	Hen egg white lysozyme		Egg white innate immunity, model folding protein	5	10 mM HEPES, pH 7.20, 200 mM NaCl	129	14313.004	-0.472093	0	0.24806202	0.29457364	0.18604651	Commercial - Sigma Adrich 10837059001	P00698	1DPX	
mortalin	Stress-70 protein, mitochondrial "Mortalin"		Mitochondrial chaperone	1	20 mM HEPES, pH 7.50, 100 mM NaCl	679	73679.6436	-0.4002946	24	0.25184094	0.21502209	0.27982327	Dr. Hao Shao, Gestwicki Lab, UCSF	P38646	3N8E	
MPro	SARS-CoV2 3C-like proteinase "Main protease"		Viral protease	1	50 mM Tris HCl, 150 mM NaCl, 1% DMSO, 0.01% Triton X-100, 1 mM TCEP	306	33796.295	-0.0189542	0	0.32026144	0.24836601	0.2124183	Commercial - R&D Systems E-720-050	POD0T1	622E	
nrh23_LBD	Nuclear hormone receptor 23	ligand binding domain alone	Hormone receptor involved in skin regeneration and molting of C. elegans	1	10 mM HEPES, pH 7.20, 200 mM NaCl, 1 mM TCEP	242	26998.3636	-0.0681818	2	0.30991736	0.18181818	0.31404959	Partch Lab, UCSC	P41828		
NPAS2_PAS_A	Neuronal PAS domain-containing protein 2	Pas A domain		10	10 mM HEPES, pH 7.20, 200 mM NaCl	163	18383.6534	-0.0650307	0	0.35582822	0.27607362	0.2392638	Partch Lab, UCSC	P97460		
nsp10	SARS-CoV2 nonstructural protein 10	screened in isolation (not in presence of nsp16 or 14)	Viral immune evasion complex component	1	50 mM HEPES, pH 8, 500 mM NaCl, 2 mM MgCl2, 0.5 mM TCEP	142	15062.0745	-0.0838028	4	0.22535211	0.25352113	0.21126761	Victor Long, Gross Lab, UCSF	POD0T1	6W4H	
nsp14	SARS-CoV2 nonstructural protein 14	screened in isolation (not in presence of nsp10)	Viral methyltransferase, immune evasion	1	50 mM HEPES, pH 8, 500 mM NaCl, 2 mM MgCl2, 0.5 mM TCEP	527	59815.0269	-0.1335863	1	0.33017078	0.21252372	0.19734345	Victor Long, Gross Lab, UCSF	POD0T1	7N0B_2	
nsp14_nsp10	Heterodimeric complex of SARS-CoV2 nonstructural protein 14 and 10	screened as a pre-assembled complex	Viral methyltransferase, immune evasion	1	50 mM HEPES, pH 8, 500 mM NaCl, 2 mM MgCl2, 0.5 mM TCEP	669	74859.0861	-0.1230194		0.30792227	0.22122571	0.20029895	Victor Long, Gross Lab, UCSF	POD0T1	7N0B	
nsp16	SARS-CoV2 nonstructural protein 16	screened in isolation (not in presence of nsp10)	Viral RNA capping, immune evasion	1	50 mM HEPES, pH 8, 500 mM NaCl, 2 mM MgCl2, 0.5 mM TCEP	301	33595.2116	-0.0933555	1	0.32890365	0.25249169	0.24252492	Victor Long, Gross Lab, UCSF	POD0T1	6W4H	
nsp16_nsp10	Heterodimeric complex of SARS-CoV2 nonstructural protein 16 and 10	screened as a pre-assembled complex	Viral RNA capping, immune evasion	1	50 mM HEPES, pH 8, 500 mM NaCl, 2 mM MgCl2, 0.5 mM TCEP	443	48639.2708	-0.0902935		0.29571106	0.25282167	0.23250564	Victor Long, Gross Lab, UCSF	POD0T1	6W4H	
nsp2	SARS-CoV2 nonstructural protein 2		Viral protein with putative zinc-finger	1	30 mM HEPES, pH 7.50, 60 mM NaCl, 1 mM TCEP, 10 uM ZnCl2	638	70510.6111	-0.062069	0	0.31661442	0.21630094	0.26332288	Meghna Gupta and Klim Verba, Stroud and Verba Lab, UCSF	POD0T1	7MSX_1	
nsp3_mac1	SARS-CoV2 nonstructural protein 3, macrodomain 1	macrodomain 1 of SARS-CoV2 nonstructural protein 3	Viral modulation/evasion of host innate immunity	2	50 mM Tris HCl, 150 mM NaCl, 1% DMSO, 0.01% Triton X-100, 1 mM TCEP	176	18971.4669	-0.0198864	0	0.34090909	0.26136364	0.25568182	Galen Correy, Fraser Lab, UCSF	POD0T1	6W0J	
nucleocapsid	SARS-CoV2 Nucleocapsid protein		Viral genome packing	1	50 mM Tris HCl, 150 mM NaCl, 1% DMSO, 0.01% Triton X-100, 1 mM TCEP	419	45625.1383	-0.9713604	68	0.18615752	0.31026253	0.19809069	Commercial - R&D Systems 10474-CV-050	POD0T1		
Nucleosome_all	Nucleosome	Screened as a fully assembled nucleosome with Widom 601 DNA	Key component of chromatin organization and regulation	1	TCS	495	54742.3457	-0.5757576		0.24848485	0.20202020	0.26464646	Hayden Saunders, Narlikar Lab, UCSF		2CV5_2	Gamarra, Nathan et al. 2021
Nucleosome_H2A	Histone protein 2A	Screen as a pre-assembled histone mixture	Key component of chromatin organization and regulation	1	TCS	130	14135.3121	-0.4776923	28	0.26153846	0.23076923	0.30769231	Hayden Saunders, Narlikar Lab, UCSF	Q6F113	2CV5_4	Gamarra, Nathan et al. 2021
Nucleosome_H2B	Histone protein 2B	Screen as a pre-assembled histone mixture	Key component of chromatin organization and regulation	1	TCS	126	13889.9692	-0.7174603	33	0.22222222	0.23015873	0.23809524	Hayden Saunders, Narlikar Lab, UCSF	Q93079	2CV5_5	Gamarra, Nathan et al. 2021
Nucleosome_H3	Histone protein 3	Screen as a pre-assembled histone mixture	Key component of chromatin organization and regulation	1	TCS	136	15403.9072	-0.5794118	34	0.23529412	0.19970588	0.29411765	Hayden Saunders, Narlikar Lab, UCSF	P68431	2CV5_2	Gamarra, Nathan et al. 2021
Nucleosome_H4	Histone protein 4	Screen as a pre-assembled histone mixture	Key component of chromatin organization and regulation	1	TCS	103	11367.2031	-0.5213592	20	0.2815534	0.21359223	0.2038835	Hayden Saunders, Narlikar Lab, UCSF	P62805	2CV5_3	Gamarra, Nathan et al. 2021
Nurr1_LBD	Nuclear receptor subfamily 4 group A member 2	Ligand binding domain only	Transcriptional regulator, dopaminergic neuron development	1		272	30627.0299	-0.1808824	12	0.33088235	0.25367647	0.26102941	Ward and Partch Lab, UCSC	P43354		
OGTase	UDP-N-acetylglucosamine-peptide N-acetylglucosaminyltransferase		O-linked GlcNAcylation enzyme, cellular process regulation	1	50 mM HEPES, 300 mM NaCl, pH 7.40, 1 mM DTT	1046	116923.285	-0.2424474	2	0.29923518	0.21797323	0.24945507	Matt Alteen, Volcadio Lab, U Alberta	O15294		Alteen, Matthew G., et al. 2020

p23	Prostaglandin E synthase 3		Hormone-dependent transcriptional activation	5	50 mM Tris, pH 9.0, 150 mM NaCl	160	18697.1939	-1.049375	26	0.2375	0.2375	0.2125	Andrew Ambrose, Arkin Lab, UCSF	Q15185	1E1F	
p300_KIX	Histone acetyltransferase p300	KIX domain	Histone acetylation, chromatin remodeling	1	10 mM HEPES, pH 7.20, 200 mM NaCl	80	9405.6063	-0.77	0	0.3	0.1375	0.3125	Partch Lab, UCSC	Q09472	2KWF	
PER2_PAS_ASB_wt	Period circadian protein homolog 2	wild-type tandem PAS A and B domains	Transcriptional repressor complex component in circadian rhythm	10	10 mM HEPES, pH 7.20, 200 mM NaCl	304	34637.3894	-0.2292763	7	0.32236842	0.23026316	0.22697368	Partch Lab, UCSC	O54943		Militi, Stefania, et al. 2016
PLPro	SARS-CoV2 Papain-like protease		Viral polyprotein cleavage, host immune evasion	1	50 mM Tris HCl, 150 mM NaCl, 1% DMSO, 0.01% Triton X-100, 1 mM TCEP	315	35631.1902	-0.3749206	0	0.2984127	0.19365079	0.22539683	Commercial - R&D Systems E-611	POD01	7CMD	
pp5	Serine/threonine-protein phosphatase 5		Modulation of diverse signalling pathways via dephosphorylation	1	20 mM HEPES, pH 7.50, 50 mM KCl, 1 mM TCEP	556	61859.1062	-0.1893885	5	0.30935252	0.22302158	0.28057554	Aye Thwin, Gestwicki and Southworth Lab, UCSF	P53041	1WAO	
PPIE	Peptidyl-prolyl cis-trans isomerase E		pre-mRNA splicing, spliceosome component	10	10 mM HEPES, pH 7.20, 200 mM NaCl	301	33430.4114	-0.5212625	14	0.26245847	0.23920266	0.23920266	Partch Lab, UCSC	Q9UNP9	3UCH	
PPIEdeltaN	Peptidyl-prolyl cis-trans isomerase E	isomerase domain only	pre-mRNA splicing, spliceosome component	10	10 mM HEPES, pH 7.20, 200 mM NaCl	170	18717.1834	-0.4417647	5	0.26470588	0.27058824	0.17058824	Partch Lab, UCSC	Q9UNP9	2R99	
PPIG	Peptidyl-prolyl cis-trans isomerase G		Proline isomerization, protein folding	10	10 mM HEPES, pH 7.20, 200 mM NaCl	754	88616.1421	-1.8831565	81	0.12732095	0.28647215	0.16180371	Partch Lab, UCSC	Q13427	2GW2	
PPIL1	Peptidyl-prolyl cis-trans isomerase-like 1		pre-mRNA splicing, spliceosome component	10	10 mM HEPES, pH 7.20, 200 mM NaCl	166	18236.6204	-0.3042169	5	0.28313253	0.26506024	0.19879518	Partch Lab, UCSC	Q9Y3C6	2X7K	
PPWD1	Peptidylprolyl isomerase domain and WD repeat-containing protein 1		Proline isomerization, protein folding	10	10 mM HEPES, pH 7.20, 200 mM NaCl	646	73574.0335	-0.4479876	16	0.28947368	0.21362229	0.21671827	Partch Lab, UCSC	Q96BP3	2A2N	
proteasome_alpha	Archaeal proteasome	screened as a single, pre-assembled complex with beta subunit	Non-specific protein degradation	1	50 mM Tris pH 7.5, 10 mM MgCl2, 1 mM DTT	233	25799.1671	-0.211588	0	0.32618026	0.19313305	0.27467811	Sarah Williams, Gestwicki Lab, UCSF	P25156	3IPM_1	Opoku-Nsiah, Kwadwo A., et al. 2021
proteasome_beta	Archaeal proteasome	screened as a single, pre-assembled complex with alpha subunit	Non-specific protein degradation	1	50 mM Tris pH 7.5, 10 mM MgCl2, 1 mM DTT	217	23969.4185	-0.0175115	1	0.31797235	0.18894009	0.26728111	Sarah Williams, Gestwicki Lab, UCSF	P28061	3IPM_2	Opoku-Nsiah, Kwadwo A., et al. 2021
RBD	SARS-CoV2 spike protein receptor-binding domain		Facilitates viral entry into host cells	1	50 mM Tris HCl, 150 mM NaCl, 1% DMSO, 0.01% Triton X-100, 1 mM TCEP	223	25098.1047	-0.2587444	0	0.34080717	0.29596413	0.14798206	Commercial - R&D Systems R&D 10500-CV-100	POD01	7CH5	
Retinoic_acid_receptor_alpha	Retinoic acid receptor RXR-alpha		Transcription factor	10		210	23633.3325	0.06238095	0	0.34285714	0.17619048	0.35238095	Mericka McCabe, Cuervo Lab, Albert Einstein College of Medicine	P19793	5EC9	Anguiano, Jaime, et al. 2013
snf2h	SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily A member 5		Helicase, chromatin remodeling	1	15 mM HEPES, pH 7.5, 140 mM KCl, 1 mM Mg-ADP-Befx	1052	121903.937	-0.7191065	25	0.2756654	0.20627376	0.27756654	Hayden Saunders, Narlikar Lab, UCSF	O60264	6NE3	
Widom_601_DNA	Widom 601 sequence, nucleosomal DNA	Screened as DNA alone	Classical DNA sequence for in vitro nucleosome assembly	1	TAE								Hayden Saunders, Narlikar Lab, UCSF			Gamarra, Nathan et al. 2021

Alteen, Matthew G., et al. "A direct fluorescent activity assay for glycosyltransferases enables convenient high-throughput screening: Application to O-GlcNAc Transferase." *Angewandte Chemie International Edition* 59.24 (2020): 9601-9609.

Anguiano, Jaime, et al. "Chemical modulation of chaperone-mediated autophagy by retinoic acid derivatives." *Nature chemical biology* 9.6 (2013): 374-382.

Gamarra, Nathan, and Geeta J. Narlikar. "Histone dynamics play a critical role in SNF2h-mediated nucleosome sliding." *Nature structural & molecular biology* 28.7 (2021): 548-551.

Gustafson, Chelsea L., et al. "A slow conformational switch in the BMAL1 transactivation domain modulates circadian rhythms." *Molecular cell* 66.4 (2017): 447-457.

Michael, Alicia K., et al. "Formation of a repressive complex in the mammalian circadian clock is mediated by the secondary pocket of CRY1." *Proceedings of the National Academy of Sciences* 114.7 (2017): 1560-1565.

Militi, Stefania, et al. "Early doors (Edo) mutant mouse reveals the importance of period 2 (PER2) PAS domain structure for circadian pacemaking." *Proceedings of the National Academy of Sciences* 113.10 (2016): 2756-2761.

Opoku-Nsiah, Kwadwo A., et al. "The YΦ motif defines the structure-activity relationships of human 20S proteasome activators." *Nature communications* 13.1 (2022): 1-12.

Rajagopal, Ponni, et al. "A conserved histidine modulates HSPB5 structure to trigger chaperone activity in response to stress-related acidosis." *Elife* 4 (2015): e07304.

Ravalin M, Theofilas P, Basu K, Opoku-Nsiah KA, Assimon VA, Medina-Cleghorn D, Chen YF, Bohn MF, Arkin M, Grinberg LT, Craik CS, Gestwicki JE. Specificity for latent C termini links the E3 ubiquitin ligase CHIP to caspases. *Nat Chem Biol.* 2019 Aug;15(8):786-794. doi: 10.1038/s41589-019-0322-6. Epub 2019 Jul 18. PMID: 31320752; PMCID: PMC6693490.

Thompson, Erin. Mechanism of compensation for disease variants in human glutamine synthetase. University of California, San Francisco, 2020.


Tsai JC, Miller-Vedam LE, Anand AA, Jaishankar P, Nguyen HC, Renslo AR, Frost A, Walter P. Structure of the nucleotide exchange factor eIF2B reveals mechanism of memory-enhancing molecule. *Science.* 2018 Mar 30;359(6383):eaag0939. doi: 10.1126

Tseng, Roger, et al. "Structural basis of the day-night transition in a bacterial circadian clock." *Science* 355.6330 (2017): 1174-1180.

Wasserman, Douglas R., et al. "Tissue-Restricted Inhibition of mTOR Using Chemical Genetics." *bioRxiv* (2022).

Table S2. Protein panel and biochemical properties. A complete list of all proteins which underwent dye screening for this manuscript. Dye screening conditions and calculated biochemical parameters are included, as well as construct specifications, PDB codes, protein donors and construct-specific citations where applicable. Further information, including primary sequence and amino acid composition, are available in Data S3.

paDSF dye screening procedure



Conducting a paDSF dye screen is nearly identical to a typical DSF experiment, with two modifications: (1) a paDSF dye is used in place of the commercial dye (e.g. SYPRO Orange), and (2) during heating, fluorescence readings are taken in multiple wavelengths. Analyzing data from a paDSF dye screen is different from analyzing a typical DSF experiment, because the goal is to identify which dyes produce a clear, reliable melting curve for the protein.

A general protocol for a paDSF screen testing 96 dyes against a single protein is provided below. Specific guidance for adapting this protocol, such as altering the protein concentration and number of dyes tested, is provided below.

Video protocols of experimental set up are provided in Movie S1.

Starting materials (prepared before dye screen):

- 50 μM purified protein stock
- 50 mL buffer (e.g. 10 mM HEPES pH 7.20, 200 mM NaCl, 1 mM TCEP)
- A 384-well plate, containing 1 μL of 5 mM DMSO stock of a different dye individual wells
- A white, low-volume 384-well qPCR plate
- An optically-clear plate seal

Experimental set-up

1. Dilute protein stock to 5X final screening concentration by adding 21 μL of 50 μM protein stock to 189 μL buffer and mixing thoroughly, yielding 210 μL of 5 μM protein***.
*Calculation: 2 μL of 5X protein is required for each dye tested. 2 μL per dye * 96 dyes = 192 μL protein solution needed total. Make 210 μL of protein solution, to accommodate for sample loss during pipetting. Standard final protein concentration in paDSF dye screens is 1 μM , so make 210 μL of 5 μM protein stock.*
2. Dilute dye stocks to 1.25X final screening concentration by adding 79 μL of buffer to each well of the pre-prepared dye-containing 384-well plate. Note: only 8 μL of 1.25X dye solution is required, so if nanoliter-volume liquid handling is possible, this volume can be reduced, e.g., by adding 20 μL of buffer to 250 nL of 5 mM dye.
*Calculation: 8 μL of 1.25X dye is required for each protein tested. 1.25X * 50 μM final concentration = 62.5 μM working dye concentration. 1 μL is the smallest volume that can be accurately pipetted, so 1 μL starting volume * (5 mM stock concentration / 62.5 μM desired concentration) = 80 μL resuspended volume per dye.*
3. Add 2 μL of 5X protein to 96 wells of the qPCR plate***, For the no-protein control, add 2 μL of buffer to the following 96 wells. Spin the plate down to settle droplets to the bottom of the well.
4. Transfer 8 μL of each 1.25X dye solution from the resuspended dye plate to the protein- and buffer-containing qPCR plate. Final well volumes are 10 μL , and each dye is present in two wells of the plate: one well containing 50 μM dye and 1 μM protein, and the other well containing 50 μM dye and buffer. Do not skip the no-protein control; it is crucial for accurate hit calling.
5. Seal the qPCR plate*** in optically-clear sealing film***, centrifuge briefly (e.g. 30 seconds at 1000 rcf) to remove bubbles.
6. Place the sealed qPCR plate in the qPCR instrument ***.
7. In the qPCR instrument, heat the plate from 25 $^{\circ}\text{C}$ to 95 $^{\circ}\text{C}$ at a rate of 1 $^{\circ}\text{C}$ per minute***. At every degree, measure the fluorescence in FAM, JOE, TAMRA, ROX, Cy5, and Cy5.5 channels***.

*** See the following page for specific recommendations on adapting and modifying these conditions.

paDSF dye screening procedure continued: hit validation

A general protocol for the validation of 10 paDSF hits against a single protein is provided below.

Starting materials (prepared before dye screen):

- 50 μ M purified protein stock
- 50 mL buffer (e.g. 10 mM HEPES pH 7.20, 200 mM NaCl, 1 mM TCEP)
- A 384-well plate, containing 1 μ L of 5, 2.5, 1.25, and 0.63 mM DMSO stock of each hit dye in individual wells
- A white, low-volume 384-well qPCR plate
- An optically-clear plate seal

Experimental set-up

1. Dilute protein stock to 5X final screening concentration by adding 10 μ L of 50 μ M protein stock to 90 μ L buffer and mixing thoroughly, yielding 100 μ L of 5 μ M protein***.

*Calculation: 2 μ L of 5X protein is required for each concentration of each hit dye tested. 2 μ L per dye * 10 dyes * 4 concentrations per dye = 80 μ L protein solution needed total. Make 100 μ L of protein solution, to accommodate for sample loss during pipetting. The same protein concentrations and volumes should be used for screening and validation.*

2. Dilute dye stocks to 1.25X final screening concentration by adding 79 μ L of buffer to each well of the pre-prepared dye-containing 384-well plate. For validation, this plate should contain four concentrations of each hit dye, at decreasing two-fold serial dilutions in DMSO.

Calculation: if a dye was screened at 50 μ M, validation should test 50, 25, 12.5, and 6.3 μ M dye. As before, 8 μ L of 1.25X dye is required for each protein tested.

3. Proceed with steps 3 - 7 of the dye screening protocol, unmodified.

*** See the following page for specific recommendations on adapting and modifying these conditions.

paDSF dye screening procedure continued: adapting the general procedure

Adapting protein, buffer, and dyes

1 of 3 | Protein

Purpose: combines with paDSF dyes to find matches

Protein requirements for paDSF dye screens

In the standard procedure, each dye tested requires 10 μL of 1 μM protein (prepared first as a 5 μM 5X solution). In general, however much protein you would like to use to create a paDSF assay, we recommend using about half of that for the screen. The other half of the protein sample can then be used for hit validation and replicates of the final chosen conditions.

Recommended protein concentration

- **Standard: 1 μM**
- Lowest: 250 nM
- Highest: 20 μM

Comment: like in standard DSF, protein concentrations in paDSF are flexible. Protein concentrations below 1 μM can be used for large proteins or complexes (>80 kDa or so), or to conserve precious samples. Protein concentrations above 10 μM may increase low hit rates for small proteins (<12 kDa or so), or produce marginally smoother data if protein is abundant. If the protein concentration is below 2 μM and screening data have only noisy or low-signal hits, increasing protein concentration is often the easiest fix. Increasing well volumes is less effective in our hands, and screening dyes at > 50 μM may introduce artifacts.

Recommended well volume

- **Standard: 10 μL** in low-volume 384-well plates
- Lowest: plate manufacturers minimum well volume
- Highest: plate manufacturers maximum well volume

Comment: like in standard DSF, well volumes in paDSF are flexible. Well volumes have rarely, if ever, majorly influenced results in our hands.

Recommended screen size

- **Standard: 96 dyes (1000 μL of 1 μM protein)**
- Lowest: ~10 dyes
- Full library: 3 mL of 1 μM protein

Comment: If protein is precious, fewer dyes can be tested. Most proteins get at least a couple of matches in the highest-matching 10 - 25 dyes. The 96 dyes tested in the concise sub-library have the highest protein match-rates, and most proteins get multiple matches with this library. Testing the full library (~300 dyes) will maximize the number of matches. If you want counter-selectivity against another protein (e.g. detection of one protein but not another), screening the full library can help, because the more protein-selective dyes are not included in the standard, highest-matching (and therefore, most protein-general) 96 dyes.

2 of 3 | 50 mL buffer

Purpose: resuspends protein and dyes to appropriate concentrations. Also used in place of protein for protein-free negative controls.

Dye screens are broadly buffer compatible. We recommend screening in whatever buffer you like to use for your system. Standard additives such as DMSO, Triton, EDTA, and reducing agents are fine. High concentrations of glycerol (>5%) are tolerated, but not preferred, as viscous solutions can increase background fluorescence of many dyes.

3 of 3 | Dyes

Purpose: combines with protein to find paDSF matches

To ease assembly, the dye library was designed to be largely commercially available. See Supporting Data Table 1 for detailed information on dye sources, how to assemble the dye library, and recommended subsets of the library to assemble. See Supporting Methods 1 for detailed information on dye handling and storage.

- **Standard: 1 μL of 5 mM dye in DMSO**, pre-dispensed into a well of a 384-well plate.
- If nL liquid handling possible: 250 nL of 5 mM dye in DMSO, pre-dispensed into a well of a 384-well plate.

paDSF dye screening procedure continued: adapting the general procedure

Adapting consumables, instruments, measurements

White 384-well qPCR plate Contains the final DSF experiment, run in the qTower

DSF can be plate-sensitive, so substitutions for qPCR plates are discouraged. We use Axygen PCR-384-LC480-W-NF, and have very few compatibility issues with this product code.

Optically-clear plate seals Seals the qPCR plate during the run

Any kind of optically-clear sealing film should work, so substitutions for sealing film are fine. We use and like Applied Biosystems MicroAmp Optical Adhesive Film ref 4311971.

qPCR instrument Used to run the paDSF experiment

Any qPCR instrument capable of running DSF can be used to run a paDSF dye screen. If possible, use an instrument capable of monitoring in more than one set of fluorescent wavelengths at a time. We use the Analytik Jena qTower 384G, and monitor in FAM, JOE, TAMRA, ROX, Cy5, and Cy5.5, but use of this instrument and precise channels is not necessary. Note that measuring in more wavelengths can, on some qPCR instruments, change increase the run-time, so testing the heating and measuring protocol prior to the experiment, and keeping the measured channels the same across experiments, can reduce unexpected changes to heating rate.

Like typical DSF, the heating procedure can be modified for the chosen protein. paDSF dyes do not have specific requirements for heating and measurement.

Optional: web-connected computer Used to analyze paDSF dye screen data

paDSF dye screens can produce a large amount of data. You can use the provided web app to assist data analysis. Detailed information on hit calling is described in Supplementary Note 1.