Supplementary Information

Materials and Methods

Materials. All affinity reagents used in this work are listed in the Supplementary Information (Table S2). Conjugation reagents, paramagnetic beads, and assay buffers were obtained from Quanterix Corporation. DNA oligos used in the MOSAIC assay were obtained from Integrated DNA Technologies. Antibodies used in final Simoa and MOSAIC assays (monoclonals Ab6, Ab54, 62H12, 34H7) were additionally validated by western blotting (**Figure S22**).

Preparation of capture and detector reagents. All capture antibodies and nanobodies were obtained in or dialyzed into phosphate buffered saline (PBS). For the first-generation Simoa assay, 7×10⁸ carboxvlated paramagnetic 2.7-um beads (Homebrew Singleplex Beads, Quanterix Corp.) were first washed three times with 400 μ L Bead Wash Buffer (Quanterix Corp.) and two times with 400 μ L cold Bead Conjugation Buffer (Quanterix Corp.) before being resuspended in 390 µL cold Bead Conjugation Buffer. A 1 mg vial of 1-ethyl-3-(3dimethylaminopropyl) carbodiimide hydrochloride (EDC) (Thermo Fisher Scientific) was then dissolved to 10 mg/mL in cold Bead Conjugation Buffer, and 10 μ L was added to the beads. The beads were shaken for 30 minutes at 4°C to activate the carboxyl groups on the beads, which were then washed once with 400 μ L cold Bead Conjugation Buffer and resuspended in the capture nanobody solution (10 μ g nanobody total), diluted in Bead Conjugation Buffer to a final volume of 400 μ L. The beads were shaken for two hours at 4°C, washed twice with 400 µL Bead Wash Buffer, and resuspended in 400 µL Bead Blocking Buffer (Quanterix Corp.) before shaking at room temperature for 30 minutes to block the beads. After one wash each with 400 µL Bead Wash Buffer and Bead Diluent (Quanterix Corp.), the beads were resuspended in Bead Diluent and stored at 4°C. Beads were counted with a Beckman Counter Z Series Particle Counter before using in assays. For secondgeneration Simoa assays, the following bead coupling conditions were used: 4.2×10^8 starting beads, 300 μ L wash volumes, 6 μ L EDC, and 40 μ g antibody.

For biotinylation of detector antibodies or nanobodies, a 1 mg vial of Sulfo-NHS-LC-LC-biotin was freshly dissolved in 150 μ L water and added at 80-fold molar excess to a 1 mg/mL solution of antibody or nanobody. The reaction mixture was incubated at 30 minutes at room temperature and subsequently purified with an Amicon Ultra-0.5 mL centrifugal filter (50K and 10K cutoffs for antibody and dimeric nanobody, respectively). Five centrifugation cycles of 14,000xg for five minutes were performed, with addition of 450 μ L PBS each cycle. The purified biotinylated detector reagent was recovered by inverting the filter into a new tube and centrifuging at 1000xg for two minutes. Concentration was quantified using a NanoDrop spectrophotometer.

Recombinant ORF1p protein production. ORF1p was prepared as described¹; briefly, codon optimized human ORF1p corresponding to L1RP (L1 insertion in X-linked retinitis pigmentosa locus, GenBank AF148856.1) with N-terminal His6-TEV was expressed in E. Coli, purified by Ni-NTA affinity, eluted, tag cleaved in the presence of RNaseA, and polished by size exclusion in a buffer containing 50 mM HEPES pH 7.8, 500 mM NaCl, 10 mM MgCl2, and 0.5 mM tris(2-carboxyethyl) phosphine (TCEP), resulting in monodisperse trimeric ORF1p bearing an N-terminal glycine scar.

Nanobody generation and screening. Nanobodies were generated essentially as decribed^{2,3} using mass spectrometry/lymphocyte cDNA sequencing to identify antigen-specific nanobody candidates. Briefly, a llama was immunized with monodisperse ORF1p, and serum and bone marrow were isolated. The heavy chain only IgG fraction (VHH) was isolated from serum and bound to a column of immobilized ORF1p. Bound protein was eluted in SDS and sequenced by mass spectrometry, utilizing a library derived from sequencing VHH fragments PCR-amplified from bone marrow-derived plasma cells. Candidate sequences were cloned into an E. coli expression vector with C-terminal His6 tag and expressed in 50 ml cultures in E. coli Arctic Express RP (Agilent) with 0.2 mM IPTG induction at 12°C overnight. Periplasmic extract was generated as follows: pellets were resuspended in 10 ml per L culture TES buffer (200 mM Tris-HCl, pH 8.0, 0.5 mM EDTA, and 500 mM sucrose), 20 ml/L hypotonic lysis buffer added (TES buffer diluted 1:4 with ddH2O), supplemented with 1 mM PMSF, 3 µg / ml Pepstatin A, incubated 45 min at 4°C, and centrifuged at 25,000 x g for 30 min. The supernatant (periplasmic

extract) was bound to ORF1p-conjugated Sepharose, washed 3 times, eluted with SDS at 70°C for 10 min, and periplasmic extract and elution were analyzed by SDS-PAGE to assay expression and yield. ORF1p-binding candidates were purified as below and analyzed by ELISA (**Figure S7**).

Nanobody and multimeric nanobody purification. C-terminally His6-tagged nanobody constructs were expressed and purified essentially as described². Briefly, protein was expressed in E. coli Arctic Express RP (Agilent) with 0.2 mM IPTG induction at 12°C overnight. Periplasmic extract (generated as above) was supplemented with 5 mM MgCl₂, 500 mM NaCl, and 20 mM imidazole, purified by Ni-NTA chromatography, dialyzed into 150 mM NaCl, 10 mM HEPES, pH 7.4, and concentrated to 1-3 mg/ml by ultrafiltration. "5xCys tail" constructs were purified with the addition of 5 mM TCEP-HCl in resuspension, wash, elution, and dialysis buffers.

Surface plasmon resonance (SPR) assays. Binding kinetics (k_a , k_d , and K_D) of antibody and nanobody constructs for ORF1p were obtained on a Biacore 8K instrument (Cytiva). Recombinant ORF1p was immobilized on a Series S CM5 sensor chip at 1.5 µg/ml using EDC/NHS coupling chemistry according to the manufacturer's guidelines. Nanobodies and antibodies were prepared as analytes and run in buffer containing 20 mM HEPES pH 7.4, 150 mM NaCl, and 0.05% Tween-20. Analytes were injected at 30 µl/min in single-cycle kinetics experiments at concentrations of 0.1, 0.3, 1, 3.3, and 10 nM, with association times of 120-180 sec, and a dissociation time of 1200-7200 sec, depending on observed off-rate. Residual bound protein was removed between experiments using 10 mM glycine-HCl pH 3.0. Data were analyzed using Biacore software, fitting a Langmuir 1:1 binding model to sensorgrams to calculate kinetic parameters.

For epitope binning, pairs of antibodies were sequentially flowed over immobilized ORF1p using Biacore tandem dual injections according to the manufacturer's guidelines. Antibodies were injected at concentrations of 200 nM with a flow rate of 10 μ l/min. Contact time for the first antibody was 120 sec, followed by 150 sec for the second antibody, then a 30 sec dissociation. Response signal for the second antibody was measured in a 10 sec window at the beginning of dissociation. The chip was regenerated between experiments with glycine pH 3.0 as above. Data were analyzed using the Biacore software epitope binning module.

ORF1p Simoa assays. Simoa assays were performed on an HD-X Analyzer (Quanterix Corp.), with all assay reagents and consumables loaded onto the instrument according to the manufacturer's instructions. 250,000 capture beads and 250,000 helper (non-conjugated) beads were used in each Simoa assay. A three-step assay configuration was used for the first- and second-generation assays, consisting of a 15-minute target capture step (incubation of capture beads with 100 μ L sample), 5-minute incubation with detector reagent (0.3 μ g/mL for both first- and second-generation assays), and 5-minute incubation with streptavidin- β -galactosidase (150 pM for first-generation assay; 300 pM for second-generation assays). The beads were washed with System Wash Buffer 1 (Quanterix Corp.) after each assay step. Upon the final wash cycle, the beads were loaded together with the fluorogenic enzyme substrate resorufin β -D-galactopyranoside into a 216,000-microwell array, which was subsequently sealed with oil. Automated imaging and counting of "on" and "off" wells and calculation of average enzyme per bead (AEB) were performed by the instrument. Calibration curves were fit using a 4PL fit with a 1/y² weighting factor, and the limit of detection (LOD) was determined as three standard deviations above the blank.

All plasma and serum samples were diluted four-fold in Homebrew Sample Diluent (Quanterix Corp.) with 1x Halt Protease Inhibitor Cocktail (ThermoFisher), with an additional 1% Triton-X 100 added in the second-generation assays. All recombinant ORF1p calibrators were run in triplicates, with four replicates for the blank calibrator, and all plasma and serum samples were run in duplicates. The average LOD across all sample runs was determined for each assay and depicted in each figure.

Healthy individual plasma and serum samples were obtained from the Mass General Brigham Biobank, with additional samples from the Penn Ovarian Cancer Research Center and Tomas Mustelin (University of Washington). Additional breakdown of patients within each cancer type, by demographic and clinicopathological variables, where available, is included in **Figures S2-S3**, and **S19-21**, and **Table S10**.

ORF1p large-volume MOSAIC assays. MOSAIC assays were performed as previously described, using Eppendorf tubes for the initial capture step. [insert ACS Nano MOSAIC ref] For each sample, 500 µL plasma

was diluted four-fold in Homebrew Sample Diluent with protease inhibitor and 1% Triton-X 100 to a total volume of 2 mL. Briefly, 100,000 capture beads were incubated with sample and mixed for two hours at room temperature, followed by magnetic separation and resuspended in 250 µL System Wash Buffer 1 before transferring to a 96-well plate. The beads were then washed with System Wash Buffer 1 using a Biotek 405 TS Microplate Washer before adding 100 μ L nanobody detector reagent (0.3 μ g/mL, diluted in Homebrew Sample Diluent) and shaking the plate for 10 minutes at room temperature. After washing with the microplate washer, the beads were incubated with 100 μ L streptavidin-DNA (100 pM, diluted in Homebrew Sample Diluent with 5 mM EDTA and 0.02 mg/mL heparin) with shaking for 10 minutes at room temperature, followed by another washing step. The beads were transferred to a new 96-well plate, manually washed with 180 uL System Wash Buffer 1, and resuspended in 50 μ L reaction mixture for rolling circle amplification (RCA). The RCA reaction mixture consisted of 0.33 U/uL phi29 polymerase, 1 nM ATTO647N-labeled DNA probe, 0.5 mM deoxyribonucleotide mix, 0.2 mg/mL bovine serum albumin, and 0.1% Tween-20 in 50 mM Tris-HCl (pH 7.5), 10 mM (NH₄)₂SO₄, and 10 mM MgCl₂. The beads were shaken at 37°C for one hour, followed by addition of 160 μ L PBS with 5 mM EDTA and 0.1% Tween-20. After washing the beads once with 200 μ L of the same buffer, the beads were resuspended in 140 µL buffer with 0.2% BSA. All samples were analyzed using a NovoCyte flow cytometer (Agilent) equipped with three lasers. Analysis of average molecule per bead (AMB) values was performed as previously described using FlowJo software (BD Biosciences) and Python. All code used for MOSAIC data analysis can be downloaded as part of the waltlabtools.mosaic Python module, which is available at https://github.com/tylerdougan/waltlabtools.

Classification models. Classification models were trained for (1) all healthy and all ovarian cancer patients measured by the second-generation assays; and (2) the subset of 51 ovarian cancer and 50 age-matched healthy female patients, obtained from Ronny Drapkin (University of Pennsylvania). Each dataset contained no missing values, and the measurements in the datasets were log-transformed and normalized beforehand for classification analysis of healthy and ovarian cancer subjects. Logistic regression was used for the univariate classifier and the k-nearest neighbors (KNN) and light gradient-boosting machine (LightGBM), which had the best performances among the classifiers, were used for the multivariate classifier, and implemented in Python 3.7.15 with scikit-learn version 1.0.2 package. Each classifier was given a weight optimization between classes to deal with data imbalance between healthy and cancer subjects, as well as hyperparameter tuning using grid search.

The performance of each biomarker in differentiating ovarian cancer subjects from healthy subjects was evaluated with fivefold cross validation by calculating accuracy, precision, recall, f1-value, sensitivity, specificity, and area under the receiver-operating characteristic (ROC) curve (AUC). A stratified five-fold cross-validation strategy randomly splits the positive and negative samples into five equally sized subsets. One positive subset and one negative subset were selected as the test dataset each time, and the other samples were used to train a classification model.

In the multivariate analysis, the Variance Inflation Factor (VIF) for the biomarkers was calculated, and any biomarkers with extremely high correlation with VIF greater than 10 were excluded from the classification model in advance.

Barrett's esophagus cases. A cohort of 75 esophageal biopsies with BE and varying degrees of dysplasia were assembled. Negative cases were screened to have no prior history of dysplasia. The mean age of the cohort was 67 years with a male predominance (M:F ratio = 3.7:1). All samples were re-analyzed for histological features of dysplasia by three experienced gastrointestinal pathologists (LRZ, VD, OHY) who were blinded to the original diagnosis. A consensus was reached for 72 cases and the consensus diagnosis was used as the gold standard. There was moderate agreement between pathologists (kappa 0.43-0.51).

Colon cancer tissue microarray. 178 sequential CRCs resected by a single surgeon from 2011-2013 were assembled on a 3 mm core tissue microarray. All cases were independently scored by two pathologists. The mean age of the cohort was 65 years with 49.8% males. Mean follow-up was 25 months. At resection, 23% were stage I, 33% were stage II, 44% were stage III, and 1% were stage IV.

Ovarian Cancer Samples. Age-matched ovarian cancer (n=53) and healthy control (n=50) patient plasma samples were from University of Pennsylvania Ovarian Cancer Research Center, OCRC Tumor BioTrust Collection, Research Resource Identifier (RRID): SCR_02287.

Gastroesophageal cancer treatment cohort. Nineteen patients received systemic therapy, 3 of which also underwent surgical resection. Patients were treated with concurrent chemotherapy (carboplatin/taxol) and radiation (N=3), fluorouracil/ leucovorin/ oxaliplatin/ docetaxel (FLOT, N=2), fluorouracil/ leucovorin/ irinotecan/ oxaliplatin (FOLFIRINOX, N=2), fluorouracil/ leucovorin/ oxaliplatin (FOLFOX, N=9), FOLFOX + trastuzumab (N=1), pembrolizumab (N=1) or FOLFOX then chemoradiation (1). The mean age of the cohort was 76 years. All patients were male (100%). Fifty-eight percent had locally advanced disease (stage II-III) and 42% had advanced disease (stage IV) at the time of initial diagnosis. Sixty-eight percent (N=13) were deemed responders to therapy while 32% (N=6) were deemed non-responders to standard therapy on review of re-staging imaging (CT and/or PET-CT) by investigators blinded to the assay results.

Histochemistry: ORF1p immunohistochemistry was performed as described using anti-ORF1 4H1 (Millipore)⁴ diluted 1:3000 on a Leica Bond system. Cases were scored by three experienced gastrointestinal pathologists (MST, VD, OHY). **LINE-1 in situ hybridization** was performed as described using RNAscope catalog 565098 (Advanced Cell Diagnostics) on a Leica Bond system⁵. The probe is complementary to the 5' end of L1RP (L1 insertion in X-linked retinitis pigmentosa locus). Cases were scored by three experienced gastrointestinal pathologists (MST, VD, OHY).



Supplementary Figure 1. Age and gender distributions of healthy control patients used in the first- and second-generation ORF1p Simoa assays from the (A-B) Mass General Brigham Biobank in the first-generation (A) and second-generation (B) ORF1p Simoa assays; and (C) Penn Medicine Biobank. Healthy female control patients from the Penn Medicine Biobank were age-matched with the ovarian cancer patients from the Penn Ovarian Cancer Research Center Tumor BioTrust Collection.



Supplementary Figure 2. Circulating ORF1p levels in eight stage I ovarian cancer patients, as measured by the first-generation Simoa assay. Blue dashed line indicates the assay limit of detection, accounting for the four-fold dilution factor.



Supplementary Figure 3. Circulating ORF1p levels in lung cancer patient cohort, as classified by (A) disease subtype, (B) gender, and (C) smoking status. Blue dashed line indicates the assay limit of detection, accounting for the four-fold dilution factor.



Supplementary Figure 4. Circulating ORF1p levels in patients with IPMN, measured using the first-generation Simoa assay. IPMN, intraductal papillary mucinous neoplasm. Blue dashed line denotes the assay limit of detection, accounting for four-fold dilution.



Supplementary Figure 5. Circulating ORF1p levels in patients with chronic liver disease or systemic lupus erythematous, measured using the first-generation Simoa assay. The red dot represents a chronic liver disease patient diagnosed with hepatocellular carcinoma after the time of sampling. Blue dashed line denotes the assay limit of detection, accounting for four-fold dilution.



Supplementary Figure 6. Analysis of circulating ORF1p levels in size exclusion chromatography (SEC) fractions of colorectal cancer patient plasma. 500 µL plasma was filtered through a 0.45 µm centrifugal filter and fractionated with a Sepharose CL-6B resin packed column (A) CD9 and ORF1p levels in each SEC fraction as measured by Simoa. Blue highlighted boxes denote extracellular vesicle (EV)-containing fractions. The majority of circulating ORF1p is measured in free protein fractions. (B) Determination of ORF1p levels inside EVs via proteinase K (protK) protection assays. EV-containing fractions were pooled and concentrated and treated with proteinase K, followed by the serine protease inhibitor phenylmethylsulfonyl fluoride (PMSF) and 0.5% Triton-X. Controls without proteinase K digestion (no protK) and with simultaneous proteinase K and Triton-X treatment (protK/Triton-X) were performed. Decreased ORF1p concentrations were measured after proteinase K digestion and subsequent EV lysis by Triton-X, further suggesting that only a very small amount of circulating ORF1p is in EVs.

Α																	В		ELISA	
Candid	ate Nb:	: 1		2	2	3		4		5		6	;	7	Ċ				EC50 (nM)	
		1	F	<u> </u>	F	1	F	1	F	1	F	1	F	1	F			Nb ORF1-1	19.4	
<u>kDa</u>		-	-	-	-	-	-	-	-	-	-	-	-		-			Nb ORF1-2	4.8	
198	-																	Nb ORF1-3	3.0	
								-		-								Nb ORF1-4	6.0	
98																		Nb ORF1-5	1.8	
	1																	Nb ORF1-6	11.4	
62	-	-			-	-	-	-	-	-		\ge		-				Nb ORF1-7	NA	
																		Nb ORF1-8	1.4	
49		-	_	a second	-		-	9	-	8	-	1	-	,	-			Nb ORF1-9	1.2	
38	-		-		=		_		_									Nb ORF1-10	1.2	
										2								Nb ORF1-11	1.5	
28			-	Ξ	-	\equiv	-		-	Ξ	-	18	-		-			Nb ORF1-12	1.4	
																		Nb ORF1-13	1.1	
17	1.0			-	-					-						— Nan	obody	Nb ORF1-14	NA	
14			-		-	-	-	-			-	-	۰					Nb ORF1-15	0.6	
																		Nb ORF1-16	1.0	
6	100																	Nb ORF1-17	1.0	
	L = lo	ad (peri	plasr	nic e	extra	ct)											Nb ORF1-18	5.6	
	E = El	utio	n (S	DS b	oileo	d bea	ids p	ost ii	ncub	oatio	n)							Nb ORF1-19	7.9	
																		Nb ORF1-20	1.8	
																		Nb ORF1-21	0.6	

Supplementary Figure 7. Nanobody generation. **(**A) Representative screening results for nanobodies 1-7. Periplasmic extracts from E. Coli expressing the indicated candidate nanobodies (Nb) were bound to ORF1p-conjugated beads and eluted, then analyzed by SDS-PAGE and stained by Coomassie blue. Expressed and purified nanobodies (~15 kDa) are indicated; larger/darker bands indicate higher levels of expression and/or purification. (B) EC50s from ELISAs performed against recombinant ORF1p. Highlighted nanobodies were selected for follow-up characterization, while remainder were rejected for poor affinity and/or protein expression.







Supplementary Figure 8. Representative SPR sensorgrams of ORF1 nanobodies. Single-cycle kinetics performed over immobilized ORF1, with nanobodies injected sequentially at 0.1, 0.3, 1, 3, and 10nM concentrations. Raw data (blue) and a 1:1 binding model fit (black) are plotted.

Α

Nanobody Constructs

Construct	Description				
MT991	Nb 1-1				
MT992	Nb 2-2				
MT993	Nb 5-5				
MT994	Nb 5-1				
MT995	Nb 5-2				
MT996	Nb 5-5-5				
MT997	Nb 5-5 long "LL"				







Supplementary Figure 9. Engineered Nanobody constructs. (A) Schematic of design of engineered dimeric and trimeric nanobody constructs, with flexible (GGGGS \times 4) and rigid helical (EAAAK \times 3) linkers. (B), Representative preparation of engineered nanobody constructs, Coomassie stain.





1	\mathbf{r}	

					D	ETECTO	R			
		Ab6	Nb5-5LL	33A8	34C5	34H7	64C6	61A11	62H12	55A6
	33A8	0.016								
	34C5	0.026				0.032				
	34H7	0.006	0.043							
	55A6	0.026	0.055							
ш	64C6	0.034	0.041							
- n	36D12	0.046	0.050							
APT	42D10	0.021	0.020							
U U	61A11	0.042	0.066							
	62H12	0.008	0.062							
	Nb5-5	0.050		0.133	0.022	0.047	0.038	0.042	0.120	
	4H1	0.030	0.056	0.057		0.044				0.036
	Ab6	0.102	0.044	0.063					0.113	

Limit of Detection (pg/mL)

Supplementary Figure 10. Screening of newly developed monoclonal antibody candidates (GenScript). (A) Signal:background ratios of novel rabbit monoclonal α-ORF1p antibodies (B-cell supernatants) in a modified Simoa employing the candidate mAb plus biotinylated secondary anti-Rabbit pAb as detector; these show up to 5-fold improvement vs. our prior best detection antibody, Ab6. Two different capture beads with distinct epitopes were employed. (B-C) Best performers were then synthesized and purified; mAbs were used in further screening with a dimeric nanobody and commercially available monoclonal antibodies for ORF1p detection with Simoa. (B) Signal-to-background comparisons of affinity reagents as capture/detector pairs on Simoa, using recombinant ORF1p protein. All labeled affinity reagents except Nb5-5LL are monoclonal antibodies; Nb5-5LL denotes a homodimeric form of the nanobody Nb5. (C) Limit of detection values for affinity reagent pairs selected from screening in (B). Anti-ORF1 D3W9O is from Cell Signaling Technology.

bioRxiv preprint doi: https://doi.org/10.1101/2023.01.25.525462; this version posted January 25, 2023. The copyright holder for this preprint (which was not certified by peer review) is the author/funder, who has granted bioRxiv a license to display the preprint in perpetuity. It is made available under aCC-BY-ND 4.0 International license.



Supplementary Figure 11. Representative SPR sensorgrams of ORF1 monoclonal antibodies and nanobody multimers. Single-cycle kinetics performed over immobilized ORF1, with (A) monoclonal antibodies or (B) nanobody multimers injected sequentially at 0.1, 0.3, 1, 3, and 10nM concentrations. Raw data (blue) and a 1:1 binding model fit (black) are plotted.



Supplementary Figure 12. Screening of multimeric nanobodies for ORF1p detection with Simoa. (A) Signal-tobackground (Signal:Bg) comparisons of multimeric nanobodies as detector reagents on the Simoa platform, using recombinant ORF1p protein. A four-plex assay format, with a unique fluorescent dye-encoded bead type for each of the four capture reagents tested, was used for screening. (B) Signal-to-background comparison of multimeric nanobodies as capture reagents on Simoa, using recombinant ORF1p protein, paired with the monoclonal antibody Ab5 as detector reagent. Sub-labels "A" and "B" refer to different nanobody concentrations used during conjugation to beads. (C-D) Screening of select affinity reagent pairs in small sets of plasma samples from healthy and cancer (colorectal and gastroesophageal) patients. Each assay is denoted by the capture/detector reagent pair. The first-generation assay (Nb5/Ab6) measurements are depicted for comparison. All assays were performed as three-step Simoa assays unless otherwise indicated (two-step). Blue dashed lines indicate the assay limits of detection, accounting for four-fold dilution.



Supplementary Figure 13. First round of screening of newly developed monoclonal antibody and dimeric nanobody reagent pairs in patient plasma. Affinity reagent pairs selected from initial screening with recombinant ORF1 protein (Supplementary Figure 8) were screened in small sets of plasma samples from healthy and cancer (colorectal and gastroesophageal) patients. Eight healthy and eight cancer patients were used in each set of assays in (A) and (B). Blue dashed lines indicate assay limits of detection, accounting for four-fold dilution. Each assay is denoted by the capture/detector reagent pair. The first-generation assay (Nb5/Ab6) measurements are depicted for comparison (gray box).



Supplementary Figure 14. Second round of screening of newly developed monoclonal antibody and dimeric nanobody reagent pairs in patient plasma. Affinity reagent pairs selected from a first round of screening in plasma samples (Supplementary Figure 9) were screened in 25 healthy and 25 cancer (colorectal, gastroesophageal, and breast) patient plasma samples. Each assay is denoted by the capture/detector reagent pair. The first-generation assay (Nb5/Ab6) measurements are depicted for comparison (gray box). Blue dashed lines indicate assay limits of detection, accounting for four-fold dilution. Blue boxes indicate the final three selected second-generation assays.



Supplementary Figure 15. Assay validation for second-generation ORF1p Simoa assays. (A-C) Representative calibration curve (A), dilution linearity (B), and spike and recovery results (C) for the 34H7/Nb5-5LL S16

(capture/detector) assay. (D-F) Representative calibration curve (D), dilution linearity (E), and spike and recovery results (F) for the 62H12/Nb5-5LL (capture/detector) assay. (G-I) Representative calibration curve (G), dilution linearity (H), and spike and recovery results (I) for the 62H12/Ab6 (capture/detector) assay. Error bars represent the standard deviation of triplicate measurements in each calibration curve and duplicate measurements of all plasma samples.



Supplementary Figure 16. Additional second-generation ORF1p Simoa assay measurements in healthy and cancer patients. Each assay is denoted by the capture/detector reagent pair. Blue dashed lines indicate assay limits of detection, accounting for four-fold dilution. The proportion of detectable patients within each cancer type is depicted by corresponding pie charts and early-stage patients (stage I/II) are indicated as red and orange dots, respectively.



Supplementary Figure 17. Circulating CA125 and HE4 levels, as measured by Simoa, in ovarian cancer patient cohort (obtained from Ronny Drapkin at the University of Pennsylvania). Early-stage patients (stage I/II) are indicated as red and orange dots, respectively.



Supplementary Figure 18. Calibration curve for the large-volume MOSAIC assay used in Figure 5. Blue dashed line indicates the assay limit of detection.



Supplementary Figure 19. Circulating ORF1p levels in breast cancer patients, classified by receptor status, with (A) metastatic and (B) localized disease. First- and second-generation Simoa assays were used to measure circulating ORF1p, with each assay labeled as capture/detector reagent. Blue dashed lines indicate assay limits of detection, accounting for four-fold dilution. ER, estrogen receptor; PR, progesterone receptor; HER2, human epidermal growth factor receptor 2; TNBC, triple-negative breast cancer.



Supplementary Figure 20. Circulating ORF1p (first generation assay) versus prostate-specific antigen levels in prostate cancer patient plasma. Blue dashed line indicates the assay limit of detection, accounting for four-fold dilution.



Supplementary Figure 21. Genomics analysis of lung cancer patient cohort. (A) Comparisons of aeuploidy score, amplifications, deletions, and tumor mutational burden between lung cancer patients with detectable versus undetectable circulating ORF1p levels, as measured by the first-generation ORF1p Simoa assay. (B) Mutation types and top mutated genes in lung cancer patients with detectable versus undetectable circulating ORF1p.



Supplementary Figure 22. Validation of monoclonal antibodies by western blotting. HeLa cells with stably integrated tet-On LINE-1 (L1RP)⁶ were treated with DMSO as control or with 1 μ g/ml doxycycline (Dox) for 4 days to induce L1 expression, lysed with RIPA buffer, cleared, and 20 μ g total protein (BCA) was loaded per lane and blotted with the indicated antibodies.

Platform	Assay (Capture / Detector)	Sample Volume (µL)	Incubation times (capture- detector- streptavidin minutes)	Detector (pg/mL)	Streptavidin β- galactosidase (pM)	Limit Detecti (pg/mL	of on)
Simoa	Nb5 / Ab6	100	15-5-5	0.3	150	0.056 0.032	Ħ
Simoa	34H7 / Nb5-5LL	100	15-5-5	0.3	300	0.027 0.016	±
Simoa	62H12 / Nb5-5LL	100	15-5-5	0.3	300	0.029 0.016	±
Simoa	62H12 / Ab6	100	15-5-5	0.3	300	0.016 0.007	±
MOSAIC	34H7 / Nb5-5LL	2000	120-10-10	0.3	100 (streptavidin- DNA)	0.002	

 Table S2. Affinity reagents used across all screening experiments.

Affinity reagent	Туре	Source
Ab6	Rabbit monoclonal antibody	Abcam (ab246317)
Ab54	Rabbit monoclonal antibody	Abcam (ab246320)
Nb5	Nanobody	Rockefeller University
4H1	Mouse monoclonal antibody	MilliporeSigma (MABC1152)
JH73	Rabbit monoclonal antibody	
JH74	Rabbit monoclonal antibody	
D3W9O	Rabbit monoclonal antibody	Cell Signaling Technology (88701)
Nb5-5 pMT993	Homodimeric nanobody (short linker)	This study
Nb5-5LL pMT997	Homodimeric nanobody (long linker)	This study
Nb1-1 pMT991	Homodimeric nanobody (short linker)	This study
Nb2-2 pMT992	Homodimeric nanobody (short linker)	This study
Nb5-1 pMT994	Heterodimeric nanobody (short linker)	This study
Nb5-2 pMT995	Heterodimeric nanobody (short linker)	This study
33A8	Rabbit monoclonal antibody	This study
34C5	Rabbit monoclonal antibody	This study
34H7	Rabbit monoclonal antibody	This study
36D12	Rabbit monoclonal antibody	This study
42D10	Rabbit monoclonal antibody	This study
50E9	Rabbit monoclonal antibody	This study
55A6	Rabbit monoclonal antibody	This study
64C6	Rabbit monoclonal antibody	This study
61A11	Rabbit monoclonal antibody	This study
62H12	Rabbit monoclonal antibody	This study

Construct	ka (1/Ms)	kd (1/s)	KD (M)	KD se	n
NB Orf1-1	6.1E+05	6.4E-04	1.0E-09		1
NB Orf1-2	6.4E+05	4.2E-05	7.1E-11	2.0E-11	3
NB Orf1-5	4.8E+06	2.3E-03	3.7E-10	1.9E-10	3
NB Orf1-8	9.4E+05	2.0E-04	2.0E-10	4.1E-11	3
NB Orf1-9	6.1E+05	1.3E-04	2.2E-10	2.0E-11	3
NB Orf1-10	8.0E+05	1.8E-04	2.7E-10	4.5E-11	3
NB Orf1-11	4.5E+06	3.3E-03	6.4E-10	1.6E-10	3
NB Orf1-12	1.2E+06	4.2E-04	2.9E-10	7.6E-11	3
NB Orf1-13	8.2E+05	2.5E-04	3.0E-10	2.3E-11	3
NB Orf1-15	1.1E+06	5.1E-04	4.9E-10	3.8E-11	3
NB Orf1-16	2.6E+06	3.2E-03	8.4E-10	3.4E-10	3
NB Orf1-17	1.5E+06	4.7E-04	2.9E-10	5.7E-11	3
NB Orf1-20	1.7E+06	3.8E-04	2.4E-10	4.2E-11	3
NB Orf1-21	9.6E+05	1.8E-04	1.7E-10	3.3E-11	3

Table S3. Surface plasmon resonance measurements of binding kinetics of anti-ORF1p nanobodies

Table S4. Surface plasmon resonance measurements of binding kinetics of affinity reagents used in single molecule assays. Orange-highlighted cells indicate dissociation rate constants (k_d) below the limit of detection of the instrument, resulting in upper-bound equilibrium dissociation constant (K_D) estimation. SEM, Standard Error of the Mean.

Construct	#	k _a (1/M·s)	k _d (1/s)	К _D (М)	K _D SEM	n				
Dimeric and Trimeric Nanobodies										
Nb 1-1	MT991	Weak signal								
Nb 2-2	MT992	6.4E+05	<5.0E-06	<7.9E-12	7.4E-13	3				
Nb 5-5	MT993	9.0E+05	3.4E-05	3.8E-11	1.5E-11	3				
Nb 5-1	MT994	5.2E+05	3.1E-05	6.0E-11	3.8E-11	2				
Nb 5-2	MT995	4.0E+05	1.3E-05	3.4E-11	4.7E-12	3				
Nb 5-5-5	MT996	8.1E+05	2.9E-05	3.6E-11	8.5E-12	3				
Nb 5-5 long "LL"	MT997	7.5E+05	2.1E-05	2.8E-11	1.2E-11	3				
Monoclonal Antibodie	S									
mAb	62H12	4.8E+05	<5.0E-06	<1.0E-11	4.4E-13	3				
mAb	34H7	5.2E+05	<5.0E-06	<9.6E-12	4.1E-12	3				
mAb	Ab6	2.2E+05	<5.0E-06	<2.2E-11	2.6E-12	3				
mAb	4H1	4.5E+05	<5.0E-06	<1.1E-11	2.5E-12	3				
mAb	Ab54	4.2E+05	<5.0E-06	<1.2E-11	5.5E-14	2				

Table S5. Epitope binning of ORF1p nanobodies. Fraction of 2nd antibody able to bind after saturation with 1st antibody is indicated, normalized to maximum signal bound. Overlapping epitopes are defined by <0.3 remaining binding activity (red shading). Epitope groups are labeled I-III.

		2nd Ab	2nd Ab								
	Group:	I	I								
	_	Nb-2	Nb-5	Nb-8	Nb-9	Nb-10	Nb-21				
	Nb-2	0.16	0.66	0.94	1.00	1.00	1.00				
	Nb-5	0.82	-0.06	1.00	0.99	0.93	0.90				
	Nb-8	1.00	0.53	0.07	0.13	0.26	0.24				
Ab	Nb-9	0.91	0.66	0.03	80.0	0.19	0.17				
	Nb-10	0.70	1.00	-0.08	-0.07	0.08	0.05				
1st	Nb-21	0.83	0.47	-0.06	-0.05	0.08	0.03				

Table S6. Epitope binning of ORF1p MAbs and representative nanobodies. Fraction of 2nd antibody able to bind after saturation with 1st antibody is indicated, normalized to maximum signal bound. Overlapping epitopes are defined by <0.3 remaining binding activity (red shading). Epitope groups are labeled I-V.

	2nd Ab						
	Group:	I	=	≡	IV	V	
		Nb-2	Nb-5	Nb-9	Ab6	34H7	62H12
	Nb-2	0.18	0.91	0.82	0.97	1.00	0.98
	Nb-5	0.67	0.02	0.83	0.95	0.88	0.99
	Nb-9	0.64	0.90	0.04	0.78	0.42	0.86
0	Ab6	1.00	0.94	0.87	0.07	0.71	1.00
A	34H7	0.87	0.92	0.85	0.81	0.10	0.22
1st	62H12	0.89	1.00	1.00	1.00	0.59	0.12

Table S7. Performance metrics of classification models built from all healthy and all ovarian cancer patients measured by the second-generation Simoa assays. Classification models were built using five-fold cross validation, with univariate logistic regression for each individual assay and k-nearest neighbors algorithm for the multivariate model combining two assays.

Assay	Accuracy	Sensitivity	Specificity	F1 Score
34H7 / Nb5-5LL	95.0%	85.7%	97.5%	0.876
62H12 / Nb5-5LL	94.6%	86.5%	96.9%	0.870
62H12 / Ab6	92.0%	87.3%	93.3%	0.827
34H7/Nb5-5LL and 62H12/Ab6	94.4%	84.1%	97.3%	0.865

Table S8. Performance metrics of classification models built from ovarian cancer and age-matched healthy female patients (Ronny Drapkin, University of Pennsylvania). Classification models were built using k-nearest neighbors (KNN) and light gradient-boosting machine (LightGBM) algorithms, with five-fold cross validation.

Biomarkers	Accuracy	Sensitivity	Specificity	F1 Score
CA125 and HE4	89.0%	90.0%	88.0%	0.891
CA125, HE4, and				
ORF1p	95.0%	94.0%	96.0%	0.949
(34H7/Nb5-5LL)				
CA125, HE4, and				
ORF1p	93.0%	88.2%	98.0%	0.926
(62H12/Nb5-5LL)				
CA125, HE4, and				
ORF1p	91.1%	86.2%	96.0%	0.906
(62H12/Ab6)				

 Table S9. Clinicopathological characteristics for gastroesophageal cancer cohort

Characteristic	N=19		
Median age (range)-yr	76 (37-81%)		
Male sex-no. (%)	19 (100%)		
Histology-no. (%)	Adenocarcinoma	19 (100%)	
	Esophagus	7 (37%)	
Primary tumor location	Gastroesophageal Junction	7 (37%)	
no. (%)	Stomach	5 (26%)	
Disease stage at initial diagnosis-no.	Locally advanced (II- III)	11 (58%)	
(%)	Advanced (IV)	8 (42%)	
Treatment paradigm	Neoadjuvant	9 (47%)	
during pre- and post-	Adjuvant	2 (11%)	
assay collection*	Palliative	8 (42%)	
	Surgery	3 (16%)	
Treatment response	Responder	13 (68%)	
	Non-responder	6 (32%)	

*Patients who underwent surgery also received systemic therapy thus total exceeds N=19

Table S10. Demographic and clinicopathological variables of head and neck squamous cell carcinoma patients grouped by detectability of circulating ORF1p.

	Total	ORF1p-negative	ORF1p-positive
Sex			
Male	33	10	1
Female	11	23	10
Cancer site			
Oral cavity	33	24	9
Oropharynx	3	2	1
Larynx	8	7	1
Race			
White	38	28	10
Asian	2	2	0
Black	1	0	1
Other/Unknown	3	3	0
pT Stage			
pT2	6	4	2
pT3	13	11	2
pT4a	20	15	5
Unknown	5	3	2
pN Stage			
pN0	19	15	4
pN1	8	6	2
pN2	6	6	0
pN3	6	3	3
Unknown	5	3	2
PNI			
Positive	18	14	4
Negative	22	16	6
Unknown	4	3	1
ENE			
Positive	9	6	3
Negative	12	9	3
NA	19	15	4
Unknown	4	3	1

Supplementary References

- 1. Carter, V., et al. High Prevalence and Disease Correlation of Autoantibodies Against p40 Encoded by Long Interspersed Nuclear Elements in Systemic Lupus Erythematosus. *Arthritis Rheumatol* **72**, 89-99 (2020).
- 2. Fridy, P.C., *et al.* A robust pipeline for rapid production of versatile nanobody repertoires. *Nat Methods* **11**, 1253-1260 (2014).
- 3. Mast, F.D., *et al.* Highly synergistic combinations of nanobodies that target SARS-CoV-2 and are resistant to escape. *Elife* **10**(2021).
- 4. Rodić, N., et al. Long Interspersed Element-1 Protein Expression Is a Hallmark of Many Human Cancers. *The American Journal of Pathology* **184**, 1280-1286 (2014).
- 5. Rajurkar, M., et al. Reverse Transcriptase Inhibition Disrupts Repeat Element Life Cycle in Colorectal Cancer. *Cancer Discov* (2022).
- 6. Ardeljan, D., et al. Cell fitness screens reveal a conflict between LINE-1 retrotransposition and DNA replication. *Nat Struct Mol Biol* **27**, 168-178 (2020).